



Human Epidermal Keratinocytes for iPS Cell reprogramming

**Complete Kit (Cat.# PC503hEKTN-K)
Cells Only (Cat.# PC503hEKTN-C)**

User Manual

Store kit at -80°C on receipt

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.

Contents

I. Human Epidermal Keratinocytes neonatal (HEKn).....	2
A. Description.....	2
B. Maintenance of HEKn.....	2
C. Product Information.....	4
D. References.....	5
E. Related Products.....	5
F. Technical Support.....	6
II. Licensing and Warranty Statement.....	7

List of Components

The Human Epidermal Keratinocytes, neonatal (HEKn), are available as either complete culture kit (Cat# PC503hEKTN-K) or cryopreserved cells (Cat# PC503hEKTN-C) in one vial.

HEKn, complete kit		Cat# PC503hEKTN-K
Component	Quantity	Ship
HEKn cells	1x10 ⁶ cells/vial	Cryopreserved and shipped on dry ice
HEKn, supplemented culture medium	250 ml	Store at 4° C upon receipt
Trypsin neutralization solution	50 ml TNS	Store at 4° C upon receipt
HEKn, cryopreserved		Cat# PC503hEKTN-C
HEKn cells	1x10 ⁶ cells/vial	Cryopreserved and shipped on dry ice

The cryopreserved HEKn cells are shipped on dry ice and should be immediately stored in liquid nitrogen upon receipt. Properly stored cells are stable for more than 2 years from the date received.

I. Human Epidermal Keratinocytes, neonatal

A. Description

Human epidermal keratinocyte, neonatal (HEKn) were isolated from individual neonatal foreskin (Day 1 to Day 3 of age). Each vial of this product contains 1×10^6 viable cells. In our laboratory, each lot of cells is performance tested by culturing the cells through multiple passages in the absence of antibiotics and antimycotics, no contamination was observed during this culture period. Upon thawing, the cells are guaranteed to be >75% viable and have a potential of >30 population doublings when cultured according to the instructions provided in this manual.

B. Maintenance of HEKn Cultures

Culture medium

We recommend serum-free medium supplemented with keratinocyte growth supplement. The supplemented medium included in the HEKn complete kit is complete culture medium optimized for HEKn. HEKn culture medium can also be purchased from either Invitrogen (Medium 154CF, Cat# M-154CF-500) or ATCC (Cat# PCS-200-030) and must be supplemented with calcium plus their respective keratinocyte Growth Supplement (Invitrogen: HKGS, Cat# S-001-5; ATCC Cat# PCS-200-040). We recommend a low calcium concentration (0.05 ~ 0.07 mM) to slow HEKn differentiation, the calcium concentration in the supplement medium in our kit is 0.07 mM.

Trypsin/EDTA: 0.25% trypsin-EDTA

Trypsin neutralization solution (TNS): 5% Chelexed Fetal Bovine Serum in Hanks BSS

Initiating Cultures from Cryopreserved Cells

To insure the highest level of viability, be sure to warm medium to 37°C before using it on the cells. We recommend seeding cells recovered from cryopreservation at a cell density of 4×10^3 viable cells/cm² or higher. The procedure given below is a sample protocol for establishing cultures from the contents of one vial.

1. Prepare a bottle of supplemented keratinocyte growth medium according to the instructions supplied with that product.
2. Remove the vial from liquid nitrogen storage, taking care to protect hands and eyes.
3. Dip the lower half of the vial into a 37° C water bath to thaw.
4. When the contents of the vial have thawed, wipe the outside of the vial with 75% alcohol to disinfect and move the vial to a laminar flow culture hood.
5. Open the vial and pipet the cell suspension up and transfer the cells into a 15-cm conical tube with 10 ml fresh keratinocyte growth medium.
6. Pipet up and down with a 10 ml pipette to disperse the cells and centrifuge the cells at 180 g for 5 minutes. Observe the cell pellet.
7. Remove the supernatant from the tube, being careful not to dislodge the cell pellet.
8. Dilute the cells with fresh culture medium and seed new culture vessels with 4×10^3 cells/cm².
9. Incubate the cultures in a 37°C, 5% CO₂ /95% air humidified cell culture incubator.
10. Do not disturb the culture for at least 24 hours after the culture has been initiated.

Maintaining Stock Cultures

1. Change the culture medium 24 to 36 hours after establishing the secondary culture from cryopreserved cells.
2. Change the medium every other day thereafter, until the culture is approximately 80% confluent.

Note: Do not allow stock culture to exceed 80% confluent to prevent differentiation of HEKn.

Subculturing HEKn

View the culture under the microscope to confirm that the cells are subconfluent before splitting. The following protocol is designed for the subculture of one 75 cm² culture flask.


1. Remove all of the culture medium from the flask.
2. Add 3 ml of 0.25% trypsin-EDTA solution to the flask. Rock the flask to ensure that the entire surfaced is covered.
3. Incubate the flask at room temperature until the cells have become completely round, approximately 5-10 minutes. View the culture frequently under a microscope to avoid over digestion.
4. Add 7 ml of trypsin neutralization solution to the flask and transfer the detached cells to a sterile 15 ml conical tube.
5. Centrifuge the cells at 180 g for 5 minutes.
6. Remove the supernatant from the tube, being careful not to dislodge the cell pellet.
7. Resuspend the cells pellet in 10 ml supplemented medium. Pipet the cells up and down with a 10 ml pipette to ensure a homogeneous cell suspension.
8. Determine the concentration of cells in the suspension.
9. Seed new culture vessels with 4×10^3 cells/cm², or 1 to 4 split if starting cells at a 80% confluence.
10. Incubate the cultures in a 37 °C, 5% CO₂ /95% air humidified cell culture incubator.

Cryopreserving HEKn

1. Follow steps 1-6 from the Subculturing of Cells above.
2. Resuspend the cell pellet in supplemented medium. Add approximately 1 ml for each T75 flask.
3. Count the number of cells and dilute the cell suspension to 2×10^6 cells/ml.
4. Add an equal volume of cold 2X Freezing Media to the cell suspension.
5. Aliquot 1 ml of suspension into each cryovial (1×10^6 cells/vial).
6. Place the vials in a cell-freezing container and keep it at -80°C overnight.
7. Transfer the vials to a liquid nitrogen tank for long-term storage.

C. Product Information**Organism:** *Homo sapiens* (human)**Gender:** Male**Age:** Neonatal (Day 1 to Day 3)**Morphology:** Cobblestone appearance**Comments:** The serum-free media is formulated to inhibit fibroblast growth and the low calcium concentration 50~70 μM slows differentiation. No feeder layers or extracellular matrix proteins are required.

HEKn cells grown in supplemented keratinocyte growth medium. Images of human epidermal keratinocytes in culture after thaw.

**Day 2 -HEKn****Day 4 -HEKn****Day 6 -HEKn**

D. References

Aasen T et al. 2008. Efficient and rapid generation of induced pluripotent stem cells from human keratinocytes. *Nature Biotechnology*. 26(11):1276-1284

Jiang YJ et al. 2009. Ceramide stimulates ABCA12 expression via peroxisom proliferator-activated receptor delta in human keratinocytes. *J Biol Chem*. 284:18942-52.

Carey BW et al. 2009. Reprogramming of murine and human somatic cells using a single polycistronic vector. *Proc Natl Acad Sci*. 106:157-62.

Takahashi, K. and Yamanaka, S. 2006. Induction of pluripotent stem cells from mouse embryonic and adult fibroblast cultures by defined factors. *Cell* 126: 663–676.

Takahashi K. et al. 2007. Induction of pluripotent stem cells from adult human fibroblasts by defined factors. *Cell*. 131: 861–72.

Park, IH et al. 2008. Reprogramming of human somatic cells to pluripotency with defined factors. *Nature*. 451:141–6.

E. Related Products

- **Human Foreskin Fibroblasts, neonatal as reprogramming source cell (Cat. # PC501A-HFF)**
From mesoderm origin, these neonatal HFFn are isolated from individual neonatal foreskin and therefore have a single genetic background. High purity and low passage of these neonatal HFFn (Day 1 to Day 3 of age) provides a reliable cell source for efficient reprogramming.
- **Human Foreskin Fibroblasts, neonatal as feeder cells (Cat. # PC502B-HFF).**
These pooled HFFs are optimized to provide a balanced nutrition to support your ES cell culture and reprogramming experiment. High purity and low passage of neonatal HFFn (Day 1 to Day 3 of age) are characteristics of this product.

F. Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site:

<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI)
1616 North Shoreline Blvd.
Mountain View, CA 94043

Phone: (650) 968-2200
(888) 266-5066 (Toll Free)

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com
Technical Support: tech@systembio.com
Ordering Information: orders@systembio.com

I. Licensing and Warranty Statement

Limited Use License

Use of the human epidermal keratinocytes (*i.e.*, the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for stem cell research.

SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2009 System Biosciences (SBI).