



# PrecisionX™ HR Targeting Vectors

Cat #s HRxxxPA-1, PBHR100A-1

## *User Manual*

---

**Store Kits at -20°C upon receipt**

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.

## Contents

I. Introduction.....	2
A. Overview of Genome Engineering Tools.....	2
B. Overview of PrecisionX™ Targeting Vectors for HR-Based Applications .....	5
C. Design Considerations for Homologous Recombination Arms .....	9
II. PrecisionX™ Targeting Vectors By Application: Gene Knock-Out.....	12
A. Overview of Gene Knock-Out .....	12
B. Validation Data for HR110PA-1 .....	17
C. Additional Materials Required .....	21
D. General Considerations .....	22
E. Cloning of Homology Arms into Vector .....	23
F. Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells .....	26
III. PrecisionX™ Targeting Vectors By Application: Gene Knock-In.....	28
A. Overview of Gene Knock-In .....	28
B. Additional Materials Required .....	29
C. General Considerations .....	30
D. Cloning of Homology Arms and Expression Cassette into Vector .....	31
E. Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells .....	30
IV. PrecisionX™ Targeting Vectors By Application: Gene Tagging.....	37

A.	Overview of Gene Tagging.....	32
B.	Additional Materials Required.....	43
C.	General Considerations.....	44
D.	Cloning of Homology Arms into Vector.....	44
E.	Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells.....	49
V.	PrecisionX™ Targeting Vectors By Application: Gene Editing.....	51
A.	Overview of Gene Editing.....	53
B.	Additional Materials Required (for PBHR100A-1 vector)..	53
C.	General Considerations (for PBHR100A-1 vector).....	54
D.	Design and Cloning of Homology Arms into Vector (for PBHR100A-1 vector).....	55
E.	Protocol for Gene Editing with TALEN/Cas9 system and piggyBac HR vector.....	61
VI.	References.....	63
VII.	Technical Support.....	65
VIII.	Licensing and Warranty Statement.....	65

## I. Introduction

### A. Overview of Genome Engineering Tools

Recent advances in tools available for precise genome engineering of target cells have revolutionized the field of biology. The first of these revolutionary tools, zinc finger nuclease (ZFN) technology, has demonstrated that it can be successfully employed to introduce a double-stranded DNA break at a defined location in the mammalian genome to

introduce transgenes by homologous recombination (HR) (Hockemeyer *et al.* 2009). This finding has since spurred follow-up studies using HR-based approaches to disrupt endogenous gene function and track cells that have been modified by ZFNs.

However, inherent limitations of the ZFN technology have driven researchers to look at other alternatives that would be able to 1) target any sequence of interest with efficiencies rivaling ZFNs 2) increased predictability of targeting and 3) cost-effective to implement for routine use in laboratories around the world. Recent work with transcriptional activator-like effector nucleases (TALENs), first described in the plant pathogen *Xanthomonas sp.*, have shown that researchers can efficiently target any genomic DNA sequence using a pair of custom TALEN proteins whose DNA-binding modules can recognize individual DNA nucleotides based on an elegant amino acid cipher (Boch *et al.* 2009). Based on the explosion of recent publications that have successfully demonstrated the utility of the TALEN platform in many eukaryotic model organisms ranging from zebrafish to humans, it has quickly established itself as the *de facto* platform for genome engineering.

While TALEN technology has been gaining acceptance in the research community, the continual advances being made in this field has uncovered another novel, yet potentially powerful tool for genome engineering based on a bacterial system for combating invading viral and nucleic acid sequences, termed CRISPR/Cas9. Originally discovered in the pathogenic bacterium *Streptococcus pyogenes*, this system uses an unusual series of short RNAs (collectively termed as “guide RNA”) to specifically target a complementary DNA sequence, and upon binding, leads to recruitment of an endonuclease called Cas9 to specifically induce a double-stranded break (DSB) in the DNA sequence targeted by the guide RNA

(Bhaya *et al.* 2011, Jinek *et al.* 2012). Similar to TALENs, the introduction of DSBs in the DNA leads to recruitment of the cellular repair machinery to drive non-homologous end-joining (NHEJ) processes in absence of homologous sequences, whose error-prone nature often leads to formation of indels in the targeted DNA sequences. The main advantage of the CRISPR/Cas9 system over the TALEN system is that it only requires an oligonucleotide sequence complementary to the target sequence in the format  $N_{20}NGG$  (where N=any nucleotide) to effect cleavage, which is much easier and less time-consuming than other technologies. With a simple design of an oligo sequence to target any gene or non-coding sequence, this opens up possibilities for precise, multiplexed genome engineering that would have been unthinkable even a few years ago.

However, in order to fully leverage the powerful nature of these genome engineering platforms, there needs to be a set of ancillary tools that will allow detection of cells that have been modified by ZFN, TALEN or Cas9 systems. The majority of genome editing events affect only a percentage of the total population of transfected cells, estimated at 1-80% for mono or bi-allelic modifications depending on the platform, the cell type, and the DNA target of interest. Due to the wide variation in activity, performing downstream phenotypic assays in a background of wild-type cells is challenging, especially if the phenotype in question is subtle or difficult to discern. Therefore, tools in the form of donor or targeting vectors that contain 1) fluorescent or antibiotic selection markers and 2) a gene fragment of interest to knock-in, knock-out, or correct a wild-type sequence will be extremely useful for obtaining a homogeneous population of cells whose genomes have been modified by ZFN, TALEN or Cas9 systems. These vectors, when co-transfected into target cells along with ZFNs, TALEN or Cas9 vectors will take advantage of intrinsic homologous recombination (HR) processes conserved across all organisms for accurate repair of DSBs and insertion of

exogenous gene or marker of interest. This is made possible by the fact that HR efficiencies are increased >1000-fold in the presence of a homologous donor sequence in mammalian cells (Choulika *et al.* 1995), making this a tractable approach to characterize phenotypic changes associated with gene targeting and potentially useful as a platform for gene therapy applications.

## B. Overview of PrecisionX™ Targeting Vectors for HR-Based Applications

We have generated a suite of homologous targeting vectors suitable for the purpose of engineering both protein coding and non-coding genes including microRNAs and LncRNAs (see Table 1):

**Table 1. List of Available PrecisionX™ HR Targeting Vectors By Application**

### 1) Gene Knock-Out

Cat#	Description	Size
HR110PA-1	HR Targeting Vector w/Dual Selection Markers (RFP+Puro)	10 µg
HR210PA-1	HR Targeting Vector w/Dual Selection Markers (GFP+Puro) and TK (for negative selection)	10 µg
HR410PA-1	HR Targeting Vector w/Dual Selection Markers (GFP+Puro)	10 µg
HR510PA-1	HR Targeting Vector w/Dual Selection Markers (RFP+Hygromycin)	10 µg

<b>HR700PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (GFP+Puro) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>
<b>HR710PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (RFP+Hygro) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>
<b>HR720PA-1</b>	<b>HR Targeting Vector w/Single Selection Marker (Blasticidin) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>

## 2) Gene Knock-In

<b>Cat#</b>	<b>Description</b>	<b>Size</b>
<b>HR100PA-1</b>	<b>Basic HR Targeting Vector For Gene Knock-In or Knock-Out (MCSs only)</b>	<b>10 µg</b>

## 3) Reporter Gene Tagging & Co-Expression

<b>Cat#</b>	<b>Description</b>	<b>Size</b>
<b>HR120PA-1</b>	<b>GFP-Fusion HR Targeting Vector (w/EF1a-RFP-T2A-Puro cassette)</b>	<b>10 µg</b>

<b>HR130PA-1</b>	<b>T2A-GFP Co-Expression HR Targeting Vector</b>	<b>10 µg</b>
<b>HR150PA-1</b>	<b>GFP-T2A-Luciferase HR Targeting Vector (GFP-Fusion with Co-expressed Luciferase)</b>	<b>10 µg</b>
<b>HR180PA-1</b>	<b>IRES-GFP Co-Expression HR Targeting Vector</b>	<b>10 µg</b>
<b>HR220PA-1</b>	<b>GFP-Fusion HR Targeting Vector (w/EF1a-RFP-T2A-Hygro cassette)</b>	<b>10 µg</b>

#### 4) Gene Editing

<b>Cat#</b>	<b>Description</b>	<b>Size</b>
<b>PBHR100A-1</b>	<b>piggyBac-based HR Donor Vector with GFP+Puro markers and TK (for negative selection)</b>	<b>10 µg</b>
<b>HR110PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (RFP+Puro)</b>	<b>10 µg</b>
<b>HR210PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (GFP+Puro) and TK (for negative selection)</b>	<b>10 µg</b>
<b>HR410PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (GFP+Puro)</b>	<b>10 µg</b>
<b>HR510PA-1</b>	<b>HR Targeting Vector w/Dual Selection Markers (RFP+Hygromycin)</b>	<b>10 µg</b>

Cat#	Description	Size
HR700PA-1	<b>HR Targeting Vector w/Dual Selection Markers (GFP+Puro) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>
HR710PA-1	<b>HR Targeting Vector w/Dual Selection Markers (RFP+Hygro) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>
HR720PA-1	<b>HR Targeting Vector w/Single Selection Marker (Blasticidin) and Negative Selection (TK) Against Random Integration</b>	<b>10 µg</b>

### Key Features:

- 1) Dual markers: Homology-directed recombination using an exogenously supplied vector can result in targeted gene addition. To facilitate the enrichment and selection of HR modified cells, fluorescent markers (RFP and GFP) and positive selection markers (Puromycin or Hygromycin depending on vector format) are built into the HR vectors. The EF-1 $\alpha$  alpha promoter is chosen to drive the dual marker expression because EF-1 $\alpha$  promoter tends to be active in most cell types, including primary and stem cells.
- 2) Insulation: To ensure the expression of knock-in genes, insulation sequences are placed on both sides of the expression cassette. This feature ensures the knock-in cassette minimally impacts neighboring genes, and is not silenced by the local environment where the knock-in cassette has been placed.

- 3) Built-in LoxP Sites. The built-in LoxP sites allow for conditional knock-in/knock-out. In addition, by removal of knock-in genes via CRE, the HR targeting vector can be reused to knockout another allele if double knockouts have not been achieved during the first attempt.
- 4) Extra MCS. Extra MCS are available for researchers who wish to implement more complicated applications such as gene correction, expression monitoring, and gene therapy.
- 5) **For HR7xx series vectors:** These specialized HR donor vectors contain a negative selection marker (*hsvTK*) that is located OUTSIDE of the 3' MCS to select against randomly integrated donor vector (please see Section II.A for additional details).


**Note: Above features do not apply to PBHR100A-1 vector except for insulators. PBHR100A-1 vector requires the use of the Excision-Only SuperPiggyBac transposase enzyme (Cat. no. PB220PA-1) for removal of selection cassette present in the vector.**

## C. Design Considerations for Homologous Recombination Arms

With several different genome engineering tools available to induce targeted DSBs, such as SBI's EZ-TAL™ TALEN kit (<http://www.systembio.com/tale-effector-technology/nuclease>) and Cas9 SmartNuclease system (<http://www.systembio.com/genome-engineering-cas9-crispr-smartnuclease>), there is a need for suitable vectors that will serve as the guide template to induce HR in cells that have been targeted using the aforementioned tools. A general targeting vector will contain the following features:

- 1) Homologous sequences (to the template with DSBs) at 5' and 3' ends of the DSB site.
- 2) An expression cassette bearing a promoter, insert (cDNA, microRNA, non-coding RNA, etc.), fluorescent marker, or selection agent to select cells that have undergone HR.

The 5' and 3' homologous sequences (termed “homology arms”) should be exactly homologous to the template with the DSB, and preferably adjacent to the actual DSB site. For TALEN-mediated cleavage, this would entail designing homology arms that are adjacent to the spacer region between TALEN binding sites, which spans 15-30bp and is the site of the DSB. An example of homology arm design is shown in Figure 1.




**Figure 1. Schematic Diagram of TALEN-mediated Homologous Recombination (HR) Process**

In this figure, a region of the human *AAVS1* locus (in blue) that is targeted by a pair of TALENs is used for inserting an EF1 $\alpha$ -RFP-T2A-Puro expression cassette (in red) present in the HR targeting vector. The HR targeting vector also bears homology arms at the 5' and 3' end of the expression cassette that contains ~0.8kb of sequence

homology to the genomic DNA surrounding the *AAVS1* locus. This region of homology is crucial for the success of the homologous recombination reaction, as it serves as the guide template for specifically targeting the exogenous cassette into this genomic locus. The typical range of sizes for homology arms varies by the application, but it should be anywhere from 0.5kb to 1kb for each arm for efficient recombination to occur. Please note that the actual regions of recombination (Red “Xs”) at the 5’ and 3’ of the target site can vary widely, thus it is difficult to predict the actual sites as this is determined by the cell.

For Cas9 system, the cleavage site is 2-3bp upstream of the protospacer adaptor motif (PAM) immediately following the guide RNA sequence, thus homology arms should be designed so that it is close as possible to the cut site (<10bp) at both ends for efficient HR reaction.



**Figure 2. Schematic Diagram of Cas9-mediated Homologous Recombination (HR) Process**

## II. PrecisionX™ Targeting Vectors By Application: Gene Knock-Out


### A. Overview of Gene Knock-Out

For HR-based applications needing the disruption of a particular protein-coding or non-coding gene of interest, a gene-knockout targeting vector would be desirable. This vector would include homology arms (0.5kb-1kb in length) at the 5' and 3' regions immediately flanking the DNA sequence of interest (e.g. exonic sequence) to specifically direct the targeting vector to insert an expression cassette where the targeted sequence lies, leading to the 'knock-out' of the sequence by gene addition. SBI's gene knockout targeting vectors are designed to be vectors of choice for this application. Gene-specific homology arms can be cloned into MCS1 and MCS2 to direct ZFN, TALEN or Cas9-mediated gene disruption. Correctly disrupted cell lines are selected as RFP-positive (GFP for HR210PA-1 and HR410PA-1) and puromycin-resistant colonies (hygromycin for HR510PA-1). An example of targeting by HR110PA-1 targeting vector is shown in Figure 10. In addition, the selection cassette is fully removable by expression of Cre recombinase (not included).


In addition, SBI also offers a series of three HR targeting vectors (Cat #HR700PA-1, HR710PA-1, and HR720PA-1, please see Figs. 7, 8, and 9 respectively for vector maps) containing a PGK-*hsvTK* selection cassette that is located outside of the 3' MCS homology arm cloning site. The major advantage of these vectors is to help reduce random integration events associated with plasmid transfection by using negative selection (*hsvTK*) to eliminate cells that have integrated the donor plasmid not

via homologous recombination through nuclease-mediated dsDNA breaks. When the constructs are randomly integrated into the genome, the PGK-*hsvTK* cassette will be active and these cells will be sensitized to ganciclovir or fialuridine when added to the culture medium. If the constructs are correctly targeted via HR, the PGK-*hsvTK* cassette is not incorporated into the genome because it is outside of the homology arms targeting the desired genomic locus, thus making correctly targeted cells insensitive to TK-mediated negative selection. While the frequency of random integration can be mitigated by careful titration of the donor vector, the use of these specialized vectors will assist the researcher in obtaining a higher frequency of correct clones (positive for HR without background integration) without the need for extensive optimization.


Finally, SBI also offers a customizable version of HR110PA-1 without the RFP-T2A-Puro selection cassette (Cat #HR100PA-1, see Section III). This vector contains a MCS which allows the user to insert their own cassette of choice to knock-out a particular sequence of interest.



**Figure 3. Schematic Diagram of HR110PA-1 Gene Knock-out Targeting Vector [MCS1-EF1 $\alpha$ -RFP-T2A-Puro-pA-MCS2]**



**Figure 4. Schematic Diagram of HR210PA-1 Gene Knock-out Targeting Vector [MCS1-EF1 $\alpha$ -GFP-T2A-Puro-P2A-TK-pA-MCS2]**



**Figure 5. Schematic Diagram of HR410PA-1 Gene Knock-out Targeting Vector [MCS1-EF1 $\alpha$ -GFP-T2A-Puro-pA-MCS2]**




**MCS** Multiple Cloning Sites

▶ LoxP Sites

■ Insulator Sequences

**Figure 6. Schematic Diagram of HR510PA-1 Gene Knock-out Targeting Vector [MCS1-EF1 $\alpha$ -RFP-T2A-Hygro-pA-MCS2]**




**MCS** Multiple Cloning Sites

▶ LoxP Sites

■ Insulator Sequences

**Figure 7. Schematic Diagram of HR700PA-1 Gene Knock-out Targeting Vector with Negative Selection Against Random Integration [MCS1-EF1 $\alpha$ -GFP-T2A-Puro-pA-MCS2-PGK-hsvTK]**




**MCS** Multiple Cloning Sites

▶ LoxP Sites

■ Insulator Sequences

**Figure 8. Schematic Diagram of HR710PA-1 Gene Knock-out Targeting Vector with Negative Selection Against Random Integration [MCS1-EF1 $\alpha$ -RFP-T2A-Hygro-pA-MCS2-PGK-*hsvTK*]**




**MCS** Multiple Cloning Sites

▶ LoxP Sites

■ Insulator Sequences

**Figure 9. Schematic Diagram of HR720PA-1 Gene Knock-out Targeting Vector with Negative Selection Against Random Integration [MCS1-EF1 $\alpha$ -Blasticidin-pA-MCS2-PGK-*hsvTK*]**




Note: simple, robust, and optional removal of positive selection markers

**Figure 10. Example of HR110PA-1 Gene Knock-out Targeting Vector Utilized for Gene Targeting**

## B. Validation Data for HR110PA-1

### 1) Functionality of Dual Selection Markers


The dual positive selection markers, RFP and Puromycin resistance gene have been functionally validated in human HEK293T cells. Cells become RFP and puromycin-resistant 48 hours after initial transfection (Fig. 11). EF1 $\alpha$  promoter is a robust promoter that is active most cell types, including primary and stem cells.



**Figure 11. RFP and Puromycin-Resistance in HEK293T Cells Transfected with HR110PA-1 Targeting Vector**

## **2) Robustness of Insulators in Preventing Silencing of Expression Cassette**


To ensure that integration of expression cassette will not be inactivated by epigenetic modifications or influenced by neighboring sequences, insulator sequences are implemented in HR vectors. Cells with successful HR events show prolonged RFP-positivity and persistent puromycin-resistant for more than 3 months in culture (Fig. 12).



**Figure 12. HR-modified HEK293T Cells using Targeting Vector Remain RFP-Positive and Puromycin-Resistant for More Than 3 Months**

### **3) Excision of Expression Cassette with *Cre/LoxP* System**


The HR knock-out targeting vectors contain *LoxP* sites flanking the expression cassette (“floxed”) which can be used to easily remove the expression cassette in the presence of Cre recombinase (Fig. 13) by the well-validated *LoxP*-Cre excision reaction. This feature of the donor vector allows for conditional expression of the cassette and phenotype screening in the cell line of interest.



**Figure 13. Transfection of a CRE-GFP Construct Results in Efficient Abrogation of Floxed RFP-Puro Cassette Expression in HEK293T Cells**

#### 4) Validation of HR By Genotyping

We have designed a series of PCR primer pairs for the 5' and 3' homology arms to confirm the HR event has occurred in the targeted locus of interest. By using one of the primers in the primer pair to bind to the genomic site (either 5' or 3' ends) flanking the cassette, this allows for specific detection of HR events based on the expected sizes of the PCR amplicons (Fig. 14).



**Figure 14. PCR Genotyping Analysis Demonstrating Successful HR Reaction in HEK293T cells**

### C. Additional Materials Required

- a) LB Agar and Broth containing 50 µg/ml ampicillin or Kanamycin
- b) Any high-transformation efficiency RecA- and EndA- *E.coli* competent cells
- c) Dulbecco's Modified Eagle's Medium (D-MEM) high glucose with sodium pyruvate and glutamine (Invitrogen, Cat. # 11995073)
- d) Lipofectamine 2000 transfection reagent (Invitrogen, Cat. # 11668019)
- e) Qiagen EndoFree Plasmid Maxi Kit (Qiagen, Cat. # 12362)
- f) Zyppy™ Plasmid MiniPrep Kit (Zymo Research, Cat. # D4019)
- g) Phusion High-Fidelity DNA Polymerase (NEB, Cat. # M0530S)
- h) Qiagen MinElute Gel Extraction Kit (Qiagen, Cat. #28604)

- i) Fetal Bovine Serum (Invitrogen, Cat. # 16000036)
- j) Penicillin/Streptomycin (Invitrogen, Cat. # 15070063)
- k) Surveyor Mutation Detection Kit (Transgenomics, Cat. #706020)
- l) Trypsin-EDTA (Sigma, Cat. # T3924)
- m) 6-well Tissue Culture Plates and Related Tissue Culture Supplies
- n) Other specific media and additives specific for cell type of interest
- o) **\*\*Optional\*\*** - For difficult-to-transfect cells, the use of an electroporation system (e.g. Lonza's NucleoFector or Invitrogen's Neon system) is highly recommended

## D. General Considerations

We recommend propagation of the HR targeting vector plasmid prior to starting the experiments. The plasmids can be transformed using standard conditions suitable in any high-quality RecA- and EndA- *E.coli* competent cell.

For cells transformed with any of the targeting vectors, we suggest plating 50-200  $\mu$ l of transformed cells on fresh ampicillin or kanamycin plates (50  $\mu$ g/ml) depending on the HR vector. Incubate the plates at 37°C overnight. Colonies picked from the transformation can be grown at 37°C overnight in ~200 ml of LB media containing 50  $\mu$ g/ml of the appropriate antibiotic. After overnight growth, plasmid DNA can be harvested from culture using an endotoxin-free DNA plasmid maxiprep kit.

For confirmation of the plasmid, we recommend performing restriction digestion analysis or direct sequencing to confirm integrity of the amplified plasmids.

## E. Cloning of Homology Arms into Vector

For fast and efficient cloning of homology arms into the vector, we suggest the use of restriction enzyme, ligase-free methods such as SBI's Cold Fusion Kit. Cloning both of the arms into the targeting vector will require sequential cloning reactions, one for each arm. Please refer to the user manual for the Cold Fusion Kit for details on the cloning process (<http://www.systembio.com/molecular-tools/cold-fusion-cloning/literature>)

For standard cloning strategies, please refer to the following protocol for more details. Since there are two homology arms that need to be cloned, the cloning will need to be repeated for each arm.

### 1. PCR of homology arms

- a) PCR amplify homology arms using proofreading PCR enzymes (e.g. Phusion) with primers containing compatible restriction sites to MCS1 (for 5' homology arm) and MCS2 (for 3' homology arm).
- b) Run out PCR products on 1.5-2% agarose gel and extract correct size bands using a suitable gel extraction kit

### 2. Ligation of insert into vector

- a) Digest the targeting vector using appropriate restriction enzymes for the homology arm being cloned and gel-purify the digestion product
- b) Dilute gel-purified, digested vector to 10 ng/μl
- c) Set up 10 μl ligation reactions for each control and test samples as below:
- d)

<u>Volume</u>	<u>Item</u>
1.0 $\mu$ l	HR Targeting vector
7.0 $\mu$ l	5' or 3' homology arm PCR product (~70 ng) or water control
1.0 $\mu$ l	10X T4 DNA ligase buffer
1.0 $\mu$ l	T4 DNA Ligase (40 U/ $\mu$ l)
<b>10.0 <math>\mu</math>l</b>	<b>Total Reaction Volume</b>

- e) Incubate reactions at 25°C for 1-2 hours (sticky-end ligation) or O/N at 16°C (for blunt-end ligation)

### 3. Transform *E. coli* with the ligation product

Transform competent cells (with a transformation efficiency of at least  $1 \times 10^9$  colonies/ $\mu$ g pUC19) with the whole ligation reaction (10  $\mu$ l) following the protocol provided with the competent cells. Plate the transformed bacteria on 50  $\mu$ g/ml of appropriate antibiotic in agar plates.

### 4. Identify clones with the correct insert

- a) Depending on the ratio of colony numbers for the test sample vs. the negative control sample, randomly pick 2 or 3 colonies and grow each clone in 100  $\mu$ l of LB Broth with 100  $\mu$ g/ml of appropriate antibiotic at 37°C for 2 hours with shaking.
- b) Use 1  $\mu$ l of each bacterial culture for screening homology arm inserts by PCR and continue to grow the culture for another 4 hours. Store the culture at 4°C.
- c) Prepare a PCR Master Mix with PCR primers flanking the insert:

1 rxn		10 rxn		Composition
0.5	μl	5	μl	PCR primer 1 (10 μM)
0.5	μl	5	μl	PCR primer 2 (10 μM)
0.5	μl	5	μl	50X dNTP mix (10 mM of each)
2.5	μl	25	μl	10X PCR Reaction Buffer
19.5	μl	195	μl	Nuclease-free water
0.5	μl	5	μl	Taq DNA polymerase (approx. 5 U/μl)
<b>24.0</b>	<b>μl</b>	<b>240</b>	<b>μl</b>	<b>Total volume</b>

d) Mix the master mix very well and aliquot 24 μl into each well of 96-well PCR plate or individual tubes.

e) Add 1 μl of each bacterial culture from step (b) into each well (or tube).

f) Proceed with PCR using the following program:

94°C, 4 min	1 cycle
94°C, 0.5 min, then 68°C, 1 min/1 kb*	25 cycles
68°C, 3 min	1 cycle

\* Depending on the size of final PCR product, use a shorter or longer time.

g) Take 5 μl of the PCR reaction and run it on a 1.5-2% agarose/EtBr gel in 1X TAE buffer to identify clones with correct insert.

h) Grow a positive clone containing insert in an appropriate amount of LB broth containing antibiotics, and purify the construct using an endotoxin-free plasmid purification kit. Sequence verification of the insert is optional.

i) Repeat the entire process with remaining homology arm to build the final targeting construct

## **F. Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells**

### Notes:

- I. Depending on the cell type being transfected, please choose a transfection protocol which results in maximal transfection efficiencies. For adherent cell lines such as HEK293T, passive transfection methods using cationic lipid-based methods (e.g. Lipofectamine 2000, FuGene HD) work very well in transfection of targeting and TALEN/Cas9 plasmids. For other types of cells such as primary, stem, or suspension cells, we suggest transfection using electroporation methods (NucleoFection or Neon) for optimal results.
  - II. The plasmids should be mixed well in minimal serum/no antibiotic media + cationic lipid transfection reagent or electroporation buffer to maximize efficiency of delivery.
  - III. For selection of target cells, we strongly recommend testing different concentrations of puromycin on untransfected cells to determine the optimal concentration of selection agent – which is kills ~90-100% of cells within 48-72 hours after administration.
- 1) Plate 200,000 to 300,000 cells (e.g. 293T cells) into a single well of a 6-well plate in 2 ml of appropriate growth medium. Include a single well of cells as

negative control (which can be non-relevant plasmid DNA)

- 2) Next day, or when cells are 50-60% confluent, co-transfect target cells with either TALEN or Cas9 plasmid(s) and the HR targeting vector using a suitable transfection reagent following the manufacturer's recommended protocol for 6-well plates. The use of reduced or serum-free media containing no antibiotics to dilute the vector/transfection complex is highly recommended.

**Note: For 293T cells, we suggest 0.5 µg of each TALEN (1 µg total for pair) or 0.5 µg of SBI's Cas9 SmartNuclease vector in conjunction with 1 µg of the HR targeting vector into cells for efficient cleavage and HR reaction. For other cell types, we suggest optimizing the amounts and ratios of Cas9 or TALEN to targeting vector for optimal results in a target cell line.**

- 3) Allow at least 12 hours before changing transfection media to complete growth media
- 4) Assay for positive HR events 48-96 hours after co-transfection. Select cells with insertion of targeting vector using FACS-based sorting of RFP/GFP marker or Puromycin/Hygromycin/Blasticidin selection. If using antibiotic selection, select cells for a minimum of 5 days prior to further characterization. Cutting efficiency of TALEN/Cas9 can be measured by Surveyor Nuclease Assay and HR efficiency by FACS based on % of RFP or GFP positive cells
- 5) Genotyping of selected cells can be done by PCR of genomic DNA-insert junctions at 5' and 3' ends of HR


site (see Fig. 14 for details). In brief, design PCR primer pairs with one of the primers in the genomic DNA region and the other located in the HR vector. Both 5' and 3' integration ends can be confirmed for complete verification of insertion.

- 6) Selected cells can be clonally expanded for further phenotypic characterization

### **III. PrecisionX™ Targeting Vectors By Application: Gene Knock-In**

#### **A. Overview of Gene Knock-In**

For targeted insertion of any expression cassette into a specific location in the genome such as a “safe harbor” locus (Sadelain *et al.* 2011) for stable expression with minimal context-dependent effects, a HR-based vector encoding for homology arms at 5' and 3' ends with an expression cassette between the homology arms would be desirable. This strategy will require co-transfection of ZFN, TALEN or Cas9 plasmids with the knock-in targeting vector into a safe harbor locus (e.g. human *AAVS1* or mouse *Rosa26*). SBI's gene knock-in targeting vector (Cat #HR100PA-1) is designed to be the vector of choice for this application (Fig. 15).



**Figure 15. Schematic Diagram of HR100PA-1 Gene Knock-in Targeting Vector**

This is the most basic HR targeting vector which allows the most flexibility for the end-user. Homology arm sequences of 0.5kb to 1kb can be cloned into MCS1 (5' HA) and MCS4 (3' HA), and an expression cassette of your choice can then be cloned into MCS3. The maximal size of the expression cassette can be upwards of 6-8kb – larger sizes may impact transfection and homologous recombination efficiencies. The targeting plasmid can be co-transfected along with ZFN, TALEN or Cas9 systems to have simultaneous cutting and HR-based ‘knock-in’ at the desired genomic locus. After knock-in, the expression cassette can be excised by expression of Cre recombinase, leaving a single *LoxP* site in between the 5' and 3' homology arms.

Please note that the HR100PA-1 vector can also be utilized as a knock-out vector (see Section II.A) as the user can clone in their cassette of choice into MCS3 to target a particular gene sequence.

## B. Additional Materials Required

a) LB Agar and Broth containing 50 µg/ml ampicillin

- b) Any high-transformation efficiency RecA- and EndA- *E.coli* competent cells
- c) Dulbecco's Modified Eagle's Medium (D-MEM) high glucose with sodium pyruvate and glutamine (Invitrogen, Cat. # 11995073)
- d) Lipofectamine 2000 transfection reagent (Invitrogen, Cat. # 11668019)
- e) Qiagen EndoFree Plasmid Maxi Kit (Qiagen, Cat. # 12362)
- f) Zippy™ Plasmid MiniPrep Kit (Zymo Research, Cat. # D4019)
- g) Phusion High-Fidelity DNA Polymerase (NEB, Cat. # M0530S)
- h) Qiagen MinElute Gel Extraction Kit (Qiagen, Cat. #28604)
- i) Fetal Bovine Serum (Invitrogen, Cat. # 16000036)
- j) Penicillin/Streptomycin (Invitrogen, Cat. # 15070063)
- k) Surveyor Mutation Detection Kit (Transgenomics, Cat. #706020)
- l) Trypsin-EDTA (Sigma, Cat. # T3924)
- m) 6-well Tissue Culture Plates and Related Tissue Culture Supplies
- n) Other specific media and additives specific for cell type of interest
- o) **\*\*Optional\*\*** - For difficult-to-transfect cells, the use of an electroporation system (e.g. Lonza's NucleoFector or Invitrogen's Neon system) is highly recommended

## C. General Considerations

We recommend propagation of the HR targeting vector plasmid prior to starting the experiments. The plasmids can be transformed using standard conditions suitable in any high-quality RecA- and EndA- *E.coli* competent cell.

For cells transformed with any of the targeting vectors, we suggest plating 50-200 µl of transformed cells on fresh ampicillin plates (50 µg/ml). Incubate the plates at 37°C overnight. Colonies picked from the transformation can be grown at 37°C overnight in ~200 ml of LB media containing 50 µg/ml of ampicillin. After overnight growth, plasmid DNA can be harvested from culture using an endotoxin-free DNA plasmid maxiprep kit.

For confirmation of the plasmid, we recommend performing restriction digestion analysis or direct sequencing to confirm integrity of the amplified plasmids.

## **D. Cloning of Homology Arms and Expression Cassette into Vector**

For fast and efficient cloning of homology arms and expression cassette into the vector, we suggest the use of restriction enzyme, ligase-free methods such as SBI's Cold Fusion Kit. Cloning both of the arms and the expression cassette into the targeting vector will require a total of three cloning reactions. Please refer to the user manual for the Cold Fusion Kit for details on the cloning process (<http://www.systembio.com/molecular-tools/cold-fusion-cloning/literature>)

For standard cloning strategies, please refer to the following protocol for more details. Since there are two homology arms and an expression cassette that need to be cloned, the cloning will need be done in a sequential manner except for the initial PCR reaction.

### **1. PCR of homology arms and expression cassette**

- a) PCR amplify homology arms and expression cassette using proofreading PCR enzymes (e.g. Phusion) with primers containing compatible restriction sites to MCS1 (for 5' homology arm), MCS3 for expression cassette, and MCS4 (for 3' homology arm). Depending on the size of the expression cassette, it may be preassembled by prior cloning through fusion assembly methods (e.g. SBI's Cold Fusion) or long-range PCR
- b) Run out PCR products on 1-2% agarose gel and extract correct size bands using a suitable gel extraction kit

## 2. Ligation of insert into vector

- a) Digest the targeting vector using appropriate restriction enzymes for the homology arm(s) or expression cassette being cloned and gel-purify the digestion product
- b) Dilute gel-purified, digested vector to 10 ng/ $\mu$ l
- c) Set up 10  $\mu$ l ligation reactions for each control and test samples as below:

<u>Volume</u>	<u>Item</u>
1.0 $\mu$ l	HR Targeting vector
7.0 $\mu$ l	5' or 3' homology arm or expression cassette PCR product (~70 ng) or water control
1.0 $\mu$ l	10X T4 DNA ligase buffer
1.0 $\mu$ l	T4 DNA Ligase (40 U/ $\mu$ l)
<b>10.0 <math>\mu</math>l</b>	<b>Total Reaction Volume</b>

- d) Incubate reactions at 25°C for 1-2 hours (sticky-end ligation) or O/N at 16°C (for blunt-end ligation)

## 3. Transform *E. coli* with the ligation product

Transform competent cells (with a transformation efficiency of at least  $1 \times 10^9$  colonies/ $\mu\text{g}$  pUC19) with the whole ligation reaction (10  $\mu\text{l}$ ) following the protocol provided with the competent cells. Plate the transformed bacteria on 50  $\mu\text{g}/\text{ml}$  ampicillin agar plates

#### 4. Identify clones with the correct insert

- a) Depending on the ratio of colony numbers for the test sample vs. the negative control sample, randomly pick 2 or 3 colonies and grow each clone in 100  $\mu\text{l}$  of LB Broth with 100  $\mu\text{g}/\text{ml}$  ampicillin at 37°C for 2 hours with shaking.
- b) Use 1  $\mu\text{l}$  of each bacterial culture for screening homology arm inserts by PCR and continue to grow the culture for another 4 hours. Store the culture at 4°C.
- c) Prepare a PCR Master Mix with PCR primers flanking the insert:

1 rxn		10 rxn		Composition
0.5	$\mu\text{l}$	5	$\mu\text{l}$	PCR primer 1 (10 $\mu\text{M}$ )
0.5	$\mu\text{l}$	5	$\mu\text{l}$	PCR primer 2 (10 $\mu\text{M}$ )
0.5	$\mu\text{l}$	5	$\mu\text{l}$	50X dNTP mix (10 mM of each)
2.5	$\mu\text{l}$	25	$\mu\text{l}$	10X PCR Reaction Buffer
19.5	$\mu\text{l}$	195	$\mu\text{l}$	Nuclease-free water
0.5	$\mu\text{l}$	5	$\mu\text{l}$	Taq DNA polymerase (approx. 5 U/ $\mu\text{l}$ )
<b>24.0</b>	<b><math>\mu\text{l}</math></b>	<b>240</b>	<b><math>\mu\text{l}</math></b>	<b>Total volume</b>

- d) Mix the master mix very well and aliquot 24  $\mu\text{l}$  into each well of 96-well PCR plate or individual tubes.
- e) Add 1  $\mu\text{l}$  of each bacterial culture from step (b) into each well (or tube).
- f) Proceed with PCR using the following program:

94°C, 4 min	1 cycle
94°C, 0.5 min, then 68°C, 1 min/1 kb*	25 cycles
68°C, 3 min	1 cycle

\* Depending on the size of final PCR product, use a shorter or longer time.

g) Take 5  $\mu$ l of the PCR reaction and run it on a 1.5-2% agarose/EtBr gel in 1X TAE buffer to identify clones with correct insert.

h) Grow a positive clone containing insert in an appropriate amount of LB-ampicillin broth, and purify the construct using an endotoxin-free plasmid purification kit. Sequence verification of the insert is optional.

i) Repeat the entire process with remaining homology arm or expression cassette to build the final targeting construct

## **E. Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells**

### Notes:

- I. Depending on the cell type being transfected, please choose a transfection protocol which results in maximal transfection efficiencies. For adherent cell lines such as HEK293T, passive transfection methods using cationic lipid-based methods (e.g. Lipofectamine 2000, FuGene HD) work very well in transfection of targeting and TALEN/Cas9 plasmids. For other types of cells such as primary, stem, or suspension cells, we suggest transfection using electroporation methods (NucleoFection or Neon) for optimal results.
  - II. The plasmids should be mixed well in minimal serum/no antibiotic media + cationic lipid transfection reagent or electroporation buffer to maximize efficiency of delivery.
  - III. For selection of target cells, we strongly recommend testing different concentrations of puromycin on untransfected cells to determine the optimal concentration of selection agent – which is kills ~90-100% of cells within 48-72 hours after drug administration.
- 1) Plate 200,000 to 300,000 cells (e.g. 293T cells) into a single well of a 6-well plate in 2 ml of appropriate growth medium. Include a single well of cells as negative control (which can be non-relevant plasmid

DNA)

- 2) Next day, or when cells are 50-60% confluent, co-transfect target cells with either TALEN or Cas9 plasmid(s) and the HR targeting vector using a suitable transfection reagent following the manufacturer's recommended protocol for 6-well plates. The use of reduced or serum-free media containing no antibiotics to dilute the vector/transfection complex is highly recommended.

**Note: For 293T cells, we suggest 0.5 µg of each TALEN (1 µg total for pair) or 0.5 µg of SBI's Cas9 SmartNuclease vector in conjunction with 1 µg of the HR targeting vector into cells for efficient cleavage and HR reaction. For other cell types, we suggest optimizing the amounts and ratios of Cas9 or TALEN to targeting vector for optimal results in a target cell line.**

- 3) Allow at least 12 hours before changing transfection media to complete growth media
- 4) Assay for positive HR events 48-96 hours after co-transfection. Select cells with insertion of targeting vector using fluorescent or antibiotic selection (if applicable). If using antibiotic selection, select cells for a minimum of 5 days prior to further characterization. Cutting efficiency of TALEN/Cas9 can be measured by Surveyor Nuclease Assay and HR efficiency by % of fluorescence signal (if applicable) via FACS. If there is no available marker, the HR efficiency can be determined by performing dilution of transfected cells and individual characterization of clones by direct sequencing
- 5) Genotyping of selected cells can be done by PCR of

genomic DNA and homology arm sites MCS1 and MCS4 as well as the expression cassette (see Section II. Fig. 14 for details). In brief, design PCR primer pairs with one of the primers in the genomic DNA region and the other located in the HR vector. For the expression cassette, one of the primers can be internal to the expression cassette and the other primer immediately outside of the cassette.


- 6) Selected cells can be clonally expanded for further phenotypic characterization

## **IV. PrecisionX™ Targeting Vectors By Application: Gene Tagging**

### **A. Overview of Gene Tagging**


If the researcher is interested in tagging an endogenous sequence of interest with a marker to create a fusion or to co-express a tag with the endogenous sequence (e.g. tracking spatial dynamics of a protein, lineage tracking, etc.), SBI provides a collection of four tagging and co-expression vectors for HR-based applications (Figs. 16-21). Any type of DNA sequence, from protein-coding ORFs to non-coding sequences can be efficiently tagged or co-expressed with GFP and/or luciferase using these vectors. For protein-coding genes, the 3' end (typically the last exon) can be targeted using these vectors for efficient generation of C-terminal fusion proteins or co-expression of the gene and the tag via T2A or IRES elements if a fusion protein is not suitable.

## 1) GFP Fusion Vector




**Figure 16. Schematic Diagram of HR120PA-1 GFP Fusion Targeting Vector [GFP-pA-EF1 $\alpha$ -RFP-T2A-Puro-MCS]**

The HR120PA-1 targeting vector is designed to direct the in-frame fusion of GFP to the 3' end of your gene of choice. The 5' homology arm/GFP fusion segment is created using seamless fusion cloning technologies (e.g. SBI's Cold Fusion Kit) at the provided EcoRI site. The 3' homology arm can be cloned into the MCS using Cold Fusion or traditional cloning using the provided restriction sites. This vector offers a Puromycin resistance marker for selection of HR-positive cells. An example of utilizing the HR120PA-1 vector to generate a C-terminal GFP fusion construct is illustrated in Fig. 18. The EF1 $\alpha$ -RFP-T2A-Puro cassette can be excised by expression of Cre recombinase, leaving only the GFP fusion protein and a single *LoxP* site.




**Figure 17. Schematic Diagram of HR220PA-1 GFP Fusion Targeting Vector [GFP-pA-EF1 $\alpha$ -RFP-T2A-Puro-MCS]**

The HR220PA-1 targeting vector is a variant of HR120PA-1 designed to direct the in-frame fusion of GFP to the 3' end of your gene of choice. The 5' homology arm/GFP fusion segment is created using seamless fusion cloning technologies (e.g. SBI's Cold Fusion Kit) at the provided restriction sites. The 3' homology arm can be cloned into the MCS using Cold Fusion or traditional cloning using the provided restriction sites. This vector offers a Hygromycin-selection marker for selection of HR-positive cells. The EF1 $\alpha$ -RFP-T2A-Hygro cassette can be excised by expression of Cre recombinase, leaving only the GFP fusion protein and a single *LoxP* site.



**Figure 18. Utilization of HR120PA-1 GFP Fusion Targeting Vector to Specifically Tag C-terminal Region of Gene of Interest**

## 2) GFP Co-Expression Vector




**Note:** GFP is co-expressed with endogenous gene. However, GFP protein is not fused with endogenous protein. This format is used when fusion protein will interfere with the function of endogenous protein.

**Figure 19. Schematic Diagram of HR130PA-1 GFP Co-Expression Targeting Vector [T2A-GFP-pA-EF1 $\alpha$ -RFP-T2A-Puro-MCS]**

The HR130PA-1 vector is designed to direct the T2A-peptide/GFP linkage to the 3' end of your gene of choice. This vector will be used in cases where a GFP fusion will

affect the function of the endogenous gene of interest. The 5' homology arm/T2A-GFP linkage segment is created using seamless fusion cloning technologies (e.g. SBI's Cold Fusion Kit) at the provided EcoRI site. The 3' homology arm can be cloned into the MCS using Cold Fusion or by traditional cloning using the provided restriction sites. The EF1 $\alpha$ -RFP-T2A-Puro cassette can be excised by expression of Cre recombinase, leaving only the GOI-T2A-GFP construct and a single *LoxP* site.

### 3) GFP Fusion and Luciferase Co-Expression Vector




**Note:** Endogenous protein is tagged with GFP. Co-expression of luciferase allows quantification of gene expression and in vivo imaging via BLI. Therefore, gene expression, localization and expression levels can be monitored and quantified.

### Figure 20. Schematic Diagram of HR150PA-1 GFP Fusion and Luciferase Co-Expression Targeting Vector [GFP-T2A-Luc-pA-EF1 $\alpha$ -RFP-T2A-Puro-MCS]

The HR150PA-1 vector is designed to direct the in-frame fusion of GFP-T2A-Luciferase to the 3' end to your gene of your choice. The gene of interest will be fused to GFP and will express luciferase as a separate protein for applications where quantitative measurement of endogenous gene expression is necessary. The 5' homology arm/GFP-T2A-Luciferase fusion segment is created using seamless fusion cloning technologies (e.g.

SBI's Cold Fusion Kit) at the provided EcoRI site. The 3' homology arm can be cloned into the MCS using Cold Fusion or by traditional cloning using the provided restriction sites. The EF1 $\alpha$ -RFP-T2A-Puro cassette can be excised by expression of Cre recombinase, leaving only the GFP Fusion-T2A-luciferase construct and a single *LoxP* site.

#### 4) GFP Co-Expression Vector via Internal Ribosomal Entry Site (IRES)



**Note:** Endogenous transcripts will be transcribed along with IRES-GFP sequences. This format simplify the design and allow easy implementation. Useful for non-coding RNA studies when fusion protein is not an option.

#### Figure 21. Schematic Diagram of HR180PA-1 IRES-GFP Co-Expression Targeting Vector [IRES-GFP-pa-MCS1-EF1 $\alpha$ -RFP-T2A-Puro-MCS2]

The HR180PA-1 vector is designed for coupled expression of the gene of interest to an eGFP marker by the use of an IRES element, which aids in translation of the eGFP protein. This vector is especially useful for monitoring expression of non-coding transcripts (e.g. microRNA, lincRNAs) that are targeted by TALEN or Cas9 systems for functional studies. The 5' homology arm can be cloned into the vector by traditional restriction enzyme cloning or by fusion cloning methods to the 3' end of your gene of your choice using the provided restriction sites directly upstream of the IRES element. The 3' homology

arm can be cloned into MCS2 using the same cloning methods. Since the IRES controls the expression of the eGFP, the gene of interest does not have to be in-frame with eGFP and in cases of cDNA inserts, the presence of a stop codon in the sequence does not affect eGFP expression. The EF1 $\alpha$ -RFP-T2A-Puro cassette can be excised by expression of Cre recombinase, leaving only the GOI-IRES-eGFP construct and a single *LoxP* site.

## B. Additional Materials Required

- a) LB Agar and Broth containing 50  $\mu$ g/ml ampicillin
- b) Any high-transformation efficiency RecA- and EndA-*E.coli* competent cells
- c) Dulbecco's Modified Eagle's Medium (D-MEM) high glucose with sodium pyruvate and glutamine (Invitrogen, Cat. # 11995073)
- d) Lipofectamine 2000 transfection reagent (Invitrogen, Cat. # 11668019)
- e) Qiagen EndoFree Plasmid Maxi Kit (Qiagen, Cat. # 12362)
- f) Zyppy™ Plasmid MiniPrep Kit (Zymo Research, Cat. # D4019)
- g) Phusion High-Fidelity DNA Polymerase (NEB, Cat. # M0530S)
- h) Qiagen MinElute Gel Extraction Kit (Qiagen, Cat. #28604)
- i) Fetal Bovine Serum (Invitrogen, Cat. # 16000036)
- j) Penicillin/Streptomycin (Invitrogen, Cat. # 15070063)
- k) Surveyor Mutation Detection Kit (Transgenomics, Cat. #706020)
- l) Trypsin-EDTA (Sigma, Cat. # T3924)
- m) Promega Luciferase Assay System (Promega, Cat. # E4030)

- n) 6-well Tissue Culture Plates and Related Tissue Culture Supplies
- o) Other specific media and additives specific for cell type of interest
- p) **\*\*Optional\*\*** - For difficult-to-transfect cells, the use of an electroporation system (e.g. Lonza's NucleoFector or Invitrogen's Neon system) is highly recommended

## C. General Considerations

We recommend propagation of the HR targeting vector plasmid prior to starting the experiments. The plasmids can be transformed using standard conditions suitable in any high-quality RecA- and EndA- *E.coli* competent cell.

For cells transformed with any of the targeting vectors, we suggest plating 50-200  $\mu$ l of transformed cells on fresh ampicillin plates (50  $\mu$ g/ml). Incubate the plates at 37°C overnight. Colonies picked from the transformation can be grown at 37°C overnight in ~200 ml of LB media containing 50  $\mu$ g/ml of ampicillin. After overnight growth, plasmid DNA can be harvested from culture using an endotoxin-free DNA plasmid maxiprep kit.

For confirmation of the plasmid, we recommend performing restriction digestion analysis or direct sequencing to confirm integrity of the amplified plasmids.

## D. Cloning of Homology Arms into Vector

In general, for all of the tagging and co-expression vectors with a single EcoRI site (excluding HR180PA-1), the cloning of the 5' homology arm is best achieved through fusion cloning strategies such as SBI's Cold Fusion Kit. Please refer to the user manual for the Cold Fusion Kit for



Reverse Primer:

5' – CTCGTCGCTCTCCATGAATTCNNNNNNNNNNNNNNNNNN – 3'

N = 5' Homology Arm sequence (18-20bp in length)

Cloning of both homology arms into the targeting vector will require sequential Cold Fusion cloning reactions, one for each homology arm.

For standard cloning strategies, please refer to the following protocol for more details. Since there are two homology arms that need to be cloned, the cloning will need be done in a sequential manner except for the initial PCR reaction.

## 1. PCR of homology arms

- a) PCR amplify homology arms using proofreading PCR enzymes (e.g. Phusion) with primers containing compatible restriction site(s) for 5' homology arm and 3' homology arm depending on the chosen targeting vector
- b) Run out PCR products on 1.5-2% agarose gel and extract correct size bands using a suitable gel extraction kit

## 2. Ligation of insert into vector

- a) Digest the targeting vector using appropriate restriction enzymes for the homology arm being cloned and gel-purify the digestion product
- b) Dilute gel-purified, digested vector to 10 ng/ $\mu$ l
- c) Set up 10  $\mu$ l ligation reactions for each control and test samples as below:

<u>Volume</u>	<u>Item</u>
1.0 µl	HR Targeting vector
7.0 µl	5' or 3' homology arm PCR product (~70 ng) or water control
1.0 µl	10X T4 DNA ligase buffer
1.0 µl	T4 DNA Ligase (40 U/µl)
<b>10.0 µl</b>	<b>Total Reaction Volume</b>

- d) Incubate reactions at 25°C for 1-2 hours (sticky-end ligation) or O/N at 16°C (for blunt-end ligation)

### 3. Transform *E. coli* with the ligation product

Transform competent cells (with a transformation efficiency of at least  $1 \times 10^9$  colonies/µg pUC19) with the whole ligation reaction (10 µl) following the protocol provided with the competent cells. Plate the transformed bacteria on 50 µg/ml ampicillin agar plates

### 4. Identify clones with the correct insert

- a) Depending on the ratio of colony numbers for the test sample vs. the negative control sample, randomly pick 2 or 3 colonies and grow each clone in 100 µl of LB Broth with 100 µg/ml ampicillin at 37°C for 2 hours with shaking.
- b) Use 1 µl of each bacterial culture for screening homology arm inserts by PCR and continue to grow the culture for another 4 hours. Store the culture at 4°C.
- c) Prepare a PCR Master Mix with PCR primers flanking the insert:

1 rxn		10 rxn		Composition
0.5	μl	5	μl	PCR primer 1 (10 μM)
0.5	μl	5	μl	PCR primer 2 (10 μM)
0.5	μl	5	μl	50X dNTP mix (10 mM of each)
2.5	μl	25	μl	10X PCR Reaction Buffer
19.5	μl	195	μl	Nuclease-free water
0.5	μl	5	μl	Taq DNA polymerase (approx. 5 U/μl)
<b>24.0</b>	<b>μl</b>	<b>240</b>	<b>μl</b>	<b>Total volume</b>

d) Mix the master mix very well and aliquot 24 μl into each well of 96-well PCR plate or individual tubes.

e) Add 1 μl of each bacterial culture from step (b) into each well (or tube).

f) Proceed with PCR using the following program:

94°C, 4 min	1 cycle
94°C, 0.5 min, then 68°C, 1 min/1 kb*	25 cycles
68°C, 3 min	1 cycle

\* Depending on the size of final PCR product, use a shorter or longer time.

g) Take 5 μl of the PCR reaction and run it on a 1.5-2% agarose/EtBr gel in 1X TAE buffer to identify clones with correct insert.

h) Grow a positive clone containing insert in an appropriate amount of LB-ampicillin broth, and purify the construct using an endotoxin-free plasmid purification kit. Sequence verification of the insert is optional.

i) Repeat the entire process with remaining homology arm to build the final targeting construct

## **E. Protocol for Co-Transfection of Targeting Vector with TALEN or Cas9 plasmids and Characterization of Cells**

### Notes:

- I. Depending on the cell type being transfected, please choose a transfection protocol which results in maximal transfection efficiencies. For adherent cell lines such as HEK293T, passive transfection methods using cationic lipid-based methods (e.g. Lipofectamine 2000, FuGene HD) work very well in transfection of targeting and TALEN/Cas9 plasmids. For other types of cells such as primary, stem, or suspension cells, we suggest transfection using electroporation methods (NucleoFection or Neon) for optimal results.
  - II. The plasmids should be mixed well in minimal serum/no antibiotic media + cationic lipid transfection reagent or electroporation buffer to maximize efficiency of delivery.
  - III. For selection of target cells, we strongly recommend testing different concentrations of puromycin on untransfected cells to determine the optimal concentration of selection agent – which is kills ~90-100% of cells within 48-72 hours after drug administration.
- 1) Plate 200,000 to 300,000 cells (e.g. 293T cells) into a single well of a 6-well plate in 2 ml of appropriate growth medium. Include a single well of cells as negative control (which can be non-relevant plasmid

DNA)

- 2) Next day, or when cells are 50-60% confluent, co-transfect target cells with either TALEN or Cas9 plasmid(s) and the HR targeting vector using a suitable transfection reagent following the manufacturer's recommended protocol for 6-well plates. The use of reduced or serum-free media containing no antibiotics to dilute the vector/transfection complex is highly recommended.

**Note: For 293T cells, we suggest 0.5 µg of each TALEN (1 µg total for pair) or 0.5 µg of SBI's Cas9 SmartNuclease vector in conjunction with 1 µg of the HR targeting vector into cells for efficient cleavage and HR reaction. For other cell types, we suggest optimizing the amounts and ratios of Cas9 or TALEN to targeting vector for optimal results in a target cell line.**

- 3) Allow at least 12 hours before changing transfection media to complete growth media
- 4) Assay for positive HR events 48-96 hours after co-transfection. Select cells with insertion of targeting vector using FACS-based sorting of RFP or antibiotic selection markers Puromycin or Hygromycin depending on the HR donor vector. If using antibiotic selection markers, select cells for a minimum of 5 days prior to further characterization. Cutting efficiency of TALEN/Cas9 can be measured by Surveyor Nuclease Assay and HR efficiency by FACS based on % of RFP positive cells.
- 5) Genotyping of selected cells can be done by PCR of genomic DNA and vector for 5' and 3' homology arms. In brief, design PCR primer pairs with one of the

primers in the genomic DNA region and the other located in the HR vector. Both 5' and 3' integration ends can be confirmed for complete verification of insertion. In addition, the junction between the targeted exon and the ORF for GFP or T2A can be sequenced to confirm successful in-frame fusion.


- 6) Selected cells can be clonally expanded for further phenotypic characterization

## V. PrecisionX™ Targeting Vectors By Application: Gene Editing

### A. Overview of Gene Editing

One of the most challenging applications for genome engineering using HR-mediated approaches is gene editing, or the ability to correct or insert a single base (e.g. SNP) in presence of excess genomic background. While such applications have been demonstrated to correct single basepair mutations, they often leave residual sequences behind such as *loxP/FRT* sites (van der Weyden *et al.* 2002) or even selection markers, with potentially unknown effects on transcription of neighboring genetic elements (Meier *et al.* 2010). The ideal system would be amenable for HR-based correction of the mutation followed by a footprint-free removal of selection cassettes or any recombinase artifacts. Previous work on the *piggyBac* transposon system, originally derived from baculovirus infecting moths, has demonstrated that it can efficiently act as a transposon to insert and remove transgenes in mammalian cells by specific recognition of inverted terminal repeats flanking the transgene (Lacoste *et al.* 2010). By exploiting this aspect of the *piggyBac*

system, researchers have designed an HR targeting vector containing *piggyBac* inverted terminal repeats, in conjunction with ZFNs, to specifically target the *A1AT* locus in human iPSC cells to correct a single mutation responsible for  $\alpha$ -antitrypsin deficiency in a footprint-free manner (Yusa *et al.* 2011).



**Figure 22. Schematic Diagram of PBHR100PA-1 HR Targeting Vector [MCS1-5'PB ITR-EF1 $\alpha$ -GFP-T2A-Puro-T2A-*hsvTK*-pA-MCS2-3' PB ITR]**

SBI's PBHR100A-1 HR-based gene editing vector (Fig. 22) has been designed to facilitate such precision correction of target genomes in a footprint-free manner, leaving only the native sequence behind after excision of the selection cassette when utilized in conjunction with ZFN, TALEN, or Cas9 systems.

By cloning of homology arms at the 5' and 3' end of the cutting site, and one of the homology arms containing the desired correction, this vector can be co-transfected into target cells with ZFN, TALEN, or Cas9 vectors for HR-based placement into the target genome. The existing GFP and Puro positive selection markers driven by a robust EF1 $\alpha$  promoter can be used to select for cells that are positive for HR, followed by the use of the Excision – only SuperPiggybac transposase enzyme (Cat.

#PB220PA-1) to remove the selection cassette and any residual elements between the inverted terminal repeat elements flanking the selection cassette. In addition, a *hsvTK* negative selection marker is present in the selection cassette to select out any cells that have not excised the cassette when ganciclovir or fialuridine is added to the cells.

Alternatively, SBI's Gene Knock-Out vectors can also be utilized for gene editing approaches by using either the 5' or 3' homology arm for single nucleotide modifications and inserting the selection cassette into a nearby intronic sequence. This approach, while technically not "footprint-free", may be a more convenient method to generate lines harboring the desired change(s) with the flexibility of having a functional selection cassette located in a region of the genome that is not expected to have deleterious effect(s) on cell biology. After generating clonal lines, the selection cassette can also subsequently be excised using Cre recombinase, leaving behind only a single LoxP site. SBI's Knock-Out vectors are a good choice for many gene editing applications, and can be used following the guidelines for transfection of target cells in Section II.F (pg. 26).

## **B. Additional Materials Required (for PBHR100A-1 vector)**

- a) LB Agar and Broth containing 50 µg/ml ampicillin
- b) Any high-transformation efficiency RecA- and EndA- *E.coli* competent cells
- c) Dulbecco's Modified Eagle's Medium (D-MEM) high glucose with sodium pyruvate and glutamine (Invitrogen, Cat. # 11995073)
- d) Lipofectamine 2000 transfection reagent (Invitrogen, Cat. # 11668019)

- e) Qiagen EndoFree Plasmid Maxi Kit (Qiagen, Cat. # 12362)
- f) Zippy™ Plasmid MiniPrep Kit (Zymo Research, Cat. # D4019)
- g) Phusion High-Fidelity DNA Polymerase (NEB, Cat. # M0530S)
- h) Qiagen MinElute Gel Extraction Kit (Qiagen, Cat. #28604)
- i) Fetal Bovine Serum (Invitrogen, Cat. # 16000036)
- j) Penicillin/Streptomycin (Invitrogen, Cat. # 15070063)
- k) Surveyor Mutation Detection Kit (Transgenomics, Cat. #706020)
- l) Trypsin-EDTA (Sigma, Cat. # T3924)
- m) Excision-only SuperPiggyBac Transposase enzyme (System Biosciences, Cat. #PB220PA-1)
- n) Ganciclovir powder (Sigma, Cat. #G2536)
- o) 6 and 12-well Tissue Culture Plates and Related Tissue Culture Supplies
- p) Other specific media and additives specific for cell type of interest
- q) **\*\*Optional\*\*** - For difficult-to-transfect cells, the use of an electroporation system (e.g. Lonza's NucleoFector or Invitrogen's Neon system) is highly recommended

### **C. General Considerations (for PBHR100A-1 vector)**

We recommend propagation of the PBHR100A-1 targeting vector plasmid prior to starting the experiments. The plasmids can be transformed using standard conditions suitable in any high-quality RecA- and EndA- *E.coli* competent cell.

For cells transformed with any of the targeting vectors, we suggest plating 50-200 µl of transformed cells on fresh


ampicillin plates (50 µg/ml). Incubate the plates at 37°C overnight. Colonies picked from the transformation can be grown at 37°C overnight in ~200 ml of LB media containing 50 µg/ml of ampicillin. After overnight growth, plasmid DNA can be harvested from culture using an endotoxin-free DNA plasmid maxiprep kit.

For confirmation of the plasmid, we recommend performing restriction digestion analysis or direct sequencing to confirm integrity of the amplified plasmids.

## **D. Design and Cloning of Homology Arms into Vector (for PBHR100A-1 vector)**

In general, cloning of the 5' and 3' homology arms is best achieved through fusion cloning strategies such as SBI's Cold Fusion Kit. Please refer to the user manual for the Cold Fusion Kit for details on the cloning process (<http://www.systembio.com/molecular-tools/cold-fusion-cloning/literature>). The cloning of the homology arms into the PBHR100A-1 vector has some special considerations that need to be followed for successful cloning. The design necessitates that the 5' homology arm terminates with TTAA and the 3' homology arm begins with TTAA. This is purposely designed to mimic the insertion of a transgene cassette into the genome by the *piggyBac* transposase, which recognizes a TTAA site in the genome and inserts the cassette flanked by inverted tandem repeats into this site. The inserted cassette now will be flanked by two TTAA sites, one at the 5' end and at the 3' end. When the *piggyBac* transposase is added to excise the cassette, it will leave a single TTAA site behind in the genome. If this TTAA site does not have a homologous sequence in the genome nearby the desired correction site, it will be left in the target genome as an insertion and

will likely lead to frameshifting of the coding sequence. To prevent this from occurring, the homology arm will need to be carefully designed such that the TTAA site at the 5' and 3' homology arms are homologous to another TTAA site (or sequence that can be altered to a TTAA) nearby the correction site as shown below (adapted from Yusa *et al.* 2011):



ZFN recognition site  
TTAA conversion site

5' Homology Arm

5' ....NNNNNNNNNGCTGTG**TTAA** -3'

3' Homology Arm

5' – **TTAA**CCATCGAC**GAG**AAAAGGGACTGAAGCTGCT.....3'

Correction sequence in HR arm

After co-transfection with ZFN, TALEN, or Cas9 vectors and HR of the PBHR vector, this is the expected result:

5'....NNNNNNNNNGCTGTG**TTAA****ITR**EF1a-GFP-T2A-Puro-T2A-hsvTK**ITR****TTAA**CCATCGAC**GAG**AAAAGGGACTGAAGCTGCT.... 3'

Piggybac Recognition Site  
Selection Cassette

Lastly, when the Excision-only SuperPiggyBac transposase is used to excise out the selection cassette, it will recognize the ITR elements and lift out the entire ITR and selection cassette along with a single TTAA site, and stitch the remaining sequence together as shown below:

5'-  
 NNNNNNNNNGCTGTGTTAAACCATCGACGAGAAAGGGACTGAAGCT  
 GCT - 3'

Corrected sequence

In this particular example, the TTAA in the 5' homology arm is inserted into a site ~10bp upstream of the correction site whose nucleotide sequence can be altered from CTG to TTA as a silent mutation. The TTAA in 3' the homology arm shares the homology as well, and contains the corrected base downstream. The net result is seamless, footprint-free correction of the genomic sequence.

**\*\*IMPORTANT NOTE\*\***

When using ZFN, TALE-Nuclease, or Cas9-mediated cleavage of the target genome, please make sure that the recognition sequence is **NOT PRESENT** in its entirety in one of the homology arms in the PBHR vector. If this occurs, the donor vector along with the target genome will be cut, leading to reduced or minimal HR efficiency. Please make sure that the recognition sequence spans the 5' and 3' arms of the PBHR vector or mutate the recognition site present in the donor vector to minimize any possible cutting of the donor vector.


The detailed protocol for cloning of homology arms is shown below:

**1) Preparation of 5' Homology Arm**


- a. PCR the 5' HA and add any restriction site you choose from the 5' MCS except *Bsa*I to the forward primer. Add a *Bsa*I site and the sequence GTTAA to the 3' end of the reverse HR primer. Homology arms

should be between 0.5kb to 1kb. The 5' HA cannot contain a *BsaI* site.

### 5' HA forward primer sequence



### 5' HA reverse primer sequence



- b. Run out the PCR product on an agarose gel.
- c. Cut of the band of the correct size and isolate the fragment.
- d. Cut the fragment with the 5'MCS enzyme you chose and with *BsaI*.
- e. Clean up the reaction and quantify the DNA.

## 2) Preparation of 3' Homology Arm

1. PCR the 3'HA and add the *BspQI* site to the forward primer. Add any site you wish from the 3' MCS except *BspQI* to the reverse primer. Homology arms should be between 1 and 0.5 Kb. The 5' HA cannot contain a *BspQI* site.



1. Cut the PBHR vector with the 5' MCS enzyme you chose and with BsaI.
2. Run out the reaction on an agarose gel.
3. Cut out the band of the correct size and purify.
4. Quantitate the DNA.
5. Ligate in the 5'HA fragment. Use a recombination deficient competent bacterial strain for the transformation.
6. Screen colonies for the clone of interest.
7. Cut the PBHR vector containing the 5'HA with the enzyme you chose in the 3' MCS and with BspQI.
8. Run out the reaction on an agarose gel.
9. Cut out the band of the correct size and purify.
10. Quantitate the DNA.
11. Ligate in the 3'HA fragment. Use a recombination deficient competent bacterial strain for the transformation.
12. Screen colonies for the clone of interest.

## **E. Protocol for Gene Editing with TALEN/Cas9 system and piggyBac HR vector**

### Notes:

- I. Depending on the cell type being transfected, please choose a transfection protocol which results in maximal transfection efficiencies. For adherent cell lines such as HEK293T, passive transfection methods using cationic lipid-based methods (e.g. Lipofectamine 2000, FuGene HD) work very well in transfection of targeting and TALEN/Cas9 plasmids. For other types of cells such as primary, stem, or suspension cells, we suggest transfection using electroporation methods (NucleoFection or Neon) for optimal results.
  - II. The plasmids should be mixed well in minimal serum/no antibiotic media + cationic lipid transfection reagent or electroporation buffer to maximize efficiency of delivery.
  - III. For selection of target cells, we strongly recommend testing different concentrations of puromycin on untransfected cells to determine the optimal concentration of selection agent – which is kills ~90-100% of cells within 48-72 hours after drug administration.
- 1) Plate 200,000 to 300,000 cells (e.g. 293T cells) into a single well of a 12-well plate in 1 ml of appropriate growth medium. Include a single well of cells as negative control (which can be non-relevant plasmid

DNA)

- 2) Next day, or when cells are 50-60% confluent, co-transfect target cells with either TALEN or Cas9 plasmid(s) and the HR targeting vector using a suitable transfection reagent following the manufacturer's recommended protocol for 12-well plates. The use of reduced or serum-free media containing no antibiotics to dilute the vector/transfection complex is highly recommended.

**Note: For 293T cells, we suggest 0.5 µg of each TALEN (1 µg total for pair) or 0.5 µg of SBI's Cas9 SmartNuclease vector in conjunction with 1 µg of the HR targeting vector into cells for efficient cleavage and HR reaction. For other cell types, we suggest optimizing the amounts and ratios of Cas9 or TALEN to targeting vector for optimal results in a target cell line.**

- 3) Allow at least 12 hours before changing transfection media to complete growth media
- 4) Assay for HR events 48-96 hrs after co-transfection. Cutting efficiency of TALEN/Cas9 can be measured by Surveyor Nuclease Assay and HR efficiency by FACS based on % of GFP positive cells. Puromycin selection can start 5 days after co-transfection.
- 5) Continue growing cells in puromycin until distinct colonies begin to form (5-7 days depending on cells)
- 6) When large enough, pick colonies to individual well of a 6-well plate under puromycin selection
- 7) Isolate genomic DNA from each line and screen by PCR for clones that have the PBHR vector integrated at the correct location. Design one genomic primer

outside each homology arm and on in the PBHR vector.

- 8) Transfect lines with correct integrations with the Excision-only SuperPiggyBac transposase using recommended amount for transfection (~0.2ug/6 well plate).

**Note: To increase the chance for excision, it is recommended to switch to media w/o puromycin for 3-4 days which decreases selection pressure on cells to retain the PB cassette**

- 9) Split cells 48-72 hours after transfection with the transposase
- 10) 24 hours later, switch to media containing 2-10µg/ml of ganciclovir to select out cells containing the expression cassette
- 11) Continue growing cells in ganciclovir media until distinct colonies begin to form (~5-10 days depending on cells)
- 12) Remaining colonies should be GFP<sup>-</sup>, puro<sup>s</sup> and ganciclovir<sup>r</sup>.
- 13) Isolate genomic DNA from cell lines and use PCR to confirm the absence of the PBHR vector and to amplify the region containing the correction.
- 14) Confirm correction by DNA sequencing

## VI. References

Hockemeyer D *et al.* Efficient targeting of expressed and silent genes in human ESCs and iPSCs using zinc-finger nucleases. *Nat Biotechnol.* 2009; 27:851–857.

Boch J *et al.* Breaking the code of DNA binding specificity of TAL-type III effectors. *Science* 2009 Dec 11;326(5959):1509-12.

Bhaya D *et al.* CRISPR-Cas systems in bacteria and archaea: Versatile small RNAs for adaptive defense and regulation. *Annu Rev Genet.* 2011;45:273-97.

Jinek M *et al.* A programmable dual-RNA-guided DNA endonuclease in adaptive bacterial immunity. *Science* 2012 Aug 17;337(6096):816-21.

Choulika A *et al.* Induction of homologous recombination in mammalian chromosomes by using the *I-SceI* system of *Saccharomyces cerevisiae*. *Mol Cell Biol.* 1995 Apr;15(4):1968-73.

Sadelain M *et al.* Safe harbours for the integration of new DNA in the human genome. *Nat Rev Cancer* 2011 Dec 1;12(1):51-8.

Van der Weyden L *et al.* Tools for targeted manipulation of the mouse genome. *Physiol Genomics* 2002 Dec 3;11(3):133-64.

Meier ID *et al.* Short DNA sequences inserted for gene targeting can accidentally interfere with off-target gene expression. *FASEB J.* 2010; 24:1714–1724.

Lacoste A *et al.* An efficient and reversible transposable system for gene delivery and lineage-specific differentiation in human embryonic stem cells. *Cell Stem Cell* 2009 Sep 4;5(3):332-42.

Yusa K *et al.* Targeted gene correction of  $\alpha$ 1-antitrypsin deficiency in induced pluripotent stem cells. *Nature* 2011 Oct 12;478(7369):391-4.

## VII. Technical Support

For more information about SBI products or to download manuals in PDF format, please visit our website:

<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

[tech@systembio.com](mailto:tech@systembio.com)

650-968-2200

## VIII. Licensing and Warranty Statement

### Limited Use License

Use of the PrecisionX™ HR Targeting Vectors (*i.e.*, the “Product”) is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for stem cell research.

SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

### **Limited Warranty**

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2014 System Biosciences (SBI).