



Exosome ELISA Complete Kits

EXOEL-CD9A-1, EXOEL-CD63A-1, EXOEL-CD81A-1

User Manual

See PAC for Storage Conditions for Individual Components

Version 12
4/17/2017

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the License and Warranty Statement contained in this user manual.

Contents

Product Description	1
List of Components	1
Equipment to be supplied by user	1
Storage	2
General Information	2
Protocol	2
For exosomes precipitated with one of SBI's ExoQuick reagents	2
For exosomes precipitated by UltraCentrifugation	2
Exosome protein standard curve	2
ELISA procedures	3
Before starting	3
ELISA assay	4
Next Steps and Related Products	5
Example Data and Applications	5
Technical Support	6
Licensing and Warranty Statement	6

Product Description

The ExoELISA kit is designed as a direct Enzyme-Linked ImmunoSorbent Assay (ELISA). The lysed exosome particles and their proteins are directly immobilized onto the wells of the microtiter plate. After binding, the wells are coated with a blocking agent to prevent non-specific binding of the primary detection antibody. The detection (primary) antibody is added to the wells for binding to a specific antigen (e.g. CD63) protein present in the exosome lysates. A Horseradish Peroxidase enzyme linked secondary antibody (Goat anti-Rabbit) is used for signal amplification and to increase assay sensitivity. A colorimetric substrate (Extra-sensitive TMB) is used for the assay read-out. The accumulation of the colored product is proportional to the specific antigen present in each well. The results are quantitated by a microtiter plate reader at 450 nm absorbance.

List of Components

ExoELISA kit Components	Amount
Exosome binding buffer	20 ml
20X Wash buffer	18 ml
Blocking buffer	30 ml
ExoELISA protein standard (Check the PAC for Lot-specific information)	400 µl
Exosome specific primary antibody (CD63) (Check the PAC for Lot-specific information)	2 x 25 µl
Exosome validated secondary antibody (Goat anti-Rabbit HRP)	10 µl
ELISA substrate (Super-sensitive TMB)	6 ml
Stop buffer	6 ml
96 well ExoELISA plate (12x8 well strips)	1 plate

Equipment to be supplied by user

1. Microtiter plate sealing film/cover
2. 37°C incubator
3. Microtiter plate shaker
4. Microtiter plate spectrophotometer with 450 nm absorbance capability
5. Multichannel pipets (recommended)

Storage

The kits are shipped on blue ice and individual components should be stored according to the information below upon receipt. Properly stored kits are stable for 3 months from the date received.

Item	Storage Temperature	Notes
*Protein Standards	-20°C	Make single use aliquots. AVOID FREEZE-THAW CYCLES!
Other Kit Components	+4°C	

General Information

The standards included in this kit are extremely temperature sensitive, so that fluctuations or freeze-thaw cycles will negatively impact their function. Please take extreme care to make single use aliquots.

Information on how to dilute the standards is contained in the Product Analysis Certificate (PAC) and is lot-specific. Please refer to the PAC for each set of standards being used in this kit. Standards should be thawed and diluted immediately before using them.

Protocol

For exosomes precipitated with one of SBI's ExoQuick reagents

1. For exosomes precipitated with ExoQuick, resuspend the pellet in 200 μ L Exosome Binding Buffer, and vortex for 15 sec.
2. Incubate at 37 °C temperature for 20 minutes to liberate exosome proteins
3. Centrifuge at 1500 \times g for 5 minutes to remove all residual precipitation solution
4. Transfer supernatant to new centrifuge tube on ice
5. Exosome protein is now ready for immobilization onto micro-titer plate

For exosomes precipitated by UltraCentrifugation


1. Carefully remove all traces of supernatant, keep exosome pellet
2. Resuspend exosome pellet in 200 μ l of Exosome binding buffer and vortex 15 seconds
3. Incubate at room temperature for 10 minutes and keep on ice.
4. Exosome protein is now ready for immobilization onto micro-titer plate

Exosome protein standard curve

A standard curve should be prepared each time the assay is performed

1. Thaw ExoELISA protein standard on ice
2. Dilute ExoELISA protein standard by performing serial dilutions with Exosome Binding buffer in microcentrifuge tubes
3. Below is an example of how to create a standard curve. The number of exosomes that the standard curve corresponds to is lot-specific and will be contained in the PAC that came with the kit you purchased. Example dilutions for standard curve:

Tube	# Exosomes	Dilution factor	ExoELISA protein standard	Exosome binding buffer
0	0	Blank	-	60
1	1.35×10^{10}	1	60 μ l	-
2	6.75×10^9	1:2	60 μ l	60 μ l
3	3.37×10^9	1:4	60 μ l	60 μ l
4	1.68×10^9	1:8	60 μ l	60 μ l
5	8.44×10^8	1:16	60 μ l	60 μ l
6	4.21×10^8	1:32	60 μ l	60 μ l
7	2.10×10^8	1:64	60 μ l	60 μ l



ELISA procedures

Before starting

1. The ExoELISA micro-titer plate is provided in a convenient 8 well X 12 strip format. We recommend using at least one strip for the standard curve and additional strips depending on the number of samples tested. Unused 8-well strips can be removed and stored at room temperature for later use.
2. Make sure to warm the Super-sensitive TMB ELISA substrate to room temperature before adding to the ELISA plate wells in step #10.

3. Dilute stock 20X Washing buffer into 1X working Wash buffer with purified water (each 8-well strip requires approximately 10 ml of 1X Washing solution)

ELISA assay

1. Add 50 μ l of prepared protein standards and exosome protein sample to the appropriate well of the micro-titer plate

2. Cover plate with sealing film/cover

3. Incubate the plate at 37°C from 2 hours to overnight (recommended)

4. After incubation step, invert the plate to empty all contents.

5. Wash the plate 3 times for 5 minutes each with 100 μ l 1X Wash buffer

- A micro-titer plate shaker is recommend for all subsequent the washing and incubation steps
- Residual liquid should be removed by hard-tapping the plate on fresh paper towels, while taking care not to let the wells dry out completely

6. Dilute Exosome specific primary antibody (CD63, CD9 or CD81) 1:100 in 1X blocking buffer and add 50 μ l of to each well and incubate at room temperature for 1 hour with shaking

7. Wash the plate 3 times for 5 minutes each with 100 μ l 1X Wash buffer

8. Dilute Exosome validated secondary antibody 1:5,000 1X blocking buffer and add 50 μ l to each well and incubate at room temperature for 1 hour with shaking

9. Wash the plate 3 times for 5 minutes each with 100 μ l 1X Wash buffer

10. Add 50 μ l of Super-sensitive TMB ELISA substrate and incubate at room temperature for 15 to 45 minutes with shaking

- 15 to 45 minutes substrate incubation time is optimized for the recommended exosome protein standard curve
- Further optimization maybe required by the user for individual sample.

11. Add 50 μ l of Stop buffer to provide a fixed endpoint for the assay


Note that the initial color of a positive sample is blue and the color changes to yellow when Stop buffer is added

12. Quantitate results with a spectrophotometric plate reader at 450 nm absorbance.

Next Steps and Related Products

Application	Related Products	Website links
Protein Characterization of Exosomes		
Western blotting	Exosome antibodies	https://www.systembio.com/micrna-research/exosome-antibody/exoab
Antibody Arrays	ExoCheck Assays	https://www.systembio.com/micrna-research/exosome-antibody-arrays
Quantification of Exosomes		
Quantification of exosomes	EXOCET Assays	https://www.systembio.com/micrna-research/exosome-antibody/exocet-assay
RNA extraction from Exosomes		
RNA extraction and profiling	SeraMir kits	https://www.systembio.com/micrna-research/seramir-exosome-rna-profiling/overview

Example Data and Applications



Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site: <http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200
(888) 266-5066 (Toll Free)

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com

Technical Support: tech@systembio.com

Ordering Information: orders@systembio.com

Licensing and Warranty Statement

Limited Use License

Use of the Exo-ELISA kits (*i.e.*, the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This

limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2017 System Biosciences (SBI), All Rights Reserved



System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200
Toll Free (888) 266-5066
Fax: (650) 968-2277

E-mail:
General Information: info@systembio.com
Technical Support: tech@systembio.com
Ordering Information: orders@systembio.com