

# **ExoQuick-TC<sup>™</sup>** Exosome precipitation

Exosomes are 60 –120 nm membrane vesicles secreted by most cell types *in vivo* and *in vitro*. Exosomes are found in blood, urine, amniotic fluid, malignant ascite

fluid, tissue culture media and contain distinct subsets of microRNAs and proteins

depending upon the cell-type from which they are secreted. SBI's ExoQuick-TC exosome precipitation reagent is a special polymer formulation distinct from the

original ExoQuick reagent for serum. ExoQuick-TC has been optimized for exosome

isolation from media and urine samples. This technology makes microRNA and

Time saving, cost-effective

solution for studying

exosomes from culture media

Optimized one-step solution for rapidly isolating exosomes from tissue culture media and urine for biomarker analysis



# Highlights


- Precipitate exosomes from cell culture media or urine samples easily
- Exosomes are released from tumors in high abundance
- Exosome cargo reflects the origin and physiological state of the source cells
- MicroRNAs are found in high abundance in circulating exosomes
- Discover novel disease-specific biomarkers

#### Isolate exosomes with ease

- No time-consuming ultracentrifugation
- Less expensive than costly DynaBeads
- More effective than any other method

0.09

• Use as little as 5ml of cell culture media or urine



protein biomarker discoveries simple, reliable and quantitative.


#### 0.08 0.07 **αCD9** 0.06 **OD450** blot 0.05 Western 0.04 0.03 0.02 CD9 (~22 kDa) 0.01 0 0 5 10 15 20 25 30 35 40 45 Amount of Urine exosomes added (µl)

### **ExoQuick-TC precipitates urine exosomes**

Ten milliliters of normal human urine was combined with 2ml ExoQuick-TC to precipitate urine exosomes. The exosome pellet was resuspended 175µl buffer and increasing amounts of the exosome suspension was loaded onto an ELISA-ready plate. The CD9 protein was detected using SBI's rabbit anti-CD9 primary antibody and SBI's HRP-conjugated secondary goat anti-rabbit antibody. The size of urine CD9 proteins was determined using Western blot analysis with the same set of antibodies (see inset).





#### NanoSight Analyses on Exosomes from Tissue Culture Media and Urine



# Media from HT1080 cells


Human HT1080 lung sarcoma cell line was cultured in conditioned media (serum-free) for 72 hours. Ten milliliters of the media was combined with 2ml ExoQuick-TC to pellet the exosomes overnight. The exosome pellet was resuspended in 1ml PBS, diluted 1:40 and visualized on the NanoSight LM10 instrument. The analysis shows that ExoQuick isolated 133nm exosomes with a recovery of 1.74 x 10^9 particles/ml.






# Media from HEK293 cells


Human embryonic kidney (HEK) cell line was cultured in conditioned media (serum-free) for 72 hours. Ten milliliters of the media was combined with 2ml ExoQuick-TC to pellet the exosomes overnight. The exosome pellet was resuspended in 1ml PBS and visualized on the NanoSight LM10 instrument undiluted. The analysis shows that ExoQuick isolated 137nm exosomes with a recovery of 3.33 x 10^8 particles/ml.





### **Human Urine**

A Normal human urine sample was used. Five milliliters was combined with 2.5ml ExoQuick-TC to pellet the exosomes overnight. The exosome pellet was resuspended in 1ml PBS, diluted 1:50 and visualized on the NanoSight LM10 instrument. The analysis shows that ExoQuick isolated 107nm exosomes with a recovery of 4.80 x 10^9 particles/ml.



### We Also Offer Custom Services

System Biosciences offers a wide-range of custom services to support your research, allowing you to spend less time making tools, and more time making discoveries. To learn more, visit our website at www.systembio.com/service or call us at 888-266-5066.



**System Biosciences, Inc.** 265 North Whisman Rd. Mountain View, CA 94043 Toll Free: 888.266.5066 Fax: 650-968-2277 Email: info@systembio.com www.systembio.com © 2011 System Biosciences, Inc. All rights reserved. System Biosciences and the System Biosciences logo are trademarks of System Biosciences, Inc.