

Human and Mouse LncProfilers™ qPCR Array Kits

Quantitate long non-coding RNAs (lncRNAs) by real-time qPCR

Cat # RA900/910A-1, RA930A-1

User Manual

Store kit at -20°C on receipt

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.

Contents

I. Introduction and Background	
A. Overview	2
B. LncProfiler qPCR profiler workflow	6
C. How LncProfiler cDNA synthesis works	7
D. List of components	8
II. Protocol	
A. LncProfiler cDNA reaction set up	9
B. Real-time qPCR Reaction setup	10
C. How the lncRNA specific primers are designed	11
D. LncProfiler qPCR array contents	13
E. Endogenous reference controls	15
III. Sample Data and Quality Control	
A. LncProfiler qPCR sample data	16
B. Specificity Tests	19
IV. Troubleshooting	21
V. lncRNA Selected References	22
VI. Technical Support	23
IX. Licensing and Warranty Statement	24

I. Introduction and Background

A. Overview

For the last few decades of the 20th century, the underlying dogma of molecular biology has been that the purpose of RNA is to direct the assembly of proteins from amino acids. A few exceptions to this paradigm were known (for example, ribosomal RNA and transfer RNA, which are functional RNA macromolecules that do not code for protein, or viral genomes that exist as or pass through an RNA phase as part of total genome replication).

Non-coding RNAs (ncRNAs) include the familiar “housekeeping” RNAs (ribosomal, transfer, small nuclear, and small nucleolar RNAs) and the thousands of regulatory RNAs that are the subject of recent intense exploration. Regulatory ncRNAs are arbitrarily classified by size: small ncRNAs (sncRNA) being less than 200 bp, and long ncRNAs (lncRNA) greater than 200 bp. The sncRNAs include other sub-classifications: microRNA (miRNA), endogenous small inhibitor RNA (endo-siRNA), and PIWI-associated RNA (piRNA).

The roles of lncRNAs in the regulation of gene expression and organismal development are diverse and just beginning to be discovered. Biological processes dependent upon lncRNAs include imprinting and gene dosage regulation, stem cell pluripotency, embryonic development and segmentation, hematopoiesis, and neural cell fate determination. LncRNAs may employ a number of mechanisms to impact gene expression via *cis* and *trans* processes.

Gene imprinting: While the function of parental gene imprinting is still unclear, lncRNAs have been found to participate in imprinting processes. Imprinting Control Regions (ICRs) are DNA regions that are differentially methylated depending on their parental origins. Unmethylated ICRs cause specific expression of nearby lncRNAs, which then suppress neighboring genes in *cis*. *Air* and *Kcnq1ot1* are examples of lncRNAs that cause suppression of paternally inherited genes.

Gene dosage and X chromosome inactivation: The discovery of *Xist* was one of the defining moments in the realization that ncRNAs can have profound roles in the control of gene expression. *Xist* is an lncRNA that suppresses the inactive, non-coding X chromosome (Xi) in female cells. In all, 7 ncRNAs are found as part of the X inactivation center on the X chromosome, including *Xist*. Initially, *Xist* and its antisense transcript *Tsix* are expressed on both X chromosomes. However, *Tsix* expression continues on the X that will remain active (Xa) and this activity recruits DNMT3A to suppress *Xist* from being transcribed on Xa. On Xi, it is *Tsix* that is suppressed, potentially via another lncRNA that is part of the X inactivation center, *Jpx*. With *Tsix* suppressed, the protein PRC2 is recruited to induce histone modification marks at the 5' end of *Xist*. This upregulates *Xist* expression on Xi and causes further propagation of these silencing marks throughout Xi, which are maintained across the lifetime of the organism.

Embryonic development and segmentation: The expression of HOX genes is also regulated by lncRNAs. Some HOX-related lncRNAs operate in *cis*, having either enhancing or repressive effects. However, some like the human HOTAIR work in *trans*, and may function as scaffolds for histone-modifying complexes. It is not clear if *trans*-acting lncRNAs like HOTAIR are involved in the process of identifying the DNA sites to which the complexes will be recruited, or if that function is retained by the protein elements of the complex.

Stem cell pluripotency: The promoters of more than 100 lncRNAs are bound by stem cell factors. Disruption of these lncRNAs can alter cell differentiation. lincRNA-RoR is involved in the reprogramming of fibroblasts back to a pluripotent state. Thus, lncRNAs are likely to play important roles in both normal development and processes that require maintenance of adult stem cell pools.

Cell fate determination: LncRNAs are implicated in cell fate determination events in multiple cell lineages, including the nervous system. *TUG1* is an lncRNA that may enhance rod gene expression and suppress cone gene expression in the developing eye. *Evf2* is a mouse lncRNA that appears to have both *cis* and *trans* effects to repress *Dlx5*, *Dlx6*, and *Gad1* during forebrain development.

Dysregulated expression of lncRNAs has been shown to be associated with a broad range of diseases such as Alzheimer's, psoriasis and many cancers. Studying the expression patterns of lncRNAs will be a crucial method to understanding the roles they play in many model systems. SBI has built a sensitive, accurate and robust qPCR array to enable researchers to closely profile the expression changes in the top lncRNAs known to date.

This manual provides details and information necessary to use the LncProfiler™ Kit to tag and convert small non-coding RNAs into detectable and quantifiable cDNAs. The system allows for the ability to quantitate dynamic fold differences of lncRNAs across 20 separate experimental RNA samples. The array plate also includes 5 endogenous RNA assays as normalization signals. To ensure optimal results, please read the entire manual before using the reagents and material supplied with this kit.

These LncProfiler qPCR Array comes with all the reagents necessary to tag all long non-coding as well as small RNAs from 20 different total RNA samples into quantifiable cDNA. The kits include assays in pre-formatted plates for well-annotated human, lncRNAs with three endogenous reference RNA controls on each plate. All of the lncRNAs on the qPCR array have validated primer sets for well-annotated lncRNAs that are registered in the lncRNA database created by Dr. John Mattick (www.lncrnadb.org).

Potential functions of lncRNAs

To date, lncRNAs have been found to exhibit a wide range of functions ranging from signaling, serving as molecular decoys, guiding ribonucleoprotein complexes to specific chromatin sites and also participating as scaffolds in the formation of complexes.

I. Signaling

The transcription of certain lncRNAs is very tissue and temporal specific. Their expression can be in response to certain stimuli, such as cellular stress and temperature. Thus, lncRNAs can serve as molecular signals and can act as markers of functionally significant biological events. Examples include imprinting lncRNAs **XIST**, **AIR** and **Kcnq1ot1**.

II. Decoys

The molecular decoy type of activity takes place when specific lncRNAs are transcribed and then bind to and titrate away protein factors. Decoy lncRNAs can "sponge" protein factors such as transcription factors and chromatin modifiers. This leads to broad changes in the cell's transcriptome. Example is **MALAT1**.


III. Guides

lncRNAs can be molecular guides by localizing particular ribonucleoprotein complexes to specific chromatin targets. This activity can cause changes in gene expression either in *cis* (on neighboring genes) or in *trans* (distantly located genes) that cannot be easily predicted by just the lncRNA sequence itself. Some example lncRNAs that act as guides are **XIST** and **HOTTIP**.


IV. Scaffolds

Assembly of complex protein complexes can be supported by lncRNAs, linking factors to together to form new functions. Some lncRNAs possess different domains that bind distinct protein factors that altogether, may impact transcriptional activation or repression. Some examples of scaffold lncRNAs are **HOTAIR** and **ANRIL**.


I. Signaling




II. Guides




III. Decoy




IV. Scaffolds





B. LncProfiler qPCR profiler workflow



C. How the LncProfiler cDNA synthesis works



The initial polyadenylation step greatly enhances cDNA synthesis yields of lncRNAs (over 100-fold) and enables the usage of small RNAs (like U6 and RNU43) to be included as reference controls on the qPCR array.



D. List of components

40 µl	5X PolyA Polymerase Buffer	10 µl Poly A Reaction
10 µl	PolyA Polymerase	(enough for 20 reactions)
20 µl	25 mM MnCl ₂	
30 µl	5 mM ATP	
10 µl	Oligo dT Adaptor	20 µl RT Reaction
80 µl	5X Reverse Transcriptase Buffer + Random Primer Mix	(enough for 20 reactions)
20 µl	Reverse Transcriptase	
30 µl	0.1 M Dithiothreitol (DTT)	
50 µl	dNTP Mix	
3500 µl	2X SYBR Green master Mix *	Enough for first 2 profiles using 2X Maxima SYBR Green from Fermentas
	Array Primers, dried down in Primer plate (400 pmoles); resuspend in 44µl RNase-free Water	(enough for 20 profiles)
1.2 ml	RNase-free Water	

*Catalog# RA910A-1 contains all of the components listed above.

*Catalog# RA900A-1 has all of the above components except for the 2X SYBR Green reagent for first 2 profiles.

The kit is shipped on blue ice and should be stored at -20°C upon arrival. Properly stored kits are stable for 1 year from the date received. The oligonucleotides for the specific lncRNAs are dried-down in the wells of the optical qPCR plates. Resuspend in 10µl RNase-free water. SBI recommends using the LncProfiler qPCR array with the following SYBR Green reagents:

- 2X Maxima® SYBR Green with Rox (Cat# K0223) from Fermentas (highly recommend).
- Power SYBR Master Mix® (Cat. #s 4368577, 4367650, 4367659, 4368706, 4368702, 4368708, 4367660) from Applied Biosystems.

II. Protocol

A. LncProfiler cDNA reaction setup


(for 1 RNA sample to be assayed on qPCR 96-well plates)



It is important to start with total RNA that includes the lncRNA fraction. RNA input can be as low as 1-2 µg total. For optimum signals, perform the following.

➡ Dilute your RNA to ~200-400 ng/µl

Start: In a thin-walled PCR tube or PCR-compatible plate well combine:



Done!

* The lncRNA cDNAs can be stored at -20°C. For more sensitive applications, a single phenol:chloroform extraction with ethanol precipitation can be performed on the cDNA to remove proteins, unused dNTPs and primers, typically this is not necessary.

B. Mastermix qPCR Reaction Setup for 1 entire 96-well qPCR plate

To determine the expression profile for the lncRNAs under study, mix the following for 1 entire qPCR plate:

For 1 entire plate:

	1,750	µl	2X SYBR Green qPCR Mastermix buffer
+	20	µl	LncRNA cDNA (from Step A)
	1,730	µl	RNase-free water
	3,500	µl	Total

Aliquot 28µl of Mastermix per well in your qPCR Plate.

Resuspend Primers in Primer plate with **44µl RNase-free water** per well before use. (the primers are dried-down in the stock primer plate)

Then :

Load 2µl per well of each of the Primers from the Primer plate into your qPCR plate (well A1 into qPCR plate A1, etc.)

The Mastermix contents can be scaled up or down depending upon your experimental needs. If you want to perform the reactions in triplicate, scale up the cDNA synthesis reactions by 3-fold and add 3X the RNA input. Or, simply follow the above recipe three times for each of the qPCR plates you want to run as replicates. Once reagents are loaded into the wells, cover the plate with an optical adhesive cover and spin briefly in a centrifuge to bring contents to bottom of wells. Place plate in the correct orientation (well A1, upper left) into the Real-time qPCR instrument and perform analysis run.



*** Use a Multichannel pipette to load the qPCR plate with MasterMix and Primers:** Pour the Mastermix into a reservoir trough and use a 8 or 12 channel pipette to load the entire 96-well qPCR plate with the Mastermix. Then load the primers from the primer plate to the qPCR plate using a separate multichannel pipette.



2. Real-time qPCR instrument parameters

Follow the guidelines as detailed for your specific Real-time instrumentation. The following parameters tested by SBI were performed on an Applied Biosystems 7300/7500 Real-time PCR System but can also apply to any other 96-well systems. The details of the thermal cycling conditions used in testing at SBI are below. A screenshot from the Real-time instrument setup is shown below also. Default conditions are used throughout.


Create a detector:

1. Create a new Detector

2. Name the Detector (any name will do)

3. Select Reporter Dye as "SYBR Green"

4. Select Quencher Dye as (none)



Instrument setup:

qPCR cycling and data accumulation conditions:

Standard Protocol

1. 50°C 2 min.
 2. 95°C 10 min.
 3. 95°C 15 sec.
 4. 60°C 1 min.
- (40 cycles of Stage 3), data read at 60°C 1 min. Step.



An additional recommendation is to include a **Dissociation Stage** after the qPCR run to assess the T_m of the PCR amplicon to verify the specificity of the amplification reaction. Refer to the User Manual for your specific instrument to conduct the melt analysis and the data analyses of the amplification plots and Cycle Threshold (Ct) calculations. In general, Cycle thresholds should be set within the exponential phase of the amplification plots with software automatic baseline settings.

C. How the IncRNA-specific primers are designed for detection and quantitation in the qPCR array

SBI's LncProfiler is complete cDNA synthesis kit combined with a 96-well based qPCR assay set. The qPCR assays have been validated across numerous cell types for robust and specific performance. Some lncRNAs have endogenous polyA tails, while other lncRNAs do not. To enhance qPCR assay performance, the cDNA synthesis kit includes reagents to polyadenylate all lncRNAs before cDNA conversion with the oligo dT adaptor and random primers. SBI's lncRNA qPCR assays are derived from published primer sequences and others were designed in-house for robust performance and all amplicons are below 200 bp in size.

All qPCR assays are designed to detect Human lncRNAs that are annotated in Dr. John Mattick's lncRNA database.

<http://www.lncrnadb.org/>

Below is a screenshot of the lncRNA database interface, the example is for the bc200 lncRNA database entry where you can find useful information about its discovery, expression, function, conservation as well as some citations for the lncRNA.

lncrna db

Home Search Submit Help

bc200

aliases
BCYRN1 (brain cytoplasmic RNA 1)

annotation

Section	Description
Characteristics	200 nucleotide ncRNA (Tiedge (1993)) exapted from an Alu element (Watson (1997)). Transcribed by RNA polymerase III (Mataricelli (1993)).
Expression	Three structural domains, 5' that shares homology with Alu elements, a central A rich region and a 3' unique region (Tiedge (1993)). Expressed predominantly in different regions of the brain. Also shows low level expression in testis but not in other normal tissues examined (Watson (1997), Tiedge (1993), Kurnobov (2001)). RNA sequencing of 11 human tissues confirmed up-regulation of expression in brain (hypothalamus) and low or no expression elsewhere (Cottle (2010)). Deregulated in cancer: expressed in a number of human tumours but not in corresponding normal tissue (Chen (1997)). Link with aging and Alzheimer disease: BC200 expression decreases with aging but is upregulated in Alzheimer's disease (AD). In AD affected brain regions, expression increased with disease severity. RNA localisation showed intense perikaryal staining, showing build up of RNA in cell body (Ma (2007)). Like BC1 ncRNA found in rodents, BC200 also suggested to localise to dendrites (Tiedge (1993)).
Function	Binds several proteins including the signal recognition particle SRP9/14 heterodimer (Gostreva et al (1998)), eukaryotic initiation factor 4A helicase (eIF4A) (Lin (2008)) and Poly(A)-binding protein (PABP), binding to PABP requires the central A rich region (Mukherjee (2002)). Inhibits translation in-vitro and in cultured cells similar to BC1. BC200 binding to eIF4A inhibits it by uncoupling eIF4A ATPase activity from its helicase duplex unwinding activity (Lin (2008)). Translational inhibition also involves BC200 binding to PABP (Gostreva et al (1998)). Anthropoid primates (monkeys, apes and humans) (Shaykin (1998)).
Conservation	BC200, BC1 and G22 likely form a family of independently exapted repetitive elements which have evolved to carry out similar functions in different mammalian species (Mukherjee (2002), Khanam (2007)).
Misc.	More than 200 pseudogenes reported (Kurnobov (2001)).

literature

PUB Med ID	Author	Title	Year
20660672	Cottle	Digital genome-wide ncRNA expression, including SnoRNAs, across 11 human tissues using polyA-neutral amplification.	2010
18316401	Lin	Translational control by a small RNA: dendritic BC1 RNA targets the eukaryotic initiation factor 4A helicase mechanism.	2008
17175535	Khanam	Two primate-specific small non-protein-coding RNAs in transgenic mice: neuronal expression, subcellular localisation and binding partners.	2007

D. LncProfiler qPCR array contents

The qPCR array plate contains assays for 90 lncRNAs and also includes 5 endogenous reference RNAs as normalization signals. Please see the SBI website to download the qPCR array arrangement and $\Delta\Delta CT$ analysis software. www.systembio.com/LncRNA

Human LncProfiler qPCR Array

	1	2	3	4	5	6	7	8	9	10	11	12
A	21A	75K	75L	Alr	AK023948	Alpha 280	Alpha 250	ANRIL	anti-NO524	antiPeg11	BACE1AS	BC200
B	CAR Intergenic	DHFR upstream	Dio3os	DISC2	DLG2AS	E2F4 antisense	EgoA	EGOB	Emx2os	Evf1 and EVF2	GASS	Gomafu
C	H19	H19 antisense	H19 upstream	HAR1A	HAR1B	HOTAIR	HOTAIRM1	HOTTIP	Hoxa11as	HOXA3as	HOXA6as	HULC
D	IGF2AS	IPW	Jpx	Kcnq1ot1	KRAS P1	L1PA16	p21	RoR	SFMBT2	VLDLR	LOC 285194	LUST
E	Malat1	masCRNA	MEG3	MEG9	MER11C	ncR-uPAR	NDM29	NEAT1	Nespas	NRON	NTT	p53 mRNA
F	PCGEM1	PR antisense	PRINS	PSF inhibiting	PTENP1	RINC3	SAF	SCA8	snaR	SNHG1	SNHG3	SNHG4
G	SNHG5	SNHG6	Sox2ot	SRA	ST7OT	TEA ncRNAs	Tmevpg1	TncRNA	Tsix	TUG1	UCA1	UM9-5
H	WT1-AS	Xist	Y RNA-1	Zeb2NAT	Zfas1	Zfx2as	18S rRNA	RNU43	GAPDH	LAMIN A/C	U6	No assay control

Mouse LncProfiler qPCR Array

	1	2	3	4	5	6	7	8	9	10	11	12
A	Adapt33	Alr	AK007836-upstream of	AK141265-Nanog	AK028326-Oct4	AK062072	ATIA	antiPeg11	B2 SINE RNA	BACE1AS	BC1	BGn-As
B	BORG	CDR1-antisense	Dio3os	Dtx1as	Emx2os	Evf2	Foxn2-as	GASS	Gomafu	Gt2-as	H19	H19 antisense
C	m HOTAIR	HOTTIP	Hoxa11as	IGF2AS	Jpx	Kcnq1ot1	linc1242 LINC-Enah	LINC1331	linc1368	linc1612	linc1547	linc1582
D	linc1609-long	linc1609-short	linc1610-long	linc1610-medium	linc1610-short	Linc 1623	Linc1633	lincENC1	lincRNA-Cox2	lincRNA-p21	lincRNA-Sox2	LINC-MD1
E	LXRBSV	Malat1	masCRNA	MEG3	MEG9	MSUR1	Msx1as	Neat1 v1/MEN	Neat1 v2/Men beta	Nespas	Nkx2.2AS	NRON
F	Otx2os	PINC	PINC 1kb isoform	Pldi	Recombinant ion hot spot	RepA transcript	Rlan	Rmst	RNCR3	SCA8 (KLHL1-AS)	Stx3os	Stx3os-clone9
G	SNHG1	SNHG3	SNHG4	SNHG5	SNHG6	Sox2ot	SRA	Tsix	TUG1	Vax2os1	VL30 RNAs	WT1-AS
H	Xist	Y RNAs	Zeb2NAT	Zfas1	Zfx2as	Mistral	18S rRNA	RNU43 (snoRNA)	GAPDH	Beta Actin	U6 snRNA	No assay control

Performing $\Delta\Delta$ Ct data analysis

Test sample plate data				Control sample plate data				Normalized LncRNA Expression levels			
	Well	LncRNA	Average Ct value		Well	LncRNA	Average Ct value	ACT 1	LncRNA	ACT 1	Test expression level compared to Control
1	A1	21A	27.099	27.0947	A1	21A	16.744	3.3657	21A	952.708	852.708
2	A2	73K	27.714	27.7157	A2	73K	33.621	33.6216	73K	230.883	320.883
3	A3	73L	23.82	23.8225	A3	73L	34.411	34.4129	73L	80508.587	80508.587
4	A4	AP	23.557	23.5571	A4	AP	27.548	27.5484	AP	53124.855	53124.855
5	A5	AKO22943	22.442	22.4427	A5	AKO22943	37.711	37.7161	AKO22943	112678.647	112678.647
6	A6	Alpha 210	22.402	22.4021	A6	Alpha 210	37.438	37.4382	Alpha 210	106889.512	106889.512
7	A7	Alpha 258	37.275	37.2794	A7	Alpha 258	35.011	35.0162	Alpha 258	1.320	1.000
8	A8	APRS	29.034	29.0345	A8	APRS	20.452	20.4554	APRS	0.128	0.128
9	A9	ant-RO52A	38.702	38.7024	A9	ant-RO52A	35.302	35.3027	ant-RO52A	0.166	0.166
10	A10	ant-Pog1	38.170	38.1701	A10	ant-Pog1	36.484	36.4852	ant-Pog1	0.984	0.984
11	A11	ant-11A.5 (family)	34.005	34.0051	A11	ant-11A.5 (family)	22.308	22.308	ant-11A.5 (family)	0.001	0.001
12	A12	BC-104	34.724	34.7241	A12	BC-100	22.745	22.7447	BC-100	0.001	0.001

Paste raw Ct
Data for Test
Sample in
Column C.

Paste raw Ct
Data for Control
Sample in
Column H.


Your Data is analyzed automatically
with geometric mean Normalization.
The Fold change levels are in Column M.

E. Endogenous reference controls


Five Subcellular reference controls

LncRNAs can localize and function in the nucleolus, nucleus and in the cytoplasm. The LncProfiler qPCR array includes RNA reference controls to allow for subcellular fractionation studies to identify and profile three separate subcellular compartments.

- **Nucleolus: SnoRNA RNU43 (and some 18S rRNA)**
- **Nucleus: Small Nuclear splicing snRNA U6B**
- **Cytoplasm: GAPDH, Lamin A/C and 18S rRNA**



**Profile your lncRNAs
wherever they are
located within the cell.**




**Example of reference control
amplification plots. Results may
vary depending upon the cell
types analyzed.**

III. Sample Data

A. LncProfiler qPCR Array sample data


Comparison of lncRNA expression in Human normal Fibroblasts and Fibrosarcoma cells

Total RNA was isolated from approximately 2×10^6 cells using standard Trizol/Ethanol precipitation methods. Two micrograms of the RNA samples were used as input into the lncRNA cDNA synthesis kit. The entire cDNA synthesized was then profiled across the 90 lncRNAs and 5 reference controls on the LncProfiler qPCR array. lncRNAs known to be upregulated in cancer, such as Neat1, Anril, Hotair and Malat1 were all detected at highly elevated levels in the fibrosarcoma cells (HT1080) when compared to the normal fibroblast cells (HFF). The data are plotted as normalized expression levels. Selected data are shown in the bar graph.



Profiling Mouse lncRNA expression profiles across various tissues

Two micrograms of the various tissue RNA samples were used as input into the lncRNA cDNA synthesis kit. The entire cDNA synthesized was then profiled across the 90 lncRNAs and 5 reference controls on the Mouse LncProfiler qPCR array. Mouse lncRNAs are differentially expressed across the tissue sets. Representative data are shown in the bar graph.



Protocol per one well of 6-well plate


1. Confluent cells in a 6-well, remove media
2. Add 1ml Trizol directly to cells on plate
3. Incubate at Room temperature for 5 minutes for complete lysis
4. Collect Trizol/cell mixture and transfer to 1.5ml tube
5. Add 200 μ l Chloroform, vortex 15 seconds
6. Centrifuge mixture for 15 minutes at 4°C
7. Collect aqueous layer and transfer to fresh 1.5 ml tube
8. Add equal volume (~250 μ l) Isopropanol, mix by inversion
9. Precipitate RNA overnight at -20°C
10. Centrifuge at 13,000 rpm for 20 minutes
11. Remove supernatant
12. Wash 1X with 500 μ l 80% Ethanol
13. Centrifuge again for 5 minutes at 13,000 rpm
14. Remove supernatant and let air dry 5 minutes
15. Resuspend RNA pellet in 50 μ l water (RNase-free)
16. Use 5 μ l of RNA per cDNA synthesis

LncRNAs are present in serum exosomes


There is high interest in discovering and developing useful RNA-based biofluid markers. The RNAs in patient fluids are present in circulating exosomes. Exosomes are 40 - 100 nm membrane vesicles secreted by most cell types in vivo and in vitro.

Exosomes are found in blood, urine, amniotic fluid, malignant ascite fluids and contain distinct subsets of microRNAs depending upon the tumor from which they are secreted. We wanted to test whether lncRNAs may be present in circulation exosomes as well by using the LncProfiler qPCR array. We precipitated exosomes from a human pooled serum sample (1ml) using SBI's ExoQuick™ exosome precipitation reagent (cat#EXOQ5A-1). The exosome vesicles were then lysed and the exosomal RNA purified using SBI's SeraMir™ kit (cat#RA806A-1). The resulting exoRNA was converted to cDNA using the LncProfiler cDNA synthesis kit. The cDNA was tested in duplicate across all of the LncProfiler qPCR array assay set.

Purifying exosome RNAs and profiling lncRNAs




Discover new biofluid lncRNA biomarkers using the LncProfiler qPCR array kit



B. Specificity Tests


Dissociation analysis

Sample dissociation analyses for reference controls and some lncRNA qPCR assays performed in duplicate.




Clean assay design with no background

Total RNA was prepared from human HT1080 cells in culture. As a control, 2 μ g of this RNA was checked in a mock cDNA synthesis reaction where the reverse transcriptase (RT) was left out. The sample was then tested across the LncRNA Profiler qPCR assays. Separately, we spiked in 10ng of human genomic DNA and tested this sample with the qPCR assays as well. There is ZERO background in the minus RT controls and the LncRNA Profiler assays do not show any amplification signals even with spiked-in genomic DNA. Profile with confidence and only detect lncRNAs.



Gel analysis

Sample gel analyses for selected lncRNA qPCR assay amplicon products from HFF cDNA and separated on a 1.5% agarose gel and stained with ethidium bromide.



IV. Troubleshooting

Problem	Possible Solution
Too much background in qPCR signals	Use much less cDNA in the SYBR Green Mastermix.
No qPCR signals	Did you select SYBR Green as the Detector's Reporter Dye? Did the controls work? Use more cDNA in Mastermix. Check Mastermix contents and try a subset with the controls as a positive control. Also try lowering the Annealing Temperature to 55°C.
How do I select the Threshold level for Ct analysis ?	Typically place the threshold setting in the center of the exponential phase of the amplification curve. Also see the User Manual for your specific instrument or phone their technical support team for guidance.

V. LncRNA Technical References (selected)

Wang KC, Chang HY. **Molecular Mechanisms of Long Noncoding RNAs**. Mol Cell. 2011 Sep 16;43(6):904-14.

Sotillo E, Thomas-Tikhonenko A. **The long reach of noncoding RNAs**. Nat Genet. 2011 Jun 28;43(7):616-7. doi: 10.1038/ng.870.

Guttman, M., J. Donaghey, B.W. Carey, M. Garber, J.K. Grenier, G. Munson, G. Young, A.B. Lucas, et al. 2011. **lincRNAs act in the circuitry controlling pluripotency and differentiation**. Nature. 2011 Aug 28. doi: 10.1038/nature10398.

Cabili, M.N., C. Trapnell, L. Goff, M. Koziol, B. Tazon-Vega, A. Regev, and J.L. Rinn. 2011. **Integrative annotation of human large intergenic noncoding RNAs reveals global properties and specific subclasses**. Genes Dev. 2011 Sep 2.

Cunnington MS, Santibanez Koref M, Mayosi BM, Burn J, Keavney B. **Chromosome 9p21 SNPs Associated with Multiple Disease Phenotypes Correlate with ANRIL Expression**. PLoS Genet. 2010 Apr 8;6(4):e1000899.

Kogo R, Shimamura T, Mimori K, Kawahara K, Imoto S, Sudo T, Tanaka F, Shibata K, Suzuki A, Komune S, Miyano S, Mori M.

Long non-coding RNA HOTAIR regulates Polycomb-dependent chromatin modification and is associated with poor prognosis in colorectal cancers. Cancer Res. 2011 Aug 23.

Bellucci M, Agostini F, Masin M, Tartaglia GG. **Predicting protein associations with long noncoding RNAs**. Nat Methods. 2011 Jun;8(6):444-5.

Perez DS, Hoage TR, Pritchett JR, Ducharme-Smith AL, Halling ML, Ganapathiraju SC, Streng PS, Smith DI. **Long, abundantly expressed non-coding transcripts are altered in cancer**. Hum Mol Genet. 2008 Mar 1;17(5):642-55.

Mus E, Hof PR, Tiedge H. **Dendritic BC200 RNA in aging and in Alzheimer's disease**. Proc Natl Acad Sci U S A. 2007 Jun 19;104(25):10679-84.

Peterlin BM, Brogie JE, Price DH. **7SK snRNA: a noncoding RNA that plays a major role in regulating eukaryotic transcription**. Wiley Interdiscip Rev RNA. 2011 Aug 18. doi: 10.1002/wrna.106.

VIII. Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site:

<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI)
265 North Whisman Rd.
Mountain View, CA 94043

Phone: (650) 968-2200
(888) 266-5066 (Toll Free)

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com

Technical Support: tech@systembio.com

Ordering Information: orders@systembio.com

IX. Licensing and Warranty Statement

Limited Use License

Use of the LncProfilers (*i.e.*, the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Limited Warranty

SBI warrants that the Product meets the specifications described in the accompanying Product Analysis Certificate. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2011 System Biosciences (SBI).