

ExoGlow[™]-Membrane EV Labeling Kit

Cat# EXOGM600A-1

User Manual

Store kit at +4°C

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the License and Warranty Statement contained in this user manual.

Version 2 4/26/2018

Contents

Product Description	2
ist of Components	2
Storage	. 2
General Information	2
Protocol for Labeling	3
Example Data and Applications	4
Fechnical Support	4
icensing and Warranty Statement	
Protocol for Labeling Example Data and Applications Fechnical Support	3 4 4

Product Description

SBI's ExoGlow-Membrane[™] EV Labeling Kit (Cat #EXOGM600A-1) is the latest generation of fluorescent labeling reagent to robustly and specifically label EV membranes. Based on a highly specific membrane sensor, the reagent specifically labels intact EVs for fluorescence-based detection, while showing minimal background in an unbound state. As a result, the kit can be used for most applications that require visualization of labeled EVs for tracking studies. With very low intrinsic background levels, the dye is well suited for demanding EV labeling studies that require high signal-to-noise ratios and specificity.

List of Components

Item	Volume	Storage Temperature
Reaction Buffer	300 µl	4ºC
Labeling Dye	50 µl	4ºC

*The kit is for 25 labeling reactions

Storage

The ExoGlow-Membrane kit is shipped on ice and should be **stored** at +4°C. Properly stored kits are stable for 12 months from the date received.

General Information

1. The reaction size is based on using 50-100 μg of exosomes.

2. The ExoGlow-Membrane excitation and emission values in the bound state are 465nm and 635nm, respectively.

3. We recommend mixing Reaction Buffer well before the usage.

4. Protect labeling dye from light.

5. We recommend to remove free dye before applying labeled exosomes on the cells with PD MiniTrap or PD SpinTrap G-25 buffer exchange column (GE, Cat# 28-9180-07 or 28-9180-04).

6. Due to the nature of the sensor dye and its broad emission spectrum, the signal can be also detected in green channel. Therefore, we don't recommend this product for multiplexing experiments.

Protocol for Labeling:

- 1. Into 12 μl of Reaction Buffer add 2 μl of Labeling Dye and mix well until dye is dissolved completely to make Labeling Reaction buffer.
- 2. Add 50-100 μg exosomes to the Labeling reaction buffer from step 1.
- 3. Mix well the sample and incubate for 30 minutes at RT. The tubes do not need to be rotated during the incubation period.

Protect the tubes from light.

- 4. Remove free unlabeled dye. We recommend buffer exchange column PD MiniTrap or PD SpinTrap G-25 from GE (not provided). A single 1 min spin at the recommended speed will be sufficient for dye removal.
- 5. Alternative protocol for free dye removal. If you have high background after PD SpinTrap G-25, the reprecipitation protocol with ExoQuick-TC is a better choice:
 - Add 35 μL of ExoQuick-TC to 100 μL of Sample and incubate at 4°C from 30 min to overnight.
 - Spin the Eppendorf tube at 10,000 rpm for 10 minutes.
 - Carefully aspirate the supernatant from the corner of the tube.
 - Resuspend the EV pellet in PBS and proceed with downstream applications.
- 6. Labeled exosomes can be added to the cells or used for downstream applications.

Example Data and Applications

FluorescenceBright-fieldMergedImage: DescenceImage: DescenceIm

Figure 1. ExoGlow-Membrane enables clear visualization of labeled EVs being internalized by target cells. We labeled 50 μ g of HEK293T EVs with ExoGlow-Membrane and followed uptake by HEK293T cells. The evenly distributed fluorescence signal shows internalization of labeled EVs by the target cells and the distribution of EV membranes to cellular membranes.

Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site: <u>http://www.systembio.com</u>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI) 2438 Embarcadero Way Palo Alto, CA 94303

Phone:	(650) 968-2200
Toll-Free:	(888) 266-5066
Fax	(650) 968-2277

E-mail:General Information:info@systembio.comTechnical Support:tech@systembio.comOrdering Information:orders@systembio.com

Licensing and Warranty Statement

Limited Use License

Use of the ExoGlow-Membrane EV Labeling Kit (*i.e.*, the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2018 System Biosciences (SBI), All Rights Reserved

This page intentionally left blank.

This page intentionally left blank.

This page intentionally left blank.



System Biosciences (SBI) 2438 Embarcadero Way Palo Alto, CA 94303

 Phone:
 (650) 968-2200

 Toll-Free:
 (888) 266-5066

 Fax:
 (650) 968-2277

E-mail:

General Information:info@systembio.comTechnical Support:tech@systembio.comOrdering Information:orders@systembio.com