



cfMAX™ cfDNA Isolation System

Cat # CFMAX100A-1

User Manual

**Storage: Magnetic Beads, 2-8°C (Do NOT freeze!);
Buffers, room temperature**

Version 2
11/5/2019

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the License and Warranty Statement contained in this user manual.

Contents

Product Description.....	2
List of Components.....	3
Additional Required and Optional Equipment Not Included in Kit	3
Prior to Initial Use.....	4
Protocol	5
Example Data and Applications.....	8
Technical Support.....	11
Licensing and Warranty Statement	11

Product Description

Reliable insights from cfDNA start with high-quality, high-yield isolation

While the existence of cell-free DNA (cfDNA) has been known for almost a century¹, it's only in the past decade or so that researchers have been able to recognize and fully exploit the valuable insights to be gained from analyzing cfDNA. In oncology, cfDNA is promising to reduce the need for invasive tissue biopsies by providing insight into the genetic status of solid tumors through a simple blood draw, and in prenatal testing, cfDNA is a powerful screening tool for genetic abnormalities that is reducing the need for invasive amniocentesis procedures.

cfDNA is released into the bloodstream from normal and cancerous apoptotic or necrotic cells and is highly fragmented, with an average length of around 170 bp². By detecting and analyzing cfDNA, as well as extracellular vesicles (including exosomes) and circulating tumor cells, researchers can obtain rich information to aid in the development of tools for the diagnosis, treatment, prognosis and monitoring of cancer.

Isolating cfDNA from plasma and serum can be challenging, as the amount of cfDNA present in individual samples can vary widely³ and analysis of cfDNA can be complicated by the presence of even low amounts of longer contaminating genomic DNA. With the cfMAX™ cfDNA Isolation System (Cat.# CFMAX100A-1) from SBI, you can overcome these challenges and consistently obtain higher yields of cfDNA than other kits can deliver with a kit optimized for isolation of short DNA fragments. The result is more reliable downstream analysis and greater insights into disease states.

With cfMAX you can:

- Reliably achieve high yields of cfDNA
- Maximize isolation of short DNA fragments while reducing contamination from longer genomic DNA
- Get superior performance compared to competitor kits
- Maximize productivity with the quick and easy magnetic bead isolation workflow

The cfMAX cfDNA Isolation System maximizes your flexibility and accommodates your lab's capabilities:

- Suitable for a range of input volumes
- Compatible with both manual and automated isolation workflows
- Interchangeable with MagMAX™ cfDNA isolation kit in Kingfisher systems

A faster, easier workflow than competitor kits

The cfMAX cfDNA Isolation System uses a fast and simple magnetic bead isolation workflow—lyse your serum or plasma sample with the included cfMAX buffer, add beads, immobilize beads with the magnet, wash the beads, and elute cfDNA.




Figure 1. The fast cfMAX workflow.

References

1. Mandel P, Metais P. Les acides nucleiques du plasma sanguin chez l'homme [in French]. *C R Seances Soc Biol Fil.* 1948;**142**:241–243.
2. Johann DJ Jr, et al. Liquid biopsy and its role in an advanced clinical trial for lung cancer. *Exp Biol Med (Maywood)*. 2018 Feb;**243**(3):262–271. PMID: PMC5813874.
3. Babji D, Nayak R, Bhat K, and Kotrashetti V. Cell-free tumor DNA: Emerging reality in oral squamous cell carcinoma. *J Oral Maxillofac Pathol.* 2019 May-Aug; **23**(2): 273–279. PMID: PMC6714275.

List of Components

Table 1. cfMAX cfDNA Isolation System Components*		
Component	Qty/Volume	Storage Temperature
Magnetic Beads	2 ml	2-8°C (Do NOT freeze!)
Lysis/Binding Buffer	85 ml	RT
Wash Buffer	40 ml	
Elution Buffer	3 ml	

*Each kit contains enough reagents to process 80 mL of serum or plasma. To purchase bulk volumes, please inquire at services@systembio.com.

Additional Required and Optional Equipment Not Included in Kit

- Pipettes
- Fisher Scientific Vortex Mixer or similar vortexing mixer
- Fisherbrand™ Vortex Mixer accessories: Insert Retainer (cat# 02-215-391)
- Foam inserts set (cat# 02-215-394 or 02-215-386)
- Magnet stand for molecular applications
- 1.5 ml non-stick Eppendorf tube(s)
- Fresh absolute ethanol

Prior to Initial Use

The Lysis/Binding and Wash Buffer are shipped as a concentrate. If precipitate is present in either solution, incubate solution at 37°C for 30 minutes or more. Add absolute ethanol as indicated into the bottle of Lysis/Binding Buffer and Wash Buffer. Mix well and tightly enclose the bottle cap to prevent ethanol evaporation.

- Add 15 ml of fresh absolute ethanol to Lysis/Binding Buffer and mix by inverting gently
- Add 40 ml of fresh absolute ethanol to Wash Buffer and mix by inverting gently
- Prepare fresh 80% ethanol solution prior to each extraction

Starting Material: Both fresh and frozen plasma can be used with the Cell-Free DNA isolation protocol. Fresh plasma, however, tends to have higher yields.

Quantification: We suggest using Bioanalyzer High Sensitivity DNA assay for measuring cfDNA.

Streck Cell-Free DNA BCT Tube(s): Plasma from blood collected with Streck Cell-Free DNA BCT Tube(s) must go through a Proteinase K treatment prior to Cell-Free DNA isolation to ensure optimal yields.

Forgoing Proteinase K treatment may decrease yields by 50%. For EDTA, Lbgard tube, Proteinase K treatment can be skipped.

Protocol

Proteinase K Pre-treatment (for plasma sample collected using Streck Cell-free DNA BCT tubes)

Note: If blood was not collected with Streck Cell-Free DNA BCT tube(s) skip this step.

If samples were collected using a Streck Cell-Free DNA BCT tube(s), Proteinase K treatment is required to ensure optimal yields.

Plasma (ml)	Proteinase K	20% SDS Solution
2	30 ul	100 ul
5	75 ul	250 ul
7	105 ul	350 ul
X	0.015X ml	0.050X ml

1. Add the appropriate amount of plasma to an appropriately sized tube(s): For example, use 1.5 ml Eppendorf tubes for 0.5ml to 1ml plasma, 15 ml conical tubes for plasma 2 to 5ml, and 50ml conical tubes for plasma volume larger than 5 ml.

Note: Using a 50 ml tube(s) for 5 ml or more of plasma is recommended over a 15 ml tube(s). While a 15 ml tube(s) will work it may lead to slightly lower yields.

2. Add 15 μ l of Proteinase K (20 mg/ml) for every 1 ml of plasma used.
3. Add 50 μ l of 20% SDS solution for every 1 ml of plasma used.
4. Mix by inverting gently 5 times.
5. Incubate at 60°C for 20 minutes.
6. After incubation, place tube(s) on ice for 5 minutes to cool tube(s) to room temperature.
7. Once tube(s) are at room temperature, proceed to Step 1 of the Lysis/Binding section in Protocol.

cfDNA isolation

Before starting the protocol, determine the amount of plasma/serum to be used for extraction and calculate the amounts of buffers and beads needed. Any amount from 500 μ l to 8 ml of plasma/serum can be used. Scale buffer and bead volumes accordingly using the table below.

Small (0.2 ml to 7.9 ml) Sample Protocol

Plasma	Lysis/Binding Buffer	Beads solution	Tube size
X (x=ml of plasma)	1.25x	0.025x	Choose by plasma volume
5ml	6.25ml	125ul	15ml or 50ml *
7ml	8.75ml	175ul	50ml

*Using a 50 ml tube(s) for 5 ml or more of plasma is recommended over a 15 ml tube(s). While a 15 ml tube(s) will work it may lead to slightly lower yields.

Lysis/ Binding

8. Add the appropriate amount of plasma/serum to appropriately sized tube(s).
9. Add 1.25 ml of **Lysis/Binding Buffer** for every 1 ml of plasma/serum used. Briefly mix.
10. Add 25 μ l of **Magnetic Beads** for every 1 ml of plasma/serum.

Important: Mix beads well prior to aliquoting. There should be no visible sedimentation at the bottom of the solution after mixing. Beads will settle quickly, so be sure to mix the magnetic bead solution after adding to each sample. Failure to do so may result in inconsistent yields.

11. Vortex or shake tube(s) vigorously for 10 minutes at room temperature to bind the cfDNA to the beads.

To obtain high yields, ensure that sample/buffer solution is mixed vigorously in the tube(s). A vortexing mixer (speed 7 or higher) with a tube-holder is highly recommended.

12. Place tube(s) into a magnet stand for 2 to 5 minutes, or until solution clears.
13. While keeping the tube(s) on the magnet stand, remove supernatant. Be careful not to remove magnetic particles.
14. Keep tube(s) on magnet stand for 1 more minute, and carefully remove any residual liquid.

First Wash

15. Add 1 ml of **Wash Buffer** to the lysis/binding tube(s).
16. Resuspend beads by vortexing for 20 seconds or pipetting up and down 15 times.
17. Transfer magnetic particle suspension into a new 1.5 ml micro tube(s) on magnet stand.
18. Allow beads to attach to magnet stand for 30 seconds.
19. Pipette supernatant from the 1.5 ml tube(s) and use the supernatant to wash the lysis/binding tube(s).
20. Transfer the rest of the magnetic particles in lysis/binding tube(s) to the 1.5 ml tube(s).
21. Keep tube(s) on magnet stand for 30 seconds or until solution is clear.
22. Remove as much buffer as possible using a 1 ml pipette.
23. Tap magnet stand on bench 5 times and remove remaining wash buffer with a 200 μ l pipette.
24. Transfer tube(s) to non-magnetic rack and add 1 ml of **Wash Buffer**.
25. Resuspend beads by vortexing for 20 seconds or pipetting up and down 15 times.

26. Centrifuge tube(s) briefly.

*Centrifuge steps are needed when using vortexing to resuspend beads. Only a brief spin is needed to remove solution from tube(s) lid.

27. Place tube(s) on magnet stand for 30 seconds or until solution is clear.

28. Remove as much buffer as possible using a 1 ml pipette.

29. Tap magnet stand on bench 5 times and remove remaining wash buffer with a 200 μ l pipette.

Second Wash

30. Transfer tube(s) to non-magnetic rack and add 1 ml of **80% EtOH**.

31. Resuspend beads by vortexing for 20 seconds or pipetting up and down 15 times.

32. Centrifuge tube(s) briefly.

33. Place on magnet stand for 30 seconds or until solution clears.

34. Remove as much buffer as possible using a 1 ml pipette.

35. Tap magnet stand on bench 5 times and remove remaining EtOH with a 200 μ l pipette.

36. Transfer tube(s) to non-magnetic rack and add 1 ml of **80% EtOH**.

37. Resuspend beads by vortexing for 20 seconds or pipetting up and down 15 times.

38. Centrifuge tube(s) briefly.

39. Place on magnet stand for 20 seconds or until solution clears.

40. Remove as much EtOH as possible using a 1 ml pipette and leave cap open.

41. Tap magnet stand with tube(s) on bench 5 times and remove remaining EtOH with a 200 μ l pipette.

42. Centrifuge tube(s) briefly and place on magnet stand for 10 seconds.

43. Remove small trace of remaining EtOH with a 20 μ l pipette.

44. Leave tube(s) open on magnet stand and allow magnetic particles to air dry for 4-5 minutes.

*Be careful not to over-dry or beads may stick to tube(s).

Elution Step

45. Transfer the microtube(s) to non-magnetic rack and add desired volume of **Elution Buffer** and resuspend beads.

Important: A minimum of 12.5 μ l of **Elution Buffer** per ml of plasma is recommended to elute DNA to ensure optimal yields. Typically, 20- 30 μ l of Elution Buffer is used per ml of plasma.

46. Vortex or shake tube(s) vigorously for 5 minutes.
47. Centrifuge tube(s) briefly.
48. Place tube(s) on magnetic rack for 10 to 30 seconds.
49. Transfer elute into a new 1.5 ml tube(s).

Example Data and Applications

The cfMAX cfDNA Isolation System extracts more cfDNA than other kits




Figure 2. The cfMAX cfDNA Isolation System extracts more cfDNA than other kits. We added known amounts of a synthetic cfDNA to DNA-depleted plasma replenished with 5 ng of sheared genomic DNA and compared cfDNA isolation efficiency of the cfMAX cfDNA Isolation System to two competitor kits. The cfMAX cfDNA kit more consistently and linearly delivered the expected amount of the spiked-in DNA than the competitor kits (n=3 for each data point).

cfMAX delivers cfDNA with less contaminating genomic DNA and higher yields of cfDNA than competing kits




Figure 3. cfMAX delivers cfDNA with less contaminating genomic DNA and higher yields of cfDNA than competing kits. We isolated cfDNA from 2 mL of plasma using the cfMAX System and two competitor kits. The cfMAX system shows less contaminating genomic DNA than the other kits.

cfMAX delivers cfDNA with less protein carryover than competing kits




Figure 4. cfMAX delivers cfDNA with less protein carryover than competing kits. We isolated cfDNA from 2 mL of serum using the cfMAX System and two competitor kits. The cfMAX system shows significantly less protein carryover than Company Q’s kit, and similar low amounts of protein carryover as Company T’s kit, as determined by Qubit Protein Assay.

cfDNA isolation with cfMAX is scalable




Figure 5. cfDNA isolation with cfMAX is scalable. We isolated cfDNA from increasing volumes of the same serum sample using cfMAX. Isolation is highly linear across the input sample volumes.

cfDNA isolation with cfMAX is reproducible




Figure 6. cfDNA isolation with cfMAX is reproducible. We isolated cfDNA from 2 mL of the same serum sample in duplicate using cfMAX. The similar yields demonstrate the reproducibility of cfMAX.

cfMAX is optimized for isolating smaller DNA fragments




Figure 7. cfMAX is optimized for isolating smaller DNA fragments. Using 105- and 236 bp fragments of GAPDH, we assessed the efficiency of DNA isolation by cfMAX and the other kits from serum. The DNA isolated by the cfMAX kit shows a much larger difference in Ct values than the DNA isolated from the other kits, indicating that the cfMAX kit more preferentially isolates smaller DNA fragments.

Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site: <http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200
Toll-Free: (888) 266-5066
Fax (650) 968-2277

E-mail:

General Information: info@systembio.com
Technical Support: tech@systembio.com
Ordering Information: orders@systembio.com

Licensing and Warranty Statement

Limited Use License

Use of the cfMAX cfDNA Isolation System (*i.e.*, the “Product”) is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

This Product is covered by licensed Patents. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This

limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2019 System Biosciences (SBI), All Rights Reserved

This page intentionally left blank.

This page intentionally left blank.

This page intentionally left blank.



System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200
Toll-Free: (888) 266-5066
Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com
Technical Support: tech@systembio.com
Ordering Information: orders@systembio.com