



MicroRNA Research

Expression Profiling • Overexpression • Knockdown • Libraries • Target Sensors



Spend less time making tools, and more time making discoveries

Essential Tools for MicroRNA and LncRNA Research


MicroRNAs are a class of naturally-occurring, small non-coding RNAs that control gene expression by translational repression or mRNA degradation. They are widely expressed throughout the plant and animal kingdoms and may comprise 1-5% of animal genes. Like protein-coding genes, microRNAs are transcribed in the nucleus as long primary transcripts (pri-microRNAs). However, distinct from protein-coding genes, they are subsequently cleaved by the nuclear RNase III enzyme Drosha to produce structured, stem-loop precursor molecules (pre-microRNAs) of 70-100 nucleotides (nt) in length. The pre-microRNAs are then exported to the cytoplasm in a process involving exportin-5, where the RNase III enzyme Dicer further processes them into mature microRNAs (~22 nt). One strand of the microRNA duplex is subsequently incorporated into the RNA-induced silencing complex (RISC) that regulates target gene expression. SBI has integrated solutions for microRNA analysis from expression profiling with QuantiMir™ & RNA-Quant™ kits, to overexpression of microRNAs using constitutive LentiMir™ and inducible iMIR™ precursor clones, as well as miRZip™ microRNA knockdown constructs and the microRNA target identification system miR-Selection™. SBI has also developed new technologies to quantify long noncoding RNAs (LncRNAs) by qPCR as well as track their cellular movement using the Spinach RNA fluorophore tag technology.

Contact us today

Telephone: 650-968-2200

Email: tech@systembio.com

Web: www.systembio.com



Index

| | |
|---|---|
| » MicroRNA Profiling | 2 |
| » LncRNA Profiling | 3 |
| » Overexpress MicroRNAs | 4 |
| » Knockdown MicroRNAs | 5 |
| » MicroRNA Libraries and Target Sensor System | 6 |




Accurate and Sensitive MicroRNA Profiling

Measure microRNAs by qPCR with QuantiMir™

MicroRNA and siRNA expression analysis is made easy with the QuantiMir™ small RNA quantitation system. Generate qPCR-ready cDNA from total RNA for accurate and sensitive expression measurements. QuantiMir rapidly tags and converts all small RNAs into detectable cDNAs for qPCR—one cDNA synthesis required to quantitate any microRNA.

Capture more expression data with genome-wide miRNome qPCR Profilers


Complete Human, Mouse and Rat qPCR Arrays - 100% miRBase compatible. Characterize microRNA signatures in stem cells, cancer cell lines and FFPE samples. Custom qPCR arrays are also available.



Human complete set


Mouse complete set

Rat complete set



All miR major and minor miR* forms

miRNome Array Controls



Cancer and Stem Cell Focused qPCR Arrays

OncoMir microRNA qPCR Array

Detect microRNAs most commonly found in Cancer

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|--------------|------------|-----------|-----------|----------|----------|----------|-----------|------------|------------|---------|----------|
| A | let-7 family | miR-7 | miR-92 | miR-93 | miR-9-1 | miR-101 | miR-103 | miR-106a | miR-106b | miR-107 | miR-10b | miR-1 |
| B | miR-122a | miR-125a | miR-125b | miR-126 | miR-128b | miR-132 | miR-133a | miR-134 | miR-135b | miR-136 | miR-137 | miR-140 |
| C | miR-141 | miR-142-3p | miR-143 | miR-145 | miR-146a | miR-149 | miR-150 | miR-151 | miR-153 | miR-154 | miR-155 | miR-15a |
| D | miR-15b | miR-16 | miR-17-3p | miR-17-5p | miR-181a | miR-181b | miR-181c | miR-181d | miR-183 | miR-185 | miR-186 | miR-188 |
| E | miR-18a | miR-190 | miR-191 | miR-192 | miR-194 | miR-195 | miR-196a | miR-197 | miR-198 | miR-199ab | miR-30b | miR-19ab |
| F | miR-95 | miR-20a | miR-200a | miR-200b | miR-200c | miR-202 | miR-203 | miR-204 | miR-205 | miR-206 | miR-21 | miR-210 |
| G | miR-214 | miR-215 | miR-372 | miR-373 | miR-218 | miR-219 | miR-22 | miR-488 | miR-221 | miR-222 | miR-223 | miR-224 |
| H | miR-23a | miR-24 | miR-25 | miR-26a | miR-26b | miR-27ab | miR-30c | miR-29abc | miR-30a-3p | miR-30a-5p | miR-296 | U6 snRNA |

Stem Cell microRNA qPCR Array

Profile microRNAs in Pluripotency and Differentiation

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|------------|----------|----------|----------|------------|------------|----------|----------|----------|----------|-----------|------------|
| A | miR-18a | miR-19a | miR-19b | miR-24 | miR-25 | miR-17-5p | miR-30c | miR-34a | miR-106a | miR-106b | miR-130b | miR-141 |
| B | miR-150 | miR-199a | miR-200b | miR-200c | miR-301 | miR-302a | miR-302b | miR-302c | miR-302d | miR-367 | miR-368 | miR-369-5p |
| C | miR-369-3p | miR-370 | miR-371 | miR-372 | miR-373 | let-7a | let-7b | miR-10a | miR-10b | miR-16 | miR-17-3p | miR-20a |
| D | miR-20b | miR-23a | miR-23b | miR-26a | miR-26b | miR-30b | miR-30d | miR-32 | miR-33 | miR-92 | miR-93 | miR-99a |
| E | miR-101 | miR-107 | miR-126 | miR-130a | miR-142-5p | miR-142-3p | miR-146a | miR-146b | miR-155 | miR-181a | miR-181b | miR-181c |
| F | miR-181d | miR-191 | miR-193a | miR-193b | miR-197 | miR-221 | miR-223 | miR-339 | miR-9 | miR-103 | miR-124a | miR-125a |
| G | miR-125b | miR-127 | miR-128a | miR-128b | miR-132 | miR-134 | miR-135a | miR-135b | miR-136 | miR-138 | miR-149 | miR-153 |
| H | miR-154 | miR-183 | miR-218 | miR-219 | miR-222 | miR-1 | miR-122a | miR-133a | miR-133b | miR-195 | miR-206 | U6 snRNA |

For more information: www.systembio.com/mirnome




Quantitate Long Non-coding RNAs by qPCR


Long non-coding RNAs (lncRNAs) and large intergenic non-coding RNAs (lincRNAs) are emerging as master regulators of embryonic pluripotency, differentiation, patterning of the body axis and promoting developmental transitions. lncRNAs are larger than 200 nucleotides in length and are pervasively expressed across the genome. lncRNAs maintain the commitment to specific cellular fates through modification and remodeling of chromatin at the epigenetic level.

Dysregulated expression of lncRNAs has been shown to be associated with a broad range of diseases such as Alzheimer's, psoriasis and many cancers. Studying the expression patterns of lncRNAs will be a crucial method to understanding the roles they play in many model systems. SBI has built a sensitive, accurate and robust qPCR array that is strand-specific to enable researchers to closely profile the expression changes in the top lncRNAs known to date. All of the lncRNAs on the qPCR array have validated primer sets for well-annotated lncRNAs that are registered in the lncRNA database created by Dr. John Mattick (www.lncrnadb.org).

How the lncRNA cDNA is made



Profile lncRNAs in Skin Cancer



Human Stem Cell and Cancer LncProfiler qPCR Array

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|----------------|---------------|--------------|----------------|----------|----------------|-----------|--------|------------|---------------|------------|------------------|
| A | 21A | 7SK | 7SL | Air | AK023948 | Alpha 280 | Alpha 250 | ANRIL | anti-NOS2A | antiPeg11 | BACE1AS | BC200 |
| B | CAR intergenic | DHFR upstream | Dio3os | DISC2 | DLG2AS | E2F4 antisense | EgoA | EGOB | Emx2os | Evf1 and EVF2 | GASS | Gomafu |
| C | H19 | H19 antisense | H19 upstream | HAR1A | HAR1B | HOTAIR | HOTAIRM1 | HOTTIP | Hoxa11as | HOXA3as | HOXA6as | HULC |
| D | IGF2AS | IPW | Jpx | Kcnq1ot1 | KRASP1 | L1PA16 | p21 | RoR | SFMBT2 | VLDLR | LOC 285194 | LUST |
| E | Malat1 | masRNA | MEG3 | MEG9 | MER11C | ncR-uPAR | NDM29 | NEAT1 | Nespa | NRON | NTT | p53 mRNA |
| F | PCGEM1 | PR antisense | PRINS | PSF inhibiting | PTENP1 | RNCR3 | SAF | SCA8 | snaR | SNHG1 | SNHG3 | SNHG4 |
| G | SNHG5 | SNHG6 | Sox2ot | SRA | ST7OT | TEA ncRNAs | Tmevpg1 | TncRNA | Tsix | TUG1 | UCA1 | UM9-5 |
| H | WT1-AS | Xist | YRNA-1 | Zeb2NAT | Zfas1 | Zfx2as | 18S rRNA | RNU43 | GAPDH | LAMIN A/C | U6 | No assay control |


Mouse LncProfiler qPCR Array

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---------------|----------------|---------------------------|----------------------------|-------------------------|-----------------|-------------------|---------------|-------------------|-----------------|--------------|------------------|
| A | Adapt33 | Air | AK007836-upstream of Oct4 | AK141205-upstream of Nanog | AK028326-Oct4 | AK082072 | ATIA | antiPeg11 | B2 SINE RNA | BACE1AS | BC1 | B6n-As |
| B | BORG | CDR1-antisense | Dio3os | Dlx1as | Emx2os | Erf2 | Foxn2-as | GASS | Gomafu | Gtl2-as | H19 | H19 antisense |
| C | m HOTAIR | HOTTIP | Hoxa11as | IGF2AS | Jpx | Kcnq1ot1 | linc1242 LINC-Enh | LINC1331 | linc1368 | Linc1612 | linc1547 | linc1582 |
| D | linc1609-long | linc1609-short | linc1610-long | linc1610-medium | linc1610-short | Linc 1623 | Linc1633 | lincENC1 | lincRNA-Cox2 | lincRNA-p21 | lincRNA-Sox2 | LINC -MD1 |
| E | LXRBSV | Malat1 | masRNA | MEG3 | MEG9 | MSUR1 | Mx1as | Neat1 v1/MEN | Neat1 v2/Men beta | Nespa | Nkx2.2AS | NRON |
| F | Otx2os | PINC | PINC 1kb isoform | Pldi | Recombinat ion hot spot | RepA transcript | Rlan | Rmst | Rmst | SCA8 (KLHL1-AS) | Six3os | Six3os-clone9 |
| G | SNHG1 | SNHG3 | SNHG4 | SNHG5 | SNHG6 | Sox2ot | SRA | Tsix | TUG1 | Vax2os1 | VL30 RNAs | WT1-AS |
| H | Xist | Y RNAs | Zeb2NAT | Zfas1 | Zfx2as | Mistral | 18S rRNA | RNU43 (snRNA) | GAPDH | Beta Actin | U6 snRNA | No assay control |

Human Disease-related qPCR Array

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|------------|-----------|-----------|-----------|---------|---------------|----------------|----------|----------|---------|--------|--------------------|
| A | 21A | AAA1 | aHIF | AK023948 | ANRIL | anti-NOS2A | BACE1AS | BC017743 | BC043430 | BC200 | BCMS | BIC |
| B | CCND1 ANCR | CMPD | DD3 | DGCR5 | DISC2 | DLG2AS | EGO | GAS5 | GOMAFU | H19 | H19-AS | HAR1A |
| C | HAR1B | HOTAIR | HOTAIRM1 | HOTTIP | HOXA1AS | HOXA3AS | HOXA6AS | HOXA11AS | HULC | IPW | IGF2AS | |
| D | KRASP1 | L1PA16 | LIT | LOC285194 | LUST | LincRNA VLDLR | LincRNA SFMBT2 | MALAT1 | MEG3 | MER11C | NEAT1 | NCRMS |
| E | NDM29 | PANDA | PARS | PCAT1 | PCAT14 | PCAT29 | PCAT32 | PCAT43 | PCGEM1 | PRAT2 | PRINS | PSF inhibiting RNA |
| F | PTENP1 | RMRP | ROR | SAF | SCA8 | Sox2OT | SRA | ST7OT1 | ST7OT2 | ST7OT3 | ST7OT4 | Telomerase RNA |
| G | TMEVPG1 | TU_001769 | TUG1 | UCA1 | WT1-AS | Y1 | Y3 | Y4 | Y5 | ZEB2NAT | 7SK | Negative control |
| H | 7SLncRNA | 5.8S rRNA | U6 scaRNA | U6 snRNA | ACTB | B2M | PGK1 | GAPDH | HPRT1 | RPL1A | RPL13A | GDC |


Human Cancer and Stem Cell qPCR Assay Specificity Tests



Overexpress Human and Mouse MicroRNAs

There are expected to be over 2,000 different microRNAs encoded in the human and mouse genomes and they function by either blocking translation of, or degrading, mRNA species corresponding to specific genes. While the number of verified human and mouse microRNAs are expanding, there is an increasing need for effective functional testing. SBI offers an extensive collection of microRNA precursors in lentiviral vectors that can be used to modulate the expression of their targeted mRNAs to study microRNA function. The lentiviral vector constructs can be packaged into pseudoviral particles and delivered to primary cells, stem cells, or other hard-to-transfect cell lines and can be used *in vivo*. SBI can package any construct into high titer lentivirus as a custom service for your studies.

Each construct in SBI's collection consists of the native stem loop structure and 200-400 base pairs of upstream and downstream flanking genomic sequence. This unique feature ensures that the microRNAs expressed from SBI's constructs will be correctly processed in the cell into mature microRNA.




Constitutive Human and Mouse Lenti-miRs

The advantage of SBI's constructs is that they can be used in transient transfection experiments as well as stably expressed in a wide variety of cell types, as opposed to synthetic microRNAs that can only be transiently expressed in cells.

MCF-7 cells transduced with Lenti-miR-205



Phase Contrast

GFP Fluorescence




Inducible Human iMIRS

The dual promoter vector configuration enables the inducible expression of the microRNA precursor. The downstream EF1 α promoter will constitutively express the GFP+Puro dual markers independent of the cumate switch activation. This system is expressing the microRNA by default and you choose whether or not to add the CymR repressor into the system for inducible regulation.



Mature miRNA Expression



MCF-7 Cells infected (MOI 5:1) with Lentivirus made from PMIRH205PA-1 (pre-miR-205) construct. Images taken after 48 hours. Expression of mature miR-205 was validated and quantitated using ABI TaqMan assays.

For more information: www.systembio.com/lentimir


Block MicroRNA Functions with miRZips™

miRZip anti-sense microRNAs are stably expressed RNAi hairpins that have anti-microRNA activity. These miRZip shRNAs produce short, single-stranded anti-microRNAs that competitively bind their endogenous microRNA target and inhibit its function. The result is the derepression and elevation of the protein levels of the transcripts targeted by the microRNA being "zipped".


Permanent MicroRNA Knockdown

The miRZip short hairpin RNAs are cloned into SBI's pGreenPuro™ shRNA expression lentivector. The miRZip hairpins are designed for asymmetry such that the lower strand is preferred for producing anti-sense microRNAs that are fully complementary to a specific microRNA target, forming a stable duplex that keeps the zipped microRNA from binding target UTR sequences. Custom constructs, including multiple miRZips in a single vector, are also available.

How miRZips Work



Discover MicroRNA Phenotypes



Inhibition of miR-145 results in cMyc protein upregulation



For more information: www.systembio.com/mirzips



MicroRNA Lentivirus Screening Libraries

The Lenti-miR™, OncoMir and miRZip™ pooled virus libraries are tools that enable the study of phenotypic effects associated with the overexpression or suppression of individual microRNAs. The lentivirus preparation is pseudotyped with VSV-G that allows for broad cellular tropism. Use these libraries with hard-to-transfect mammalian cell lines, primary cells, non-dividing cells and even whole animal studies. Transduced cells exhibiting the phenotypes of interest are isolated by selection or sorting. The microRNA or microRNAs responsible for generating the phenotypes of interest may be recovered through simple genomic PCR using lentivector-specific primers followed by direct sequencing of microRNA precursor or miRZip anti-miR shRNA hairpin clones.

Perform one transduction and easily identify the microRNAs involved in your phenotypic screen

SBI's pooled lentiviral libraries allow you to perform high-throughput screening studies on a genome-wide or cancer-focused basis. Pooled lentiviral libraries enable simultaneous identification of multiple genes that alter a specific cellular phenotype in a single experiment. Lentiviral libraries are available as prepackaged virus, so you can begin transducing cells the day you receive the library.


Use microRNA libraries to find targets

Target Selection System

Identify Binding Sites through Cellular Selection and Quantitate Target Interactions

Novel cellular method to detect microRNA binding to its target mRNA using a dual reporter system featuring Luciferase (Fire) and a Cytotoxic Sensor (Ctx). The miR-Select platform captures the 3' UTR to microRNA binding event using a survival screen by modulating the reduction of the cytotoxic sensor. Validation is made simple using the built-in Luciferase reporter. This powerful and elegant technology finally enables the accurate identification of microRNA targets.


The miR-Select Fire-Ctx-UTR Lentivector



How to use Pooled Lentiviral Libraries

Overview of Library Screening Protocol

Transduction with lentivirus library



Infect Target Cells

Treatment to induce phenotype

Treat/Select Cells



Selection for phenotype

Recover Cells with Specific Phenotype

Amplification of effector sequences from selected cells



Simple Genomic PCR to Identify microRNA Effectors



Directly Sequence PCR Amplicons

MicroRNAs Affecting Phenotype Identified

For more information: www.systembio.com/microrna



265 North Whisman Rd.
Mountain View, CA 94043
Telephone: 650-968-2200
Email: info@systembio.com
www.systembio.com

