



---

# **PrecisionX™ Multiplex gRNA Cloning Kit**

**Cat. # CAS9-GRNA-KIT**

## ***User Manual***

---

**Store at -20°C upon receipt**

**A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.**

# Contents

- I. Introduction .....2
  - A. Key Features .....3
  - B. Overview of the Multiplex gRNA Cloning Kit .....4
  - C. List of Components .....7
  - D. Storage .....7
  - E. Other Reagents Needed .....7
- II. Protocol .....8
  - A. Design of Primers for PCR reaction (For Dual gRNA cloning)  
.....8
  - B. Design of Primers for PCR reaction (For “Quad-Plex” gRNA  
cloning) .....11
  - C. Setting up the PCR Reaction .....17
  - D. Fusion Reaction .....19
  - E. Transformation of Cells .....21
  - F. Confirmation of Positive Clones .....21
- III. Validation Data .....24
  - A. Gene repair by double-nicking and HR .....24
  - B. Gene Knockout by dual gRNAs targeting RFP .....25
  - C. Tandem paired-nicking (4 gRNAs in nickase vector) .....26
- IV. Appendix .....27
  - A. Primer design for dual gRNA/single promoter constructs  
27
  - B. Sequencing primer information .....29
- V. Troubleshooting .....30

VI.      References.....32

VII.     Technical Support .....32

VIII.    Licensing and Warranty .....33

I.      Introduction

The recent discovery of the CRISPR/Cas9 system has provided researchers an invaluable tool to target and modify any genomic sequence with high levels of efficacy and specificity. The system, consisting of a nuclease (Cas9) and a DNA-directed guide RNA (gRNA), allows for sequence-specific cleavage of target sequences containing a protospacer adaptor motif “NGG”. By changing the gRNA target sequence, virtually any gene sequence upstream of a PAM motif can be targeted by the CRISPR/Cas9 system, enabling the possibility of systematic targeting of sequences on a genomic scale.

In order to facilitate the cloning of gRNAs into Cas9 expression vectors, SBI has developed a revolutionary cloning kit, the PrecisionX™ Multiplex gRNA Cloning Kit (Cat# CAS9-GRNA-KIT). This system allows for the cloning of **multiple gRNAs** into any Cas9/gRNA “all-in-one” expression vector or gRNA cloning vector of the customer’s choice, including SBI’s PrecisionX Cas9/gRNA SmartNuclease plasmids and lentivectors. The Multiplex gRNA Cloning Kit is also compatible with the most popular Cas9/gRNA cloning vectors developed in labs across the world, such as pX330, pX335, pX458, and pX459. The broad applicability of this kit across different formats of Cas9/gRNA vectors will appeal to both beginners and advanced users of CRISPR/Cas9 technology alike by simplifying the process and lead time required to generate powerful CRISPR/Cas9 constructs for their experiments.

## A. Key Features

The PrecisionX Multiplex gRNA Cloning Kit provides the following advantages to the researcher:

- **Saves time and reagents - No need to clone separate constructs bearing different gRNAs**
- **Simple two-step PCR and “fusion” reaction to generate multi-cistronic constructs**
- **Compatible with all types of Cas9/gRNA expression vectors containing H1 or U6 promoters**
- **Ideal for Cas9 nickase applications requiring expression of two gRNAs simultaneously for precise targeting**
- **Single plasmid system reduces transfection variability associated with multi-plasmid transfection**
- **Facilitates the use of more advanced Cas9 targeting applications, including tandem gRNA cutting and tandem double-nicking**
- **Enables precise deletion of defined genomic segments with a single vector**

The Multiplex gRNA Cloning Kit is based on SBI's proven Cold Fusion cloning technology, which has been the standard for restriction and ligase-free cloning technology for many years, freeing customers from laborious and time-consuming traditional cloning approaches.

## B. Overview of the Multiplex gRNA Cloning Kit

### Addition of a gRNA cassette (H1-gRNA) into U6 vector




Fig. 1: General workflow of Multiplex gRNA Cloning Kit to insert H1-gRNA cassette into existing vector with U6 promoter

**Addition of a gRNA cassette (U6-gRNA) into H1 vector**




Fig. 2: General workflow of Multiplex gRNA Cloning Kit to insert U6-gRNA cassette into existing vector with H1 promoter

The multiplex cloning reaction involves two steps:

**Step 1:**

Overlapping PCR of the two primers (designed by the user) with a scaffold-promoter block (provided in the kit) to create a PCR product containing the desired amplicon containing both gRNAs.

**Step 2:**

Using the PCR amplicon created in Step 1, fuse the fragment into linearized destination vector.

This cloning process is extremely efficient, and has been demonstrated to be useful in applications such as targeting non-coding RNAs (Ho *et al.* 2014). In addition, the kit can be adapted to more than two gRNAs in single reaction. SBI has validated the cloning and fusion of 4 separate guide RNAs (Fig. 3 below) in a two-step PCR and fusion reaction into a linearized destination vector (Section II.B). The level of flexibility, ease of use, and efficiency of the Multiplex gRNA Cloning Kit is an industry-first – designed with the researcher in mind who needs freedom from limitations imposed by traditional gRNA cloning methods.




Fig. 3: Example diagram of generating a “quad-plex” multi-cistronic gRNA cassette using the Multiplex gRNA Cloning Kit

## C. List of Components

### Cat #CAS9-GRNA-KIT (10 reactions)

Component	Volume
5X Master mix	20ul
Linearized vector (positive control)	5ul
PCR product (positive control)	5ul
H1 Block	5ul
U6 Block	5ul

**Note: Primers for the PCR reaction are sequence-dependent and will need to be designed by the end-user per guidelines in the user manual (See Section II).**

## D. Storage

Store the kit at -20°C.

## E. Other Reagents Needed

- Linearized Cas9/gRNA destination vector
- Vector/gRNA-specific primers (standard desalted, non-phosphorylated oligos)
- dNTPs, Phusion (recommended) or other high fidelity polymerase, and corresponding buffers for PCR
- PCR-Grade Dimethyl Sulfoxide (Cat #D9170, Sigma)
- QIAquick PCR Purification Kit (Cat # 28106, Qiagen)
- QIAquick Gel Extraction kit (Cat # 28704, Qiagen)
- SOC or LB Broth for transformation of bacteria
- LB + 50 µg/ml Ampicillin or Kanamycin plates
- High-quality competent cells ( $>1 \times 10^9$  cfus/ug DNA)



## II. Protocol

### A. Design of Primers for PCR reaction (For Dual gRNA cloning)

Prior to the experiment, complete forward and reverse primers for generating the desired multiplex gRNA PCR amplicon will need to be designed. Once the correct sized amplicon is generated and gel-purified, it can be used as the insert in the subsequent fusion reaction with a suitable linearized destination vector.

The general format of the PCR primers for selected Cas9/gRNA vectors is shown below for generating constructs containing H1 and U6 promoters to drive gRNA expression.

**Note: For those researchers who want to use two copies of the same promoter for both gRNAs, please see the Appendix (Section IV) for primer design details.**

**For vectors with existing U6 promoter (e.g. Addgene pX series and SBI's CASLV U6-based gRNA cloning vectors):**

#### ***Forward Primer Sequence:***

5'AAAGGACGAAACACCGNNNNNNNNNNNNNNNNNNNNNNGTTT  
TAGAGCTAGAAATAGCAAG3'

#### ***Reverse Primer Sequence:***

5'TTCTAGCTCTAAAACXXXXXXXXXXXXXXXXXXXXGGATCCA  
AGGTGTCTCATAC3'

**N = Denotes gRNA1 sequence X= Denotes Reverse Complement of gRNA2 sequence**

Here is an example of two gRNAs (gRNA1 and gRNA2) that need to be cloned in to an U6-gRNA vector to make a dual gRNA vector (U6-gRNA1-H1-gRNA2)

PrecisionX™ Multiplex gRNA Cloning Kit  
Cat. # CAS9-GRNA-KIT



gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

Forward primer

5' AAAGGACGAAACACCGCCATCCGACTTCGACGATCGGTTTtagagctagaaatagcaag 3'

Reverse primer:

5' TTCTAGCTCTAAACAGAGCTCTATAGGTCGAATGGATCCAAGGTGTCTCATAC 3'

Use H1 block as template

The **H1 scaffold block** will need to be used for the PCR reaction

**For SBI's CAS8xx and 9xx vectors with H1 promoter**

***Forward Primer Sequence:***

5'CTTATAAGTTCTGTATGAGACCACTTGGATCCNNNNNNNNNN  
NNNNNNNNNNGTTTTAGAGCTAGAAATAGCAAG3'

***Reverse Primer Sequence:***

5'TTCTAGCTCTAAACXXXXXXXXXXXXXXXXXXXXCGGTGTT  
TCGTCCTTTCCAC3'

N = Denotes gRNA1 sequence X= Denotes Reverse Complement  
of gRNA2 sequence

Here is an example of two gRNAs (gRNA1 and gRNA2) that need to be cloned into these vectors to make a dual gRNA vector (H1-gRNA1-U6-gRNA2)



gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

Forward primer

5' CTTATAAGTTCTGTATGACCACTTGGATCC **CCATCCGACTTCGACGATCG** GTTTTAGAGCTAGAAATAGCAAG 3'

Reverse primer:

5' TTCTAGCTCTAAAC **AGAGCTCTATAGGTCCGAAT** CGGTGTTTCGTCCTTTCCAC 3'

Use U6 block as template

The **U6 scaffold block** will need to be used for the PCR reaction

### **For SBI's CAS7xx and CASLVxxx vectors with H1 promoter**

#### ***Forward Primer Sequence:***

5'AGACACCTTGGATCC **NNNNNNNNNNNNNNNNNNNNNN** GTTTT  
AGAGCTAGAAATAGCAAG3'

#### ***Reverse Primer Sequence:***

5'TTCTAGCTCTAAAC **XXXXXXXXXXXXXXXXXXXXX** CGGTGTT  
TCGTCCTTTCCAC3'

**N = Denotes gRNA1 sequence** **X= Denotes Reverse Complement of gRNA2 sequence**

Here is an example of two gRNAs (gRNA1 and gRNA2) that need to be cloned into these vectors to make a dual gRNA vector (H1-gRNA1-U6-gRNA2)

PrecisionX™ Multiplex gRNA Cloning Kit  
Cat. # CAS9-GRNA-KIT



gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGACCTATAGAGCTCT

Forward primer

5' AGACACCTTGGATCCCATCCGACTTCGACGATCGGTTTAGAGCTAGAAATAGCAAG 3'

Reverse primer:

5' TTCTAGCTCTAAACAGAGCTCTATAGGTCCGAATCGGTGTTTCGTCCTTCCAC 3'

Use U6 block as template

The **U6 scaffold block** will need to be used for the PCR reaction

## B. Design of Primers for PCR reaction (For “Quad-Plex” gRNA cloning)

The scalability of the Multiplex gRNA cloning kit allows for simultaneous cloning of more than 2 gRNAs at once into a single vector. This allows researchers to perform more advanced CRISPR/Cas9 techniques such as tandem double-nicking (4 gRNAs total) which allow removal of defined genomic segments using Cas9 nickase with significantly decreased chance of off-target cutting.

The cloning of four gRNAs will require the researcher to perform three separate PCR reactions with separate primer pairs and blocks. Once the correct size amplicons are generated and gel-purified, they can be mixed at equimolar ratios (1:1:1) based on their concentrations, and used as the insert in the subsequent fusion reaction with a suitable linearized destination vector.

The general format of the PCR primers for selected Cas9/gRNA vectors for this application is the following:

**For vectors with existing U6 promoter (e.g. Addgene pX series and SBI's CASLV U6-based gRNA cloning vectors):**

**Amplicon #1 (gRNA1 and gRNA2)**

***Forward Primer Sequence:***

5'AAAGGACGAAACACCGNNNNNNNNNNNNNNNNNNNNGTTT  
TAGAGCTAGAAATAGCAAG3'

***Reverse Primer Sequence:***

5'XXXXXXXXXXXXXXXXXXXXGGATCCAAGGTGTCTCATAC3'

N = Denotes gRNA1 sequence X= Denotes Reverse Complement  
of gRNA2 sequence

The **H1 scaffold block** will need to be used for the PCR reaction for generating amplicon #1.

**Amplicon #2 (gRNA2 and gRNA3)**

***Forward Primer Sequence:***

5'NNNNNNNNNNNNNNNNNNGTTTGTAGAGCTAGAAATAGCAAG3'

***Reverse Primer Sequence:***

5'XXXXXXXXXXXXXXXXXXXXCGGTGTTTCGTCCTTTCCAC3'

N = Denotes 15bp of gRNA2 sequence upstream of PAM X=  
Denotes Reverse Complement of gRNA3 sequence

The **U6 scaffold block** will need to be used for the PCR reaction for generating amplicon #2.

**Amplicon #3 (gRNA3 and gRNA4)**

***Forward Primer Sequence:***

5'NNNNNNNNNNNNNNNNNNGTTTGTAGAGCTAGAAATAGCAAG3'

### Reverse Primer Sequence:

5'TTCTAGCTCTAAAACXXXXXXXXXXXXXXXXXXXXGGATCCA  
AGGTGTCTCATAC3'

N = Denotes 15bp of gRNA3 sequence upstream of PAM X=  
Denotes Reverse Complement of gRNA4 sequence

The **H1 scaffold block** will need to be used for the PCR reaction for generating amplicon #3.

For example, given 4 gRNAs below which need to be cloned into a single destination vector, the following primers will need to be designed to generate three separate PCR amplicons:



gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

gRNA3: TACGGCAATGCAGCATTCTGA

gRNA4: CTACTGGCTAACCGTTAGCT

Forward primer-1

5' AAAGGACGAAACACCGCATCCGACTTCGACGATCGGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-1

5' AGAGCTCTATAGGTCCGAATGGATCCAAGGTGTCTCATAC 3'

H1 block as template

Forward primer-2

5' GACCTATAGAGCTCTGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-2

5' TCGAATGCTGCATTGCCGTACGGTGTTCGTCTTCCAC 3'

U6 block as template

Forward primer-3

5' CAATGCAGCATTCTGAGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-3

5' TTCTAGCTCTAAAACAGCTAACGGTTAGCCAGTAGGGATCCAAGGTGTCTCATAC 3'

H1 block as template

### For SBI's CAS8xx and 9xx vectors with H1 promoter

#### Amplicon #1 (gRNA1 and gRNA2)

#### Forward Primer Sequence:

5'CTTATAAGTTCTGTATGAGACCACTTGGATCCNNNNNNNNNN  
NNNNNNNNNNGTTTTAGAGCTAGAAATAGCAAG3'

**Reverse Primer Sequence:**

5'XXXXXXXXXXXXXXXXXXXXCGGTGTTTCGTCCTTTCCAC3'

N = Denotes gRNA1 sequence X= Denotes Reverse Complement  
of gRNA2 sequence

The **U6 scaffold block** will need to be used for the PCR reaction for generating amplicon #1.

**Amplicon #2 (gRNA2 and gRNA3)**

**Forward Primer Sequence:**

5'NNNNNNNNNNNNNNNNGTTTTAGAGCTAGAAATAGCAAG3'

**Reverse Primer Sequence:**

5'XXXXXXXXXXXXXXXXXXXXGGATCCAAGGTGTCTCATAC3'

N = Denotes 15bp of gRNA2 sequence upstream of PAM X=  
Denotes Reverse Complement of gRNA3 sequence

The **H1 scaffold block** will need to be used for the PCR reaction for generating amplicon #2.

**Amplicon #3 (gRNA3 and gRNA4)**

**Forward Primer Sequence:**

5'NNNNNNNNNNNNNNNNGTTTTAGAGCTAGAAATAGCAAG3'

**Reverse Primer Sequence:**

5'TTCTAGCTCTAAACXXXXXXXXXXXXXXXXXXXXCGGTGTT  
TCGTCCTTTCCAC3'

N = Denotes 15bp of gRNA3 sequence upstream of PAM X=  
Denotes Reverse Complement of gRNA4 sequence

## PrecisionX™ Multiplex gRNA Cloning Kit

Cat. # CAS9-GRNA-KIT

The **U6 scaffold block** will need to be used for the PCR reaction for generating amplicon #3.

For example, given 4 gRNAs below which need to be cloned into a single destination vector, the following primers will need to be designed to generate three separate PCR amplicons:



gRNA1: CCATCCGACTTCGACGATCG  
gRNA2: ATTCGGACCTATAGAGCTCT  
gRNA3: TACGGCAATGCAGCATTCTGA  
gRNA4: CTA CTGGCTAACCGTTAGCT

Forward primer-1

5' CTTATAAGTTCTGTATGAGACCACTTGGATCC**CCATCCGACTTCGACGATCG**GTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-1

5' **AGAGCTCTATAGGTCGAAT**CGGTGTTTCGTCCTTCCAC 3'

} U6 block as template

Forward primer-2

5' **GACCTATAGAGCTCT**GTTTATAGCTAGAAATAGCAAG 3'

Reverse primer-2

5' **TCGAATGCTGCATTGCCGTA**GGATCCAAGGTGTCTCATAC 3'

} H1 block as template

Forward primer-3

5' **CAATGCAGCATTCTGA**GTTTATAGCTAGAAATAGCAAG 3'

Reverse primer-3

5' TTCTAGCTCTAAAC**AGCTAACGGTTAGCCAGTAG**CGGTGTTTCGTCCTTCCAC 3'

} U6 block as template

### For SBI's CAS7xx and CASLVxxx vectors with H1 promoter

#### Amplicon #1 (gRNA1 and gRNA2)

##### **Forward Primer Sequence:**

5'AGACACCTTGGATCC**NNNNNNNNNNNNNNNNNNNNNN**GTTT  
AGAGCTAGAAATAGCAAG3'

##### **Reverse Primer Sequence:**

5'**XXXXXXXXXXXXXXXXXXXX**CGGTGTTTCGTCCTTCCAC3'



N = Denotes gRNA1 sequence X= Denotes Reverse Complement of gRNA2 sequence

The **U6 scaffold block** will need to be used for the PCR reaction for generating amplicon #1.

### **Amplicon #2 (gRNA2 and gRNA3)**

#### ***Forward Primer Sequence:***

5' NNNNNNNNNNNNNNNNNNNNN GTTTTAGAGCTAGAAATAGCAAG3'

#### ***Reverse Primer Sequence:***

5' XXXXXXXXXXXXXXXXXXXX GGATCCAAGGTGTCTCATAC3'

N = Denotes 15bp of gRNA2 sequence upstream of PAM X= Denotes Reverse Complement of gRNA3 sequence

The **H1 scaffold block** will need to be used for the PCR reaction for generating amplicon #2.

### **Amplicon #3 (gRNA3 and gRNA4)**

#### ***Forward Primer Sequence:***

5' NNNNNNNNNNNNNNNNNNNNN GTTTTAGAGCTAGAAATAGCAAG3'

#### ***Reverse Primer Sequence:***

5' TTCTAGCTCTAAAC XXXXXXXXXXXXXXXXXXXX CGGTGTT TCGTCCTTCCAC3'

N = Denotes 15bp of gRNA3 sequence upstream of PAM X= Denotes Reverse Complement of gRNA4 sequence

The **U6 scaffold block** will need to be used for the PCR reaction for generating amplicon #3.

## PrecisionX™ Multiplex gRNA Cloning Kit

Cat. # CAS9-GRNA-KIT

For example, given 4 gRNAs below which need to be cloned into a single destination vector, the following primers will need to be designed to generate three separate PCR amplicons:



gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

gRNA3: TACGGCAATGCAGCATTCTGA

gRNA4: CTACTGGCTAACCGTTAGCT

Forward primer-1

5' AGACACCTTGGATCCCATCCGACTTCGACGATCGGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-1

5' AGAGCTCTATAGGTCGGAATCGGTGTTTCGTCTTTCCAC 3'

} U6 block as template

Forward primer-2

5' GACCTATAGAGCTCTGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-2

5' TCGAATGCTGCATTGCCGTAGGATCCAAGGTGTCTCATAC 3'

} H1 block as template

Forward primer-3

5' CAATGCAGCATTCTGAGTTTATAGAGCTAGAAATAGCAAG 3'

Reverse primer-3

5' TTCTAGCTCTAAACAGCTAACGGTTAGCCAGTAGCGGTGTTTCGTCTTTCCAC 3'

} U6 block as template

## C. Setting up the PCR Reaction

1. Once the primer sets and the appropriate blocks needed for the reactions have been determined, set up the PCR reaction according to the table below (for 1 reaction):

Component	Volume
5X Phusion HF buffer	10 µl
10 mM dNTPs	1 µl
10 µM Forward Primer	2.5 µl
10 µM Reverse Primer	2.5 µl
H1 block or U6 block	0.5 µl
DMSO	1.5 µl
Phusion DNA Polymerase	0.5 µl
Nuclease-free water	To 50 µl

Note: SBI recommends the use of Phusion DNA polymerase and buffer for the PCR reaction as this combination was utilized successfully for in-house validation.

2. Run the PCR reaction according to the cycling conditions below.

Step	Temperature	Time
<u>Initial Denaturation</u>	98°C	3 min
30 cycles	98°C	5 seconds
	51°C	15 seconds
	72°C	15 seconds
Final Extension	72°C	10 minutes
Hold	4°C	

3. After completion of the PCR reaction, run out 2-5 µl of the PCR product on a 2% agarose gel to determine specificity of the PCR

reaction. If there are no additional bands (e.g. primer dimer, other non-specific bands) visible, the PCR reaction may be cleaned up using the QIAquick PCR Purification Kit (Cat # 28106, Qiagen), concentration measured by UV-Vis spectrophotometry, and directly used for the fusion reaction.

If primer dimers or non-specific products are seen, we highly recommend gel-purification of the correct band to avoid their carryover into the fusion reaction as these products will inhibit the efficiency of the reaction.

## **D. Fusion Reaction**

### **IMPORTANT!**

Complete linearization of the destination vector is critical to achieve a successful fusion reaction. Incomplete linearization of the vector will result in high background. The linearized vector can be generated by PCR (e.g. inverse PCR) or restriction enzyme digest (single or double digest) and should be purified using either a gel or PCR purification kit.

Due to the digestion efficiency, different restriction enzymes will generate different levels of background. In general, two enzyme digestion is better than a single enzyme digestion. The further the restriction sites are apart, the better the digestion efficiency. Increasing the enzyme digestion time and the digestion reaction volume will also help reduce the background. For many enzymes, we recommend incubate the digestion reaction between 3 hours and overnight in order to increase linearization and reduce background (w/exception of certain high-fidelity “HF” enzymes).

Check the background of your vector by transforming 1µl of 10-100ng/µl linearized vector into competent cells. If the background is high, continue digesting the remaining vector for a longer time after addition of more restriction enzyme(s).

We recommend digesting 2µg vector in 50µl reaction overnight. Use QIAGEN's QIAquick Spin Gel Extraction kit for gel purification and elute the DNA with 30µl dH<sub>2</sub>O.

Set up the following reaction in a 1.5 ml sterile reaction tube by mixing the following reagents gently and then spin down briefly to collect the reagents at the bottom of the tube.

### **Fusion reaction**

Linearized destination vector (10-100ng/µl)	1 µl*
PCR insert(s) (20-200ng/µl) (for each PCR Product)	1 µl*
dH <sub>2</sub> O	—µl
<u>5x master mix</u>	<u>2µl</u>
Total	10µl

### **Positive control reaction**

Linearized vector (positive control)	1 µl
500bp PCR insert (positive control)	1 µl
dH <sub>2</sub> O	6µl
<u>5x master mix</u>	<u>2µl</u>
Total	10µl

### **Negative control reaction**

Linearized destination vector (10-100ng/µl)	1 µl*
dH <sub>2</sub> O	7µl
<u>5x master mix</u>	<u>2µl</u>
Total	10µl

\* 2:1 or 1:1 molar ratio of insert: vector works well in the fusion reaction.

When using the kit for the first time, we strongly recommend that you perform the positive and negative control reactions in parallel with your fusion reaction. The positive control 500bp PCR insert and linearized vector provided in the kit have already been

PrecisionX™ Multiplex gRNA Cloning Kit  
Cat. # CAS9-GRNA-KIT

purified. There is no treatment (e.g. phosphatase, etc.) needed prior to the cloning reaction.

**Fusion Reaction Incubation**

1. 5 minutes at room temperature
2. 10 minutes on ice


**E. Transformation of Cells**

1. Add 30-50µl of competent cells ( $> 1 \times 10^9$  cfus/ug DNA) to the cloning mixture per instructions provided by manufacturer of competent cells.
2. Incubate on ice for 20 minutes
3. Heat shock at 42°C for 50 seconds
4. Transfer on ice for 2 minutes
5. Add 250µl S.O.C medium or LB broth
6. Incubate at 37°C for an hour
7. Take 100µl culture spread on pre-warmed (37°C) culture plate containing 50µg/ml Ampicillin or Kanamycin
8. Incubate the plate at 37°C overnight


**F. Confirmation of Positive Clones**

1. Pick 1 to 2 colonies, grow in LB/Kanamycin or Ampicillin medium overnight at 37°C with shaking
2. Next day, miniprep plasmid DNAs and check for proper ligation via restriction digestion using AgeI and EcoRV (for SBI's CAS7xx, CAS8xx, CAS9xx vectors), or AgeI and Stul (for SBI's CASLV5xx vectors). For non-SBI vectors, please find suitable unique restriction sites flanking gRNA cloning sites.

CAS7xx, CAS800, CAS900 digest with Agel and EcoRV



CASLV5xx digest with AgeI and StuI



Releasing Band: ~0.9kb    ~1.3kb    ~1.7kb    ~2.1kb

The correct sized restriction product should indicate successful cloning of the multi-gRNA construct into the destination vector.

3. Send out finished construct for sequencing using appropriate sequencing primer for the construct (see Appendix for primer information).



### III. Validation Data

#### A. Gene repair by double-nicking and HR







Fig. 4: Functional validation of dual gRNAs generated by the Multiplex gRNA Cloning Kit for combined double-nickase and homologous recombination to rescue EGFP in the EGIP control cell line. (EGIP cell line courtesy of Dr. Jizhong Zou, NIH CRM.)

**B. Gene Knockout by dual gRNAs targeting RFP**




Pictures were taken 10 days after transfection


Fig. 5: Functional validation of dual gRNAs cloned into SBI's CAS940A-1 vector (CMV-hspCas9-H1-gRNA) using the Multiplex gRNA Cloning Kit to targeting RFP in a cell line stably expressing a CMV-GFP-T2A-RFP cassette. PCR assay with primers flanking

RFP indicate evidence of cutting by gRNAs, and fluorescence images of targeted cells vs control show reduction in RFP expression.

C. Tandem paired-nicking (4 gRNAs in nickase vector)



Genomic DNA was extracted 5 days after transfection



Pictures were taken 10 days after transfection

Fig. 6: Functional validation of “quad-plex” gRNAs cloned into CAS840A-1 nickase vector (EF1-Nickase-H1-gRNA) using the

PrecisionX™ Multiplex gRNA Cloning Kit  
Cat. # CAS9-GRNA-KIT

Multiplex gRNA Cloning Kit for targeting GFP and RFP in a cell line stably expressing a CMV-GFP-T2A-RFP cassette. PCR assay with primers (one at 5' end of GFP and other at 3' end of RFP) indicate evidence of cutting by gRNAs. Fluorescence images of targeted cells vs control show significant reduction in both RFP and GFP expression.

## IV. Appendix

### A. Primer design for dual gRNA/single promoter constructs

For customers interested in having two copies of the same promoter driving both gRNAs, the following sets of PCR primers can be used to generate U6-U6 or H1-H1 constructs for a wide variety of different CRISPR/Cas9 expression vectors.

#### For U6-U6 constructs (e.g. Addgene pX series and SBI's CASLV U6-based gRNA cloning vectors)



Sample gRNAs to be cloned:

gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

#### **Forward Primer Sequence (gRNA1 sequence in light blue)**

5'AAAGGACGAAACACCGCCATCCGACTTCGACGATCGGTTTT  
AGAGCTAGAAATAGCAAG3'

#### **Reverse Primer Sequence (Reverse complement of gRNA2 sequence in red)**

5'TTCTAGCTCTAAAACAGAGCTCTATAGGTCCGAATCGGTGT  
TTCGTCTTTCCAC3'

**Note: Use U6 Block for PCR reaction**

**For H1-H1 constructs (e.g. SBI's CAS8xx and 9xx vectors)**



Sample gRNAs to be cloned:

gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

***Forward Primer Sequence (gRNA1 sequence in light blue)***

5'CTTATAAGTTCTGTATGAGACCACTTGGATCCCATCCGACT  
TCGACGATCGGTTTTAGAGCTAGAAATAGCAAG3'

***Reverse Primer Sequence (Reverse complement of gRNA2  
sequence in red)***

5'TTCTAGCTCTAAAACAGAGCTCTATAGGTCCGAAT  
GGATCCAAGGTGTCTCATAC3'

**Note: Use H1 Block for PCR reaction**

**For H1-H1 constructs (e.g. SBI's CAS7xx and CASLVxxx H1-  
based vectors)**



Sample gRNAs to be cloned:

gRNA1: CCATCCGACTTCGACGATCG

gRNA2: ATTCGGACCTATAGAGCTCT

***Forward Primer Sequence (gRNA1 sequence in light blue)***

PrecisionX™ Multiplex gRNA Cloning Kit  
Cat. # CAS9-GRNA-KIT

5'AGACACCTTGGATCC**CCATCCGACTTCGACGATCG**GTTTTA  
GAGCTAGAAATAGCAAG3'

***Reverse Primer Sequence (Reverse complement of gRNA2 sequence in red)***

5'TTCTAGCTCTAAAAC**AGAGCTCTATAGGTCCGAAT**  
GGATCCAAGGTGTCTCATAC3'

**Note: Use H1 Block for PCR reaction**

## **B. Sequencing primer information**

The following primer sequences can be used for sequence-verification of inserted gRNAs using the Kit.

For all U6-H1 constructs:

U6 primer: 5' GGAATATCATATGCTTACCG 3'

For all H1-U6 constructs:

H1 primer: 5' GTCATCAACCCGCTCCAAGG 3'

For non-SBI U6-U6 and H1-H1 constructs:

Please find a suitable primer sequence upstream of 5' most U6 or H1 promoter based on vector sequence.

For SBI U6-U6 or H1-H1 constructs:

GRNA-PR-F: 5' GACGAGTCGGATCTCCCTTTG 3'

GRNA-PR-R: 5' CAAGTTCCTCTCACTCTCTGAT 3'

**NOTE: The primer set above (GRNA-PR-F/R) is available from SBI (Cat# CAS9-GRNA-PR)**

## V. Troubleshooting

Problem	Probable cause	Solution
1. No or few colonies obtained from the transformation	Primer sequences are incorrect	Check primer sequences to ensure that the gRNA sequence is correct and primer ends are homologous to vector and block, respectively
	Suboptimal PCR product	Optimize your PCR amplification reactions so that you generate pure PCR products.  Use a different method to purify your PCR product.
	Low DNA concentration in reaction	It is imperative to obtain as high a DNA concentration as possible in your fusion reaction.
	Inhibitory contaminants from PCR product or linearized vector	Both the PCR product and the linearized vector should be purified.
	Transform with too much reaction mixture	Do not add more than 10 $\mu$ l of reaction mixture to 50 $\mu$ l of competent cells. Too much reaction mixture inhibits the transformation.

Problem	Probable cause	Solution
	Low quality or poor handling of competent cells	<p>Handle the competent cells gently. Do not re-freeze cells after thawing.</p> <p>Quality of competent cells may be tested by transforming a circular plasmid to determine cells' competency. Competent cells with a transformation efficiency of <math>1 \times 10^9</math> cfu/ <math>\mu</math>g are recommended.</p>
	Wrong antibiotic or too much antibiotic in the media	Choose plates with the appropriate concentration of the right antibiotic.
2. Large numbers of colonies contain no insert	Incomplete linearization of your vector	It is critical to remove any uncut vector prior to use in the fusion reaction. If necessary, re-digest your vector and gel purify.



Problem	Probable cause	Solution
2. Large numbers of colonies contain no insert	Plates are too old or contained incorrect antibiotic	Make sure that your antibiotic plates are fresh. Check the antibiotic resistance of your fragment.
3. Clones contain incorrect insert	PCR products contain non-specifically amplified artifacts	Optimize your PCR reaction to improve the specificity. Screen more colonies for the correct clones.

## VI. References

Ho *et al.* Targeting non-coding RNAs with the CRISPR/Cas9 system in human cell lines. *Nucleic Acids Research* 2014; Nov 20 [Epub ahead of print]

## VII. Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site:

<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

**Phone:** (650) 968-2200  
(888) 266-5066 (Toll Free)

**Fax:** (650) 968-2277

**E-mail:** General Information: [info@systembio.com](mailto:info@systembio.com)  
Technical Support: [tech@systembio.com](mailto:tech@systembio.com)  
Ordering Information: [orders@systembio.com](mailto:orders@systembio.com)

**Address:** System Biosciences (SBI)

265 North Whisman Rd.  
Mountain View, CA 94043

## **VIII. Licensing and Warranty**

### **Limited Use License**

Use of the PrecisionX™ Multiplex gRNA Cloning Kit (i.e., the “Product”) is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms. Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein. SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

### **Limited Warranty**

SBI warrants that the Product meets the specifications described in the accompanying Product Analysis Certificate. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2015 System Biosciences (SBI), All Rights Reserved.