

EVery EV RNA Isolation Kit

Cat # EVery100B-1

User Manual

Storage:

Store RNase-free DNase I and Glycogen at -20°C Store all other reagents at room temperature

Version 1 8/25/2022 A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the License and Warranty Statement contained in this user manual.

Contents

Product Description	2
List of Components	3
Additional Required and Optional Equipment Not Included in Kit	3
Protocol	4
Example Data and Applications	6
Technical Support	10
Licensing and Warranty Statement	10

Product Description

Discover more when you capture total EV RNA

EV samples vary drastically in their purity and concentration. Some EV samples from sera of diseased subjects and plasma bring up the challenge of heavy load of contaminating proteins that can potentially lead to column clogging thus lower yields. New column dimensions work just as smoothly for RNA isolation from all types of EV samples, irrespective of upstream EV isolation methods. Overcoming many of the other challenges with RNA isolation from extracellular vesicles (EVs), EVery EV RNA Isolation Kit is able to capture total EV RNA, including small RNAs. It is effective even with low amounts of input RNA and is capable of delivering high yields of highly pure RNA. Move quickly and confidently with exoRNA isolation that's high-yield and completed in <30 minutes.

- Find what others miss when you capture every RNA
- Achieve phenol-level yields with a safer column-based method
- Get more RNA for each downstream reaction
- Ensure delivery of highly pure RNA by using the included DNase I
- Compatible with most downstream applications, including RNA-seq and miRNA profiling

EVery EV RNA Isolation Kit can be used with EVs isolated by commonly used methods, including ExoQuick, SmartSEC, and ultracentrifugation, and comes with sufficient reagents to perform 20 purification reactions.

EVery EV RNA Isolation Kit delivers high yields of highly concentrated RNA from already isolated EVs. The columnbased workflow is easy to implement and can be completed in less than 30 minutes (Figure 1).



Figure 1. The quick and easy EVery EV RNA isolation workflow.

List of Components

Components	Qty/Volume	Storage Temperature
Lysis buffer	25 ml	RT
Lysis additive	3 ml	RT
Glycogen	100 μL	-20°C
Wash solution	18 ml	RT
Elution solution	1 ml	RT
RNA Spin column	20	RT
Collection tubes	20	RT
Elution Tubes	20	RT
RNase free DNase I	200 μL	-20°C
Enzyme buffer	1.5 ml	RT

Table 2. Components of EVery100B-1, EVery EV RNA Isolation Kit

NOTE: The table above is for the 20 reaction kit.

Additional Required and Optional Equipment Not Included in Kit

- 1. 96-100% Ethanol
- 2. β-mercaptoethanol (cat# M3148-25ML, Sigma)- optional, but highly recommended

Protocol

Before you start the protocol for exosomal RNA isolation:

- The protocol is outlined for 700 μL (SmartSEC Single) or 500 μL (ExoQuick, ExoQuick TC, ExoQuick ULTRA, ExoQuick TC ULTRA) of sample input. If processing a sample volume lower than 700 μL (500 μL), simply bring the volume of your sample up to 700 μL (500 μL) using 1X PBS and proceed as outlined below. The pellet from ExoQuick or ExoQuick TC should be resuspended in 1xPBS to 500 μL to prevent column clogging during isolation.
- All steps should be performed at room temperature and all centrifugation steps performed at ≥8000 x g (≥10,000 rpm for table top microcentrifuge).
- 3. It is highly recommended to warm up **Lysis Buffer** at 60°C for 5-10 minutes and mix well until the solution becomes clear again if precipitates are present.

! OPTIONAL (highly recommended)

The use of β -mercaptoethanol in the Lysis Buffer is highly recommended. Add 10 μ L of β -mercaptoethanol to each 1 mL of Lysis buffer.

! OPTIONAL (highly recommended)

Add 5 μ L of Glycogen to the Lysis Mix if you are expecting a low RNA yield.

- 4. Prepare a working solution of the **Wash solution** by adding 42 mL of 96-100% ethanol to the supplied bottle containing the concentrated Wash Solution.
- 5. The RNA yield may be increased by using **Elution Buffer** warmed to 60°C.

RNA isolation steps:

- 1. Add 1 mL (700 μL) of **Lysis Buffer**, 150 μL (110 μL) of **Lysis Additive** and 10 μL (7 μL) of β-mercaptoethanol (optional) to the 700 μL (500 μL) PBSx1 Buffer containing the purified exosomes.
- 2. Mix well by vortexing for 10 sec. then incubate at RT for 10 min.

! OPTIONAL (highly recommended)

Add 5 μ L of Glycogen to the lysis mix if you know that RNA yield is low.

- 3. After incubation add 1.85 mL (1.35 mL) of 96%-100% EtOH to the mix from Step3 and mix well by vortexing for 10 seconds.
- 4. Transfer 750 μL of the mixture from Step 4 into a Micro Spin column. Centrifuge for 1 minute. Discard the flowthrough and reassemble the spin column in its collection tube.
- 5. Repeat Step 5 four times to transfer the remaining mixture from Step 4 into the Micro Spin column.
- 6. Apply 400 μL of **Wash Solution** on the column and centrifuge for 2 minutes. Discard the flowthrough and reassemble the spin column in its collection tube.

! OPTIONAL (highly recommended)

On-column DNA removal:

a. For every on-column reaction prepare a mix of 7.5 µL of **RNase-free DNase I** and 50 µL of **Enzyme buffer**. Mix gently by inverting the tube a few times or flicking the tube with your fingers to mix.

! DO NOT VORTEX

- b. Apply 57 µL of DNase I mix from step a. to the column and incubate at 25°C-30°C for 15 minutes.
- 7. Apply 600 μL of **Wash solution** to the column and centrifuge for 30 seconds. Discard the flowthrough and reassemble the spin column in its collection tube.
- 8. Repeat step 9 one more time, for a total of two washes.
- 9. Centrifuge the empty column for 1 minute to completely remove any residual solution. Discard the collection tube.
- 10. Transfer the spin column to a fresh Eppendorf tube. Apply 30 μ L of **Elution solution** to the column and let it stand for 2 minutes. Centrifuge for 1 minute to elute.
- 11. To maximize the recovery of the RNA add the eluate collected from step 10 back on the column, let it stand for 2 minutes. Centrifuge for 1 minute to elute.
- 12. Exosomal RNA is now ready for downstream applications.

Example Data and Applications

EVery RNA captures Everything

To demonstrate the ability of the EVery EV RNA isolation Kit to capture the full range of RNAs, we used the kit to isolate RNA from 10,000 cells (Figure 2, lane 1), from EVs that were isolated from 250 µL of serum using <u>SmartSEC</u> <u>Single</u> (Figure 2, lane 2), and from buffer spiked with 0.1 pmol of Cel-miR-39 (Figure 2, lane 3). The high quality of the isolated RNA can be seen in lane 1, where the RNA integrity number (RIN) is 9.9 and the 28S/18S RNA ratio is 1.5. The multiple bands in lane 2 demonstrate that EVery captures RNAs of different lengths—EVerything—from EVs with no apparent bias or size preference. The strong signal from the spiked-in miRNA in lane 3 demonstrates the good recovery of even small RNAs.



Figure 2. EVery EV RNA Isolation Kit captures EVerything.

The compatibility of multiple EV isolation techniques with EVery EV RNA Isolation Kit and the excellent size distribution of RNAs isolated from those EVs is shown in Figure 3.



Figure 3. EVery EV RNA Isolation Kit is compatible with EVs isolated using ExoQuick Ultra and SmartSEC Single.

EVery EV RNA Isolation Kit delivers similar amounts of RNA as phenol-based methods

To demonstrate the excellent RNA yields and robust cDNA synthesis obtained with EVery family of products, we isolated EVs from 250 µL of serum using <u>SmartSEC Single</u>, spiked in 0.1 pmol of Cel-miR-39, and used both <u>EVery EV</u> <u>RNA Isolation Kit</u> and a phenol-based kit to isolate RNA. The isolated RNA was reverse transcribed using EVery cDNA Synthesis Kit and the copy number of Cel-miR-39 measured (Figure 4). EVery EV RNA Isolation Kit delivered similar levels of Cel-miR-39 as the phenol-based method.



Figure 4. EVery EV RNA Isolation is EVery bit as good as phenol.

EVery EV RNA Isolation Kit efficiently isolates mRNA

We used the EVery EV RNA Isolation Kit and EVeryRNA cDNA Synthesis Kit to isolate mRNA and synthesize cDNA from cells overexpressing eGFP (Figure5). Robust levels of eGFP mRNA are recovered and converted to cDNA when cells are overexpressing eGFP.



Figure 5. EVery efficiently isolates longer RNAs like mRNA

miRNA isolated with EVery100B-1 and converted to cDNA using EVery200B-1 can be used for miRNA profiling

We isolated EVs from 250 µL of normal or breast cancer patient serum using <u>SmartSEC Single</u> and used EVery EV RNA Isolation Kit and EVery cDNA Synthesis Kit to isolate and reverse transcribe EV RNAs for miRNA profiling using <u>EVery miRNome Profiler for Human Serum and Plasma (EVery500B-1)</u>. Distinct miRNA expression profiles of breast cancer patient serum compared to normal serum are shown in Figure 6. Differential expression of selected breast cancer associated miRNA markers were further confirmed by qPCR (Figure 7).



Figure 6. EVery family of products generate high-quality cDNA suitable for miRNA profiling



Figure 7. Differential expression of selected miRNAs in normal vs. breast cancer serum confirmed by qPCR. Data normalized with global mean. Expression of miR-451a, miR-21-5p, miR-16-5p and miR-103a are significantly upregulated, while expression of the tumor suppressor miR-124-3p is significantly downregulated in breast cancer serum compared to normal serum.

We were also able to show robust, successful RNA-seq runs using RNA isolated from EVs with EVery EV RNA Isolation Kit (Table 3). All three EV isolation methods tested generated high-quality RNA-seq data.

EV isolation method	Amount of RNA isolated (ng)	Number of reads	FAST-QC
ExoQuick	4.3	53,349,528	Passed
ExoQuick Ultra	5.4	107,154,128	Passed
SmartSEC Single	5.2	98,886,924	Passed

Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site: <u>http://www.systembio.com</u>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI) 2438 Embarcadero Way Palo Alto, CA 94303

Phone:	(650) 968-2200
Toll-Free:	(888) 266-5066
Fax	(650) 968-2277

E-mail:

General Information:info@systembio.comTechnical Support:tech@systembio.comOrdering Information:orders@systembio.com

Licensing and Warranty Statement

Limited Use License

Use of the EVery EV RNA Isolation Kit (*i.e.*, the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

This Product is covered by licensed Patents. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2020 System Biosciences (SBI), All Rights Reser

This page intentionally left blank.



System Biosciences (SBI) 2438 Embarcadero Way Palo Alto, CA 94303

Phone:	(650) 968-2200
Toll-Free:	(888) 266-5066
Fax:	(650) 968-2277

E-mail:

General Information:info@systembio.comTechnical Support:tech@systembio.comOrdering Information:orders@systembio.com