


Retro-Concentin™

Easily concentrate Retroviruses without ultracentrifugation or expensive columns

Retro-Concentin™ precipitates retroviruses quantitatively without ultracentrifugation or complicated column procedures. Instead, retroviruses are directly pelleted from culture medium with simple one/two steps of low speed centrifugations. In addition, the concentration solution stabilizes the retroviruses for frozen storage, which may provide another advantage over other methods. Each preparation can handle up to 200 ml of retroviral supernatant (most centrifuges) and the resulting pellet can be dissolved in a desired volume to meet your experimental requirements. Cells transduced show normal morphology and no obvious cytotoxicity is observed.




Step 1.
Add 5x Retro-Concentin to Retrovirus Supernatant

Step 2.
Centrifuge at 1500xg for 30 min. at 4°C

Step 2B (Optional).
Transfer to an Eppendorf tube and Centrifuge at 12,000 rpm for 2 min.

High titer retrovirus particles prepared with Retro-Concentin were successfully used for generating Induced Pluripotent Stem Cells (iPSCs).


Use concentrated retrovirus to make iPSCs efficiently



Step 3.
Aliquot and store at -70°C



System Biosciences
1616 North Shoreline Blvd.
Mountain View, CA 94043
Toll Free: 888.266.5066
Fax: 650-968-2277
Email: info@systembio.com



Highlights

- Greater than 80% recovery
- Quick and easy protocol
- Flexible volume options
- Rapid 15 minute protocol
- No ultracentrifugation or columns required
- Stabilizes virus for storage
- High performance particles can reprogram cells into iPSCs

Deliver more Virus using Retro-Concentin™

