



AAVanced Concentration Reagent

Cat# AAV100-SAM, AAV100A-1, AAV110A-1

User Manual

Store the reagent at room temperature

Version 8
8/4/2022

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the License and Warranty Statement contained in this user manual.

Contents

Product Description.....	1
List of Components.....	1
Storage.....	2
Equipment to be supplied by user.....	2
Protocol:	2
AAV Concentration Protocol.....	2
General rAAV Transduction Protocol	3
Example Data and Applications	4
References	5
Technical Support	6
Licensing and Warranty Statement	6

Product Description

Adeno-associated virus (AAV) is a single-strand DNA virus belonging to the *Parvoviridae* family. AAV can infect a wide range of cell types, including both dividing and quiescent cells. In humans, AAV infection is not associated with any known diseases, thus it has been widely adopted for delivery of recombinant DNA for *in vivo* applications. Recombinant adeno-associated virus (rAAV) vectors have been developed via removal of the AAV packaging signal from the virus genome. While the native virus integrates into a specific genomic locus in the host cell (*AAVS1*, located in chromosome 19q13.13), rAAV vectors persist in an extrachromosomal state with a very low frequency of random integration. Therefore, rAAV vectors have been broadly used in gene therapy and genome engineering as an alternative to other viral gene delivery methods.

Traditional rAAV production methods requires multiple steps, including cell lysis and CsCl₂ ultra-high-speed density gradient centrifugation, chromatography, or binding to affinity matrix columns. These are difficult to setup, time-consuming, and require specialized equipment for isolation of high-purity rAAV for *in vivo* experiments. Many studies have shown that several rAAV serotypes are efficiently released into the culture medium after co-transfection of the plasmids required to produce rAAV particles in human embryonic kidney 293 cells. To take full advantage of this property of rAAV, SBI has developed an innovative rAAV concentration solution called AAVanced™ Concentration Reagent – a simple, one step rAAV concentration reagent for the isolation of rAAV particles from media. AAVanced Concentration Reagent is based on a proven nanoparticle technology and is specifically optimized for the precipitation of most AAV serotypes from culture medium. The rAAV virus produced using AAVanced Concentration Reagent has been tested successfully for *in vitro* and *in vivo* applications for multiple serotypes of AAV with no observed cytotoxic effects, which provides direct evidence for the utility of this reagent for demanding applications using rAAV.

SBI's AAVanced Concentration Reagent significantly reduces the complexity and time required for rAAV virus production, allowing researchers to focus on their research experiments rather than virus production.

List of Components

Catalog Number	Description	Amount	Storage Condition
AAV100-SAM	AAVanced Concentration Reagent	5 mL	Room temperature
AAV100A-1	AAVanced Concentration Reagent	100 mL	Room temperature
AAV110A-1	AAVanced Concentration Reagent	250 mL	Room temperature

Storage

The kits are shipped at room temperature or on blue ice and should be stored at RT upon receipt. Properly stored kits are stable for 1 year from the date received.

Equipment to be supplied by user

1. HEK293TN cells (Cat #LV900A-1, System Biosciences). Other low passage HEK293T or FT cells can also be used if available.
2. AAV shuttle and packaging vectors (including plasmids encoding for Rep/Cap protein and helper factors (e.g. Adenovirus E2A, VA, and E4 proteins)
3. 150mm² tissue-culture plates
4. PureFection transfection reagent (Cat #LV750A-1, System Biosciences). Other reagents such as Lipofectamine 2000 may also be suitable.

Protocol:

AAV Concentration Protocol

AAVanced™ Concentration Reagent is provided in a 5x solution.

1. Seed 7×10^6 HEK293TN cells in 150mm² plate(s) with 20ml of complete growth medium (i.e. DMEM, 10% FBS) without antibiotics to reach 70-80% confluency within 24hrs.
2. Transfect cells with manufacturer recommended amounts of AAV shuttle, Rep/Cap, and AAV helper plasmids for 150mm² plates using PureFection™ (Cat # LV750A-1, please visit SBI website for detailed protocol)
3. Next day after transfection, change the media to complete growth media including antibiotics.
4. 72 hours after initial transfection, transfer supernatant with free rAAV virus to a sterile 50ml conical vessel. If using multiple plates, the supernatants can be combined together before addition of AAVanced™ Concentration Reagent
5. Centrifuge at $3000 \times g$ for 15 minutes to remove cells and cell debris.
6. Add 1 volume of cold AAVanced™ Concentration Reagent (4°C) to every 4 volumes of AAV particles-containing supernatant. Mix thoroughly with pipetting (For example, add 2 ml AAVanced™ Concentration Reagent for every 8 ml of supernatant)

7. Refrigerate overnight (or at least 12 hours). AAV viral particles-containing supernatants mixed with AAVanced™ Concentration Reagent are stable for up to 4-5 days at 4°C.
8. Centrifuge supernatant/ precipitation mixture at 1500 × g (or 3000 rpm) for 30 minutes at 4°C. After centrifugation, the AAV viral particles may appear as a beige or white pellet at the bottom of the vessel.
9. Remove all of the supernatant (or as much as possible) and spin the tube containing AAV pellet again 1500 × g (or 3000 rpm) for 3 minutes at 4°C. Completely remove all traces of fluid by aspiration.
10. Resuspend viral pellets in 1/1000 (for use in vivo) to 1/100 (for use in vitro) of original volume using cold, sterile Phosphate Buffered Saline (PBS) or DMEM containing 25mM HEPES buffer at 4°C. Titer the concentrated AAV with a method of your choice.
11. Aliquot into cryogenic vials and store at -80°C until ready for use.

Precipitation of AAV virus particles from large volumes can be achieved by using the Corning 250 mL polypropylene centrifuge tube (Cat. # 430776), following manufacturer's instructions.

General rAAV Transduction Protocol

A general protocol for transduction of rAAV particles is shown below. This protocol can be scaled upwards with respect to vessel size.

1. Day 1: Plate 10,000 cells per well in a 24 well plate (about 30% confluence) in culture medium.
2. Day 2: Add virus to each well at different volumes based on the experimental needs (typical amount of virus added will range from 1 – 10 ul for a single well in 24-well plate)
3. Day 3: Change to complete growth medium
4. Day 6-7: Look at the cells for reporter expression if the rAAV construct has a reporter (e.g. GFP) and/or begin appropriate antibiotic selection to establish stable cell line.

Example Data and Applications

A. AAV GFP Expression Vector



B. 6 days After Virus Infection




Fig. 1. Schematic diagram of the rAAV construct used *in vitro* and *in vivo* experiments is shown in Figure 1A. Three AAV serotypes (AAV-2, AAV-DJ and AAV-5) were precipitated using AAVanced™ Concentration Reagent. Three cell lines (HEK293T, CHO and HT1080) were used to validate the transduction efficiency of virus isolated with AAVanced™ Concentration Reagent *in vitro*. 1 ul of precipitated virus was used to infect cells, and images are taken 6 days after virus infection (Figure 1B).




Fig. 2. Representative images of hippocampal and midbrain sections from 6 week old C57 mice injected with AAV virus (ITR-PGK-GFP-ITR) concentrated using AAVanced™ Concentration Reagent. 1.5 ul of concentrated AAV virus (AAV-2 serotype) was delivered to hippocampal and midbrain regions by stereotaxic injection. Three weeks after virus injection animals are perfused by paraformaldehyde, fixed, and brains are sectioned to 40 micron slices before visualization.




Fig. 3. Representative images of mouse neurons in cortical sections of 2 month old C57 mice injected with AAV virus (ITR-PGK-GFP-ITR) concentrated using AAVanced™ Concentration Reagent. 1.5 ul of concentrated AAV virus (AAV-5 serotype) is delivered to cortical region by stereotaxic injection. Two weeks after virus injection animals are perfused by paraformaldehyde, fixed and brains are sectioned to 70 micron slices before visualization.

References

1. Vasileva A, Jessberger R. Precise hit: adeno-associated virus in gene targeting. *Nat Rev Microbiol.* 2005 Nov;3(11):837-47.
2. Petrs-Silva H, Linden R. Advances in recombinant adeno-associated viral vectors for gene delivery. *Curr Gene Ther.* 2013 Oct;13(5):335-45.
3. Vandenberghe LH et al. Efficient serotype-dependent release of functional vector into the culture medium during adeno-associated virus manufacturing. *Hum Gene Ther.* 2010 Oct;21(10):1251-7.
4. Lock M et al. Rapid, simple, and versatile manufacturing of recombinant adeno-associated viral vectors at scale. *Hum Gene Ther.* 2010 Oct;21(10):1259-71
5. Doria M et al. AAV2/8 vectors purified from culture medium with a simple and rapid protocol transduce murine liver, muscle, and retina efficiently. *Hum Gene Ther Methods.* 2013 Dec;24(6):392-8.

Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site:
<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200, (888) 266-5066 (Toll Free)

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com

Technical Support: tech@systembio.com

Ordering Information: orders@systembio.com

Licensing and Warranty Statement

Limited Use License

Use of the ExoELISA-Ultra Kits (*i.e.*, the “Product”) is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

- The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.
- The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.
- This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

Limited Warranty

SBI warrants that the Product meets the specifications described in this manual. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2018 System Biosciences (SBI), All Rights Reserved



System Biosciences (SBI)
2438 Embarcadero Way
Palo Alto, CA 94303

Phone: (650) 968-2200

Toll Free: (888) 266-5066

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com

Technical Support: tech@systembio.com

Ordering Information: orders@systembio.com