

# pCDF cDNA Cloning and Expression Lentivectors

Cat. #s CD100A-1 - CD111B-1

User Manual

Store kit at -20°C on receipt

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.

(ver. 5-061226)

## **Contents**

Ι.	Inti	roduction and Background	
	Α.		2
	В.	Advantages of the Lentivector Expression System	2
	C.	pCDF cDNA Cloning and Expression Lentivectors	
	D.	List of Components	5
	E.	List of Components Additional Required Materials	5
	F.	Safety Guidelines	6
п.		otocol	
	Α.	cDNA Amplification	.8
	В.	Preparation of Digested pCDF Vectors	8
	C.	cDNA Cloning into pCDF Vectors	9
	D.	Packaging of pCDF Expression Construct	11
III.	Tro	bubleshooting	12
	Α.	Large number of colonies on control plate	12
		No or low number of colonies on plate with cDNA sample	
		No correct cDNA inserts	
	_		
IV.	Re	ferences	14
v	۸n	pendix	
۷.	AP	Map and Features for pCDF1-MCS1 Vector	17
	А. В.	Map and Features for pCDF1-MCS1 vector	10
	ь. С.	Map and Features for pCDF1-MCS2-EF1-Puro Vector	10
		Bastriation Mana of Multimla Claning Sites (MCS)	19
	D.	Restriction Maps of Multiple Cloning Sites (MCS)	20
		Properties of copGFP Fluorescent Protein	
	F.		21
	G.	Technical Support	22
vi	Lic	ensing and Warranty Statement	23

## I. Introduction and Background

#### A. Purpose of this Manual

This manual provides details and information necessary to generate expression constructs of your gene of interest in the pCDF lentivectors. Specifically, it provides critical instructions on amplification and cloning the cDNA into the pCDF Vectors, and verifying final expression constructs. This manual does not include information on packaging the pCDF expression constructs into pseudotyped viral particles or transducing your target cells of choice with these particles. This information is available in the user manual Lentivector Expression Systems: Guide to Packaging and Transduction of Target Cells, which is available on the SBI website (www.systembio.com). Before using the reagents and material supplied with this system, please read the entire manual.

#### B. Advantages of the Lentivector Expression System

Lentiviral expression vectors are the most effective vehicles for delivering and expression of a gene of interest to almost any mammalian cell—including non-dividing cells and model organisms (C.A. Machida, 2003; M. Federico, 2003; W. C. Heiser, 2004). As with standard plasmid vectors, it is possible to introduce lentivector expression constructs in plasmid form into the cells with low-to-medium efficiency using conventional transfection protocols. However, by packaging the lentivector construct into viral particles, you can obtain highly efficient transduction of expression constructs—even with the most difficult to transfect cells, such as primary, stem, and differentiated cells. The expression construct transduced in target cells is integrated into genomic DNA and provides stable, long-term expression of the target gene.

The lentiviral cDNA expression system consists of three main components:

- (1) The lentiviral expression vector (*e.g.*, pCDF1-MCS2-EF1-Puro)
- (2) The lentiviral packaging plasmids (*e.g.*, pPACKF1<sup>™</sup> Packaging Plasmid mix)
- (3) A pseudoviral particle producer cell line (*e.g.*, 293TN cells)

The expression lentivector contains the genetic elements responsible for packaging, transduction, stable integration of the viral expression construct into genomic DNA, and expression of the target gene sequence. The packaging vector provides all the proteins essential for transcription and packaging of an RNA copy of the expression construct into recombinant viral particles. To produce a high titer of viral particles, expression and packaging vectors are transiently cotransfected into producer mammalian cells (e.g., HEK 293 cells). For a detailed description of SBI's Lentivector expression system, please refer to the Lentivector Expression Systems user manual.

SBI's novel pCDF Vectors are derived from feline immunodeficiency virus (FIV; Poeschla, 2003; for Safety Guidelines when working with these vectors, see section G). These pCDF Vectors, developed at SBI, are self-inactivating as a result of a deletion in the U3 region of 3' ΔLTR (see Appendix for Vector Features). Upon integration into the genome, the 5' LTR promoter is inactivated, which prevents formation of replication-competent viral particles.

When expressed, the hybrid CMV/FIV 5' LTR drives high level transcription of the viral construct and produces a transcript that contains all the necessary functional elements (i.e., Psi, RRE, and cPPT) for efficient packaging. When this construct is expressed in HEK 293 cells that also express viral coat proteins (i.e., a packaging cell line), the pCDF transcripts are efficiently packaged into pseudoviral particles. After isolation, these pseudoviral particles containing the RNA version of the pCDF expression cassette can be efficiently transduced into any mammalian target cells. Following transduction into the target cells, this expression cassette is reverse transcribed and integrated into the genome of the target cell. The pCDF Vectors also contain a bacterial origin of replication and ampicillin resistance (Amp<sup>R</sup>) gene for propagation and selection in E. coli. The pCDF1-MCS2-EF1-Puro Vector (Cat. # CD110B-1) contains a puromycin resistance gene, under the control of a constitutive EF1 promoter and a WPRE regulatory element, to enable selection of target cells stably expressing the cDNA template. The pCDF1-MCS2-EF1-copGFP Vector (Cat. # CD111B-1) contains a copGFP gene under the control of a EF1 promoter and WPRE CopGFP is a novel fluorescent protein ,derived from element. copepod plankton (Panalina sp.), which is similar to EGFP but has a brighter color This gene serves as a reporter for the transfected or transduced cells.

#### C. pCDF Cloning and Expression Lentivectors

The FIV derived pCDF vectors contain the following features:

- CMV promoter-promotes a high level of expression of your • gene of interest in a wide variety of cell lines.
- Multiple Cloning Site (MCS)—for cloning the gene of interest in ٠ MCS located downstream of CMV promoter.
- WPRE element—enhances stability and translation of the CMV-• driven transcripts.
- SV40 polyadenylation signal—enables efficient termination of transcription and processing of recombinant transcripts.

- Optional second expression cassette—provides expression of puromycin resistance gene or copGFP reporter under control of constitutive elongation factor 1 (EF1) promoter for selection or FACS analysis of transduced cells.
- **Hybrid CMV-5LTR promoter**—provides a high level of expression of the full-length viral transcript in producer 293 cells.
- Genetic elements (cPPT, GAG, LTRs)—necessary for packaging, transducing, and stably integrating the viral expression construct into genomic DNA.
- **SV40 origin**—for stable propagation of the pCDF plasmid in mammalian cells.
- **pUC origin**—for high copy replication and maintenance of the plasmid in *E.coli* cells.
- Ampicillin resistance gene—for selection in *E.coli* cells.

#### D. List of Components

#### pCDF cDNA Cloning and Expression Lentivectors:

Component	Cono	Amount
pCDF1-MCS2-EF1-copGFP		Cat. #: CD111B-1
pCDF1-MCS2-EF1-Puro		Cat. #: CD110B-1
pCDF1-MCS1 cDNA Expression Vector	or	Cat. #: CD100A-1

Component	Conc.	Amount
pCDF cDNA Expression Vector	0.5 μg/μl	20 µg

All plasmids are shipped at a concentration of 0.5  $\mu$ g/ $\mu$ l and an amount of 20 µg. All kits are shipped in dry ice and should be stored at -20°C upon receipt. Properly stored kits are stable for 12 months from the date received.

#### E. Additional Required Materials

#### For Cloning

- Restriction enzymes for digestion of the vectors and/or inserts (Recommended: New England BioLabs enzymes)
- High Fidelity Long-distance PCR enzymes
- T4 DNA Ligase and ligation reaction buffer (Recommended: New England BioLabs T4 DNA Ligase (400 U/µl), Cat. # M0202S. Dilute to 40 U/ul with the provided 1X reaction buffer just before use)
- High efficiency competent E. coli cells (RecA) (Recommended: Invitrogen One Shot OmniMAX 2 competent cells, Cat. # C8540-03)
- Petri plates containing LB Agar media with 50 µg/ml Ampicillin

#### For Screening Inserts and Sequencing

- Tag DNA polymerase, reaction buffer, and dNTP mix . (Recommended: Clontech Titanium™ Taq DNA polymerase, Cat. # 639208)
- PCR machine
- 2-3% 1X TAE Agarose gel

#### For Purifying cDNA Constructs after Cloning

- Plasmid purification kit (Recommended: QIAGEN Endotoxin-free Plasmid Kit. The following kit combinations can be used for Midi scale (up to 200 µg of plasmid DNA) preparation of endotoxin-free DNA:
  - QIAfilter Plasmid Midi Kit, Cat. # 12243, and EndoFree Plasmid Maxi Kit. Cat. # 12362

 QIAfilter Plasmid Midi Kit, Cat. # 12243, and EndoFree Plasmid Buffer Set, Cat. # 19048

Please visit the QIAGEN website to download the specialized protocol that is not contained in the current user manual:

http://www1.qiagen.com/literature/protocols/pdf/QP15.pdf

#### For Transfection of pCDF Constructs into Target Cells

 Transfection Reagent (Recommended: Invitrogen Lipofectamine 2000, Cat. # 11668-027)

#### For Packaging of pCDF Constructs in Pseudoviral Particles

- In order to package your pCDF cDNA constructs into VSV-G pseudotyped viral particles, you will need to purchase the pPACKF1 Lentivector Packaging Kit (Cat. # LV100A-1). The protocol for packaging and transduction of packaged pseudoviral particles is provided in the User Manual for the Lentivector Expression System.
- 293 Producer Cell Line (Recommended: SBI 293TN Cell Line, Cat. # LV900A-1 or ATCC 293 Cells, Cat. # CRL-11268)
- Transfection Reagent (Recommended: Invitrogen Lipofectamine, Cat. # 18324-111 and
- Plus Reagent, Cat. # 11514-015)

#### F. Safety Guidelines

SBI's Expression lentivectors together with the pPACK packaging plasmids comprise the third-generation lentiviral expression system. The original FIV expression system was developed by Eric M. Poeschla, David J. Looney, and Flossie Wong-Staal at UCSD (Poeschla, 1998; Poeschla 2003). The feline immunodeficiency virus (FIV) was originally isolated from cat blood. Despite common close exposure of humans to FIV through contact with domestic cats (including bites, scratches, etc.), no human infection or disease has ever been associated with FIV (Poeschla, 2003).

Both FIV-based and HIV-based lentivector systems are designed to maximize their biosafety features, which include:

- A deletion in the enhancer of the U3 region of 3'∆LTR ensures self-inactivation of the lentiviral construct after transduction and integration into genomic DNA of the target cells.
- The RSV promoter (in HIV-based vectors) and CMV promoter (in FIV-based vectors) upstream of 5'LTR in the lentivector allow efficient Tat-independent production of viral RNA, reducing the number of genes from HIV-1 that are used in this system.

- Number of lentiviral genes necessary for packaging, replication • and transduction is reduced to three (gag, pol, rev), and the corresponding proteins are expressed from different plasmids (for HIV-based packaging plasmids) lacking packaging signals and share no significant homology to any of the expression lentivectors, pVSV-G expression vector, or any other vector, to prevent generation of recombinant replication-competent virus.
- None of the HIV-1 genes (gag, pol, rev) will be present in the packaged viral genome, as they are expressed from packaging plasmids lacking packaging signal-therefore, the lentiviral particles generated are replication-incompetent.
- Pseudoviral particles will carry only a copy of your expression construct.

Despite the above safety features, use of SBI's lentivectors falls within NIH Biosafety Level 2 criteria due to the potential biohazard risk of possible recombination with endogenous viral sequences to form self-replicating virus, or the possibility of insertional mutagenesis. For a description of laboratory biosafety level criteria, consult the Centers for Disease Control Office of Health and Safety Web site at http://www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4s3.htm. It is also important to check with the health and safety guidelines at your institution regarding the use of lentiviruses and always follow standard microbiological practices, which include:

- Wear gloves and lab coat all the time when conducting the • procedure.
- Always work with pseudoviral particles in a Class II laminar flow • hood.
- All procedures are performed carefully to minimize the creation • of splashes or aerosols.
- Work surfaces are decontaminated at least once a day and after • any spill of viable material.
- cultures, stocks, and other regulated wastes All are decontaminated before disposal by an approved decontamination method such as autoclaving. Materials to be decontaminated outside of the immediate laboratory area are to be placed in a durable, leakproof, properly marked (biohazard, infectious waste) container and sealed for transportation from the laboratory.
- Please keep in mind that pCDF vectors are integrated into genomic DNA and could have a risk of insertional mutagenesis.

## II. Protocol

The following section provides general guidelines for the cloning of cDNA, amplified by PCR, into pCDF vectors.

#### A. cDNA Amplification

Full-length cDNA fragments can be recloned from another plasmid or amplified by PCR. PCR-based cloning is the most convenient way for full-length cDNA cloning in pCDF vectors. The cDNA lentivector does not contain an ATG initiation codon. A translation initiation sequence must be incorporated in the insert cDNA if the cDNA fragment to be cloned does not already have an ATG codon. We also recommend including a Kozak sequence (i.e. GCCACC) before the ATG for optimal translation. For amplification of the target cDNA fragment, design a 5'-primer (containing a Kozak sequence and ATG codon) and 3'-primer with unique restriction sites present in the MCS of the pCDF vector but not present in the cDNA sequence. Amplify the cDNA fragment by high fidelity long-distance PCR using about 200 ng of plasmid template DNA and a minimum number of cycles (usually 12-15 cycles), purify, digest the amplified product with endspecific restriction enzyme(s) and purify the digested PCR product in a 1.2% agarose gel to prevent contamination with the original plasmid used for amplification.

#### B. Preparation of Digested pCDF Vector

Digest pCDF vector with the corresponding restriction enzymes used for preparation of cDNA fragments, and verify complete digestion of the vector by agarose gel electrophoresis. We suggest that you perform only preparative gel purification of the digested vector if more than one restriction enzyme is used. If you use a single restriction enzyme, dephosphorylation and gel purification of vector is necessary to reduce the background in the vector ligation step.

#### C. Cloning of cDNA into pCDF Vector

The optimal insert-to-vector molar ratio may be different for different inserts. Always try at least two different ratios (*e.g.*, 10:1 and 30:1) for each experiment. Also make sure to include one negative control reaction, which contains only the digested vector.

#### 1. Ligation of cDNA to Vector

- a. Dilute the gel-purified digested vector to 10 ng/µl.
- b. Set up 10  $\mu l$  ligation reactions for each sample and control, as follows:

1.0 μl	Digested pCDF Vector (10 ng/µl)
7.0 μl	cDNA insert (usually 15-50 ng) or Nuclease-free water
1.0 μl	10X T4 DNA Ligase Buffer
1.0 μl	T4 DNA ligase (40 U/μl)
10.0 μl	Total volume

c. Incubate the ligation reaction at 16°C for 2-5 hrs, if it is stickyend ligation. For blunt-end ligation, an overnight incubation time is recommended.

#### 2. Transform *E. coli* with the ligation product

Transform competent cells (with a transformation efficiency of at least  $1 \times 10^9$  colonies/µg pUC19) with the whole ligation reaction (10 µl) following the protocol provided with the competent cells. Plate the transformed bacteria on LB-Ampicillin agar plates.

#### 3. Identify Clones with the cDNA Insert

- a. Depending on the ratio of colony numbers for the cDNA sample vs. the negative control sample, randomly pick 5 or more well-isolated colonies and grow each clone in 100  $\mu$ l of LB Broth with 75  $\mu$ g/ml ampicillin at 37°C for 2 hours with shaking.
- b. Use 1  $\mu$ l of each bacterial culture for screening cDNA inserts by PCR and continue to grow the culture for another 4 hours. Store the culture at 4°C.

c. Prepare a PCR Master Mix with PCR primers flanking the cDNA insert:

<u>1 rxn</u>	<u>10 rxn</u>	Composition
0.5 μl	5 μl	PCR primer 1 (10 μM)
0.5 μl	5 μl	PCR primer 2 (10 μM)
0.5 μl	5 μl	50X dNTP mix (10 mM of each)
2.5 μl	25 μl	10X PCR Reaction Buffer
19.5 μl	195 μl	Deionized water
0.5 μl	5 μl	Taq DNA polymerase (approx. 5 U/μl)
24.0 μl	240 μl	Total volume

- d. Mix the master mix very well and aliquot 24  $\mu I$  into each well of 96-well PCR plate or individual tubes.
- e. Add 1  $\mu l$  of each bacterial culture from step (b) into each well (or tube).
- f. Proceed with PCR using the following program:

94°C, 4 min	1 cycle	
94°C, 0.5 min, then 68°C, 1 min/1 kb*.	25 cycles	
68°C, 3 min	1 cycle	
* depending on the size of final PCR product, use shorter or		
longer time.		


g. Take 5  $\mu$ l of the PCR reaction and run it on a 1.2% agarose/EtBr gel in 1X TAE buffer to identify clones with correct insert.

Grow a positive clone with the cDNA insert in an appropriate amount of LB-Amp Broth, and purify the construct using an endotoxin-free plasmid purification kit (see Section I.E).

Confirm identity of the cDNA insert by sequence analysis of the construct using the one of the PCR primers. Alternatively, you may use the following sequencing primer which is located upstream of the MCS: 5'-CACGCTGTTTTGACCTCCATAGA-3'.

## D. Packaging of the pCDF expression constructs into pseudoviral particles

If you are planning to create a stably transduced cell line expressing your gene of interest, you first need to package the cDNA lentiviral construct into lenti pseudoviral particles. For this purpose, you will need to purchase the pPACKF1 Lentivector Packaging Kit from SBI (see Appendix). Figure 3 schematically shows all steps which need to be performed in order to generate pseudoviral packaged cDNA expression constructs.



**Fig. 3.** Schematic presentation of the packaging procedure for lentivector expression constructs and making of stable cell lines.

The Lentivector Expression System User Manual includes the procedural information for packaging and transducing the expression constructs. This user manual is also available on the SBI web site (<u>www.systembio.com</u>). Although you can create stable transfectants with the lentiviral construct using standard transfection and selection protocols, transduction of the lentiviral cDNA construct using packaged pseudoviral particles is the most efficient way to deliver cDNA constructs in a wide range of cells, including dividing, non-dividing, and hard-to-transfect cells.

## **III.** Troubleshooting

#### A. Large number of colonies on negative control plate

If you see that the colony number on the negative control plates (no insert) is equal or more than on the plate with the cDNA sample, there is probably undigested plasmid contamination. Check your digestion conditions, and repeat digestion with an increased concentration of restriction enzyme(s) or use a longer reaction time. For best results, gel-purify and dephosphorylate the vector after single enzyme digestion. Also, check the sequences of the PCR primers in order to be sure that the necessary restriction sites are present.

#### B. No or low number of colonies on plate with cDNA sample

The efficiency of cDNA cloning in pCDF vector depends on many factors, including size, purity, integrity, modification of insert, selection of restriction sites, etc. If your cDNA sample ligation resulted in only a few colonies, please continue with PCR screening first. If none of these few colonies has the right insert, or you did not get any colonies at all, it may be caused by:

#### 1. Inappropriate ratio of insert-to-vector

Not enough or too much insert could inhibit the ligation reaction. Try a different ratio of insert-to-vector to optimize the ligation reaction. Sometimes, the yield of the ligation reaction may also be improved by increasing both the insert and vector amounts.

#### 2. Low ligation efficiency

a.	Inactive ligase and /or ligase reaction buffer	Test your ligase and reaction buffer for activity using different vector and insert. Replace the reagents if they are proven inactive.

b. Ligation inhibitors EDTA and high salt may inhibit the are present ligation reaction.

#### 3. Low transformation efficiency

a.	Low quality or poor handling of competent cells	Handle the competent cells gently. Many cells do not allow re-freezing after thawed. Quality of competent cells may be tested by transforming a circular plasmid to determine cells' competency. Use competent cells with a transformation efficiency of at least $1 \times 10^9$ colonies/ µg of pUC19 plasmid.
b.	Wrong antibiotic or too	The plates used for cloning should

b. Wrong antibiotic or too much antibiotic in the media. The plates used for cloning should contain 50-100  $\mu$ g/ml ampicillin in the media.

#### C. No correct cDNA inserts

If the colony number for the cDNA sample is more than for the negative control sample (*i.e.* vector only), but you failed to amplify cDNA insert, it could be that:

1.	Inactive Taq polymerase or reaction buffer	Test the activity of the PCR master mix by amplifying cDNA from original template. Replace the PCR reagents if they are proven inactive.
2.	Wrong primer was used	Make sure you are using the correct primers for the specific orientation of cDNA insert.
3.	Not enough clones were screened	Pick more colonies for screening.

### IV. References

Buchschacher, G.L., and Wong-Staal, F. (2000) Development of lentiviral vectors for gene therapy for human diseases. Blood. 95:2499-2504.

Burns, J.C., Friedmann, T., Driever, W., Burrascano, M., and Yee, J.K. (1993) Vesicular stomatitis virus G glycoprotein pseudotyped retroviral vectors: concentration to a very high titer and efficient gene transfer into mammalian and non-mammalian cells. Proc. Natl. Acad. Sci. USA. 90:8033-8034.

Cann, A.J.(ed). (2000) RNA Viruses. A Practical Approach. Oxford Univ. Press.

Dull, T., Zufferey, R., Kelly, M., Mandel, R.J., Nguyen, M., Trono, D., and Naldini, L. (1998) A third-generation lentivirus vector with a conditional packaging system. J. Virol. 72:8463-8471.

Gould, D.J. and Favorov, P. (2003) Vectors for the treatment of autoimmune diseases. Gene Therapy 10:912-927.

Lee, N.S., Dohjima, T., Bauer, G., Li, H., Li, M-J., Ehsani, A., Salvaterra, P., and Rossi, J. (2002) Expression of small interfering RNAs targeted against HIV-1 rev transcripts in human cells. Nature Biotechnol. 20:500-505

Morgan, R.A., Cornetta, K. and Anderson, W.F. (1990) Application of the polymerase chain reaction in retroviral-mediated gene transfer and the analysis of gene-marked human TIL cells. Hum. Gene Ther. 1:135-149.

Pfeifer, A., Kessler, T., Yang, M., Baranov, E., Kootstra, N., Cheresh, D.A., Hoffman, R.M. and Verma, I.M. (2001) Transduction of liver cells by lentiviral vectors: Analysis in living animals by fluorescence imaging. Mol. Ther. 3:319-322.

Qin, X.F., An, D.S., Chen, I.S., and Baltimore, D. (2003) Inhibiting HIV-1 infection in human T cells by lentiviral-mediated delivery of small interfering RNA against CCR5. Proc. Natl. Acad. Sci. USA 100:183-188

Quinn, T.P., and Trevor, K.T. (1997) Rapid quantitation of recombinant retrovirus produced by packaging cell clones. Biotechniques 23:1038-1044.

Sui, G., Soohoo, C. Affar, E.B., Gay, F., Forrester, W.C., and Shi, Y. (2002) A DNA vector-based RNAi technology to suppress gene expression in mammalian cells. Proc. Natl. Acad. Sci. U.S.A 99:5515-5520

Curran MA, Nolan GP. Nonprimate lentiviral vectors. Curr Top Microbiol Immunol. 2002; 261: 75-105.

Curran MA, Nolan GP. Recombinant feline immunodeficiency virus vectors. Preparation and use. Methods Mol Med. 2002; 69: 335-50

Loewen N, Barraza R, Whitwam T, Saenz DT, Kemler I, Poeschla EM. FIV Vectors. Methods Mol Biol. 2003; 229: 251-71.

Naldini L. Lentiviruses as gene transfer agents for delivery to non-dividing cells. Curr Opin Biotechnol. 1998 Oct; 9(5): 457-63.

Sauter SL. Gasmi M. FIV vector systems. Somat Cell Mol Genet. 2001 Nov: 26(1-6): 99-129.

Alisky JM, Hughes SM, Sauter SL, Jolly D, Dubensky TW Jr, Staber PD, Chiorini JA. Davidson BL. Transduction of murine cerebellar neurons with recombinant FIV and AAV5 vectors. Neuroreport. 2000 Aug 21; 11(12): 2669-73.

Brooks AI, Stein CS, Hughes SM, Heth J, McCray PM Jr, Sauter SL, Johnston JC, Cory-Slechta DA, Federoff HJ, Davidson BL. Functional correction of established central nervous system deficits in an animal model of lysosomal storage disease with feline immunodeficiency virus-based vectors. Proc Natl Acad Sci U S A. 2002 Apr 30; 99(9): 6216-21.

Crystal RG. Bad for cats, good for humans? Modified feline immunodeficiency virus for gene therapy. J Clin Invest. 1999 Dec; 104(11): 1491-3.

Curran MA, Kaiser SM, Achacoso PL, Nolan GP. Efficient transduction of nondividing cells by optimized feline immunodeficiency virus vectors. Mol Ther. 2000 Jan; 1(1): 31-8.

Derksen TA, Sauter SL, Davidson BL. Feline immunodeficiency virus vectors. Gene transfer to mouse retina following intravitreal injection. J Gene Med. 2002 Sep-Oct; 4(5): 463-9.

Haskell RE, Hughes SM, Chiorini JA, Alisky JM, Davidson BL. Viral-mediated delivery of the late-infantile neuronal ceroid lipofuscinosis gene, TPP-I to the mouse central nervous system. Gene Ther. 2003 Jan; 10(1): 34-42.

Price MA, Case SS, Carbonaro DA, Yu XJ, Petersen D, Sabo KM, Curran MA, Engel BC, Margarian H, Abkowitz JL, Nolan GP, Kohn DB, Crooks GM. Expression from second-generation feline immunodeficiency virus vectors is impaired in human hematopoietic cells. Mol Ther. 2002 Nov; 6(5): 645-52.

Stein CS, Davidson BL. Gene transfer to the brain using feline immunodeficiency virus-based lentivirus vectors. Methods Enzymol. 2002; 346: 433-54.

Browning MT, Schmidt RD, Lew KA, Rizvi TA. Primate and feline lentivirus vector RNA packaging and propagation by heterologous lentivirus virions. J Virol. 2001 Jun; 75(11): 5129-40.

Curran MA, Kaiser SM, Achacoso PL, Nolan GP. Efficient transduction of nondividing cells by optimized feline immunodeficiency virus vectors. Mol Ther. 2000 Jan; 1(1): 31-8.

Poeschla EM, Wong-Staal F, Looney DJ. Efficient transduction of nondividing human cells by feline immunodeficiency virus lentiviral vectors. Nat Med. 1998 Mar; 4(3): 354-7.

Poeschla, E.M., Looney, D.J., and Wong-Staal, F. (2003) Lentiviral nucleic acids and uses thereof. US Patent NO. 6,555,107 B2

Dull, T., Zufferey, R., Kelly, M., Mandel, R.J., Nguyen, M, Trono, D. (1998) J. Virol., 72, 8463-8471

Miyoshi, H., Blomer, U., Takashi, M., Gage, F.N., Verma, I.M (1998), J.Virol., 72, 8150-8157.

Zufferey, R., Donello, J.E., Trono, D., Hope, T.J. (1999), J.Virol., 73, 2886-2892

Ramezani, A., Hawley, T.S., Hawley, R.G. (2000) Mol. Ther., 2, 458-469

Leung, T.H., Hoffmann, A., Baltimore, D. 2004, Cell, v. 118, 453-464

Viral vectors for gene therapy. Methods and Protocols. Eds. C.A. Machida. (2003), Humana Press.

Methods in Molecular Biology. Volume 246. Gene delivery to mammalian cells. Volume 2: Viral Gene transfer techniques. Ed. by W. C. Heiser. (2004), Humana Press.


Methods in Molecular Biology. Volume 229. Lentivirus gene engineering protocols. Ed. by M. Federico. (2003), Humana Press.

Li MJ, Rossi JJ. Lentiviral vector delivery of recombinant small interfering RNA expression cassettes. Methods Enzymol. 2005;392:218-26.

Davidson BL, Harper SQ. Viral delivery of recombinant short hairpin RNAs. Methods Enzymol. 2005;392:145-73.

## V. Appendix


#### A. Map and Features for pCDF1-MCS1 Vector



Feature	Location*	Function
CMV/5'LTR	1-415	Hybrid CMV promoter-R/U5 long terminal repeat; required for viral packaging and transcription
gag	762-1011	Packaging signal
RRE	1012-1143	Rev response element binds gag and involved in packaging of viral transcripts
cPPT	1150-1391	Central polypurine tract (includes DNA Flap region) involved in nuclear translocation and integration of transduced viral genome
CMV promoter	1407-1746	Human cytomegalovirus (CMV)constitutive promoter for transcription of cloned cDNA insert
WPRE	1817-2406	Woodchuck hepatitis virus posttranscriptional regulatory elementenhances the stability of the viral transcripts
3' ΔLTR (ΔU3)	2525-2740	Required for viral reverse transcription; self- inactivating 3' LTR with deletion in U3 region prevents formation of replication-competent viral particles after integration into genomic DNA
SV40 Poly-A	2741-2872	Transcription termination and polyadenylation
SV40 Ori	2881-3027	Allows for episomal replication of plasmid in eukaryotic cells
pUC Ori	3397-4070 (C)	Allows for high-copy replication in E. coli
AmpR	4215-5075 (C)	Ampicillin resistant gene for selection of the plasmid in <i>E. coli</i>

\* The notation (C) refers to the complementary strand.


#### B. Map and Features for pCDF1-MCS2-EF1-Puro Vector



Feature	Location*	Function
CMV/5'LTR	1-415	Hybrid CMV promoter-R/U5 long terminal repeat; required for viral packaging and transcription
gag	762-1011	Packaging signal
RRE	1012-1143	Rev response element binds gag and involved in packaging of viral transcripts
cPPT	1150-1391	Central polypurine tract (includes DNA Flap region) involved in nuclear translocation and integration of transduced viral genome
CMV promoter	1407-1746	Human cytomegalovirus (CMV)constitutive promoter for transcription of cloned cDNA insert
EF1	1807-2352	Elongation factor 1α promoterconstitutive promoter for transcription of Reporter gene (Puromycin resistance or copGFP)
Puro	2358-2957	Puromycin-resistant marker for selection of the transfected/transduced cells
WPRE	2694-3553	Woodchuck hepatitis virus posttranscriptional regulatory elementenhances the stability of the viral transcripts
3' ΔLTR (ΔU3)	3674-3889	Required for viral reverse transcription; self- inactivating 3' LTR with deletion in U3 region prevents formation of replication-competent viral particles after integration into genomic DNA
SV40 Poly-A	3890-4021	Transcription termination and polyadenylation
SV40 Ori	4030-4176	Allows for episomal replication of plasmid in eukaryotic cells
pUC Ori	4546-5219 (C)	Allows for high-copy replication in E. coli
AmpR	5364-6224 (C)	Ampicillin resistant gene for selection of the plasmid in <i>E. coli</i>

\* The notation (C) refers to the complementary strand.

#### C. Map and Features for pCDF1-MCS2-EF1-copGFP Vector



Feature	Location*	Function
CMV/5'LTR	1-415	Hybrid CMV promoter-R/U5 long terminal repeat; required for viral packaging and transcription
gag	762-1011	Packaging signal
RRE	1012-1143	Rev response element binds gag and involved in packaging of viral transcripts
cPPT	1150-1391	Central polypurine tract (includes DNA Flap region) involved in nuclear translocation and integration of transduced viral genome
CMV promoter	1407-1746	Human cytomegalovirus (CMV)constitutive promoter for transcription of cloned cDNA insert
EF1	1807-2352	Elongation factor 1α promoterconstitutive promoter for transcription of Reporter gene (Puromycin resistance or copGFP)
copGFP	2366-3124	Copepod green fluorescent protein (similar to regular EGFP, but with brighter color) as a reporter for the transfected/transduced cells
WPRE	3131-3720	Woodchuck hepatitis virus posttranscriptional regulatory elementenhances the stability of the viral transcripts
3' ΔLTR (ΔU3)	3839-4054	Required for viral reverse transcription; self- inactivating 3' LTR with deletion in U3 region prevents formation of replication-competent viral particles after integration into genomic DNA
SV40 Poly-A	4055-4186	Transcription termination and polyadenylation
SV40 Ori	4195-4341	Allows for episomal replication of plasmid in eukaryotic cells
pUC Ori	4711-5384 (C)	Allows for high-copy replication in E. coli
AmpR	5529-6389 (C)	Ampicillin resistant gene for selection of the plasmid in <i>E. coli</i>

\* The notation (C) refers to the complementary strand.

#### D. Restriction Maps of Multiple Cloning Sites (MCS)

MCS1 XbaI					
1701					TCTA
	AscI		BglII	PacI	AGAT
	SmaI			SwaI	
	XbaI ~~	BamHI	Hind	lIII	EcoRI
1751	GAGCCCGGGC GC CTCGGGCCCG CG NotI EcoRI	GCCGGATC	CAGATCTAAG	CTTAATTAAT	TTAAATGAAI
1801	~~ TCGCGGCCGC AGCGCCGGCG				
MCS2					
1701					TCTA AGAT
AscI					
	SmaI		BglII ~~~~~	Swal	NotI
	XbaI ~~	BamHI	Pa	cI	EcoRI
1751	GAGCCCGGGC GC CTCGGGCCCG CG NotI	GCCGGATC	CAGATCTTAA	TTAATTTAAA	TGAATTCGCG
1801	GCCGC CGGCG				

#### E. Properties of the copGFP Fluorescent Protein

The pCDF1-MCS2-EF1-copGFP Vector contains the full-length copGFP gene with optimized human codons for high level of expression of the fluorescent protein from the CMV promoter in mammalian cells. The copGFP marker is a novel natural green monomeric GFP-like protein from copepod (Pontellina sp.). The copGFP protein is a non-toxic, non-aggregating protein with fast protein maturation, high stability at a wide range of pH (pH 4-12), and does not require any additional cofactors or substrates. The copGFP protein has very bright fluorescence that exceeds at least 1.3 times the brightness of EGFP, the widely used Aeguorea victoria GFP mutant. The copGFP protein emits green fluorescence with the following characteristics:

emission wavelength max - 502 nm; excitation wavelength max - 482 nm; quantum yield - 0.6: extinction coefficient – 70.000 M<sup>-1</sup> cm<sup>-1</sup>

Due to its exceptional properties, copGFP is an excellent fluorescent marker which can be used instead of EGFP for monitoring delivery of lentivector constructs into cells.

#### F. Related Products

pPACKF1<sup>™</sup> Lentivector Packaging Kit (Cat. # LV100A-1) •

Unique lentiviral vectors that produce all the necessary FIV viral proteins and the VSV-G envelope glycoprotein from vesicular stomatitis virus required to make active pseudoviral particles. 293TN cells (SBI, Cat. # LV900A-1) transiently transfected with the pPACKF1 and a pCDF cDNA expression construct produce packaged viral particles containing a pCDF cDNA construct.

- HIV-Based pCDH cDNA Cloning and Expression Vectors •
  - pCDH1-MCS1 (Cat. # CD500A-1)
  - $\geq$ pCDH1-MCS2 (Cat. # CD501A-1)
  - pCDH1-MCS1-EF1-Puro (Cat. # CD510A-1)  $\geq$
  - $\geq$ pCDH1-MCS1-EF1-copGFP (Cat. # CD511A-1)

#### **RNAi Cloning and Expression Lentivectors** ٠

These FIV and HIV-based single- and double-promoter shRNA and siRNA cloning vectors allow you to clone siRNA templates and efficiently transduce these siRNA constructs in a wide range of cells. For a list of currently available vectors, please visit our website at http://www.systembio.com.

#### MicroRNA Precursor Construct Collection

FIV-based microRNA Precursor Constructs allow you to express pre-miRNA, consisting of the stem loop structure and upstream and downstream flanking genomic sequence. For a list of currently available vectors, please visit our website at http://www.systembio.com.

#### • PathNet<sup>™</sup> Transcriptional Reporter Lentivectors

FIV and HIV-based transcriptional reporter vectors, allow detection of the activation of transcriptional factors (TFs) in a natural environment (nuclei). For a list of currently available vectors, please visit our website at <u>http://www.systembio.com</u>.

#### G. Technical Support

For more information about SBI products, to download manuals in PDF format, or to obtain vector sequences, please visit our web site:

http://www.systembio.com

For additional information or technical assistance, please call or email us at:

System Biosciences (SBI) 1616 North Shoreline Blvd. Mountain View, CA 94043 Phone: (650) 968-2200 (888) 266-5066 (Toll Free) Fax: (650) 968-2277 E-mail: General Information: info@systembio.com Technical Support: icch@systembio.com Ordering Information: orders@systembio.com

## VI. Licensing and Warranty Statement

#### Limited Use License

Use of the pCDF cDNA Cloning and Expression Vector (i.e., the "Product") is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms.

The purchaser of the Product is granted a limited license to use the Product under the following terms and conditions:

The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use.

The Product may not be resold, modified for resale, or used to manufacture commercial products without prior written consent of SBI.

This Product should be used in accordance with the NIH guidelines developed for recombinant DNA and genetic research.

#### **FIV Vector System**

This Product is for non-clinical research use only. Use of this Product to produce products for sale or for any diagnostic, therapeutic, clinical (including pre-clinical), veterinary or high throughput drug discovery purpose (the screening of more than 10,000 compounds per day) is prohibited. In order to obtain a license to use this product for these commercial purposes, contact The Regents of the University of California. This Product or the use of this Product is covered by U.S. Patent No. 6,555,107 owned by The Regents of the University of California.

#### WPRE Technology

System Biosciences (SBI) has a license to sell the Product containing WPRE, under the terms described below. Any use of the WPRE outside of SBI's Product or the Products' intended use, requires a license as detailed below. Before using the Product containing WPRE, please read the following license agreement. If you do not agree to be bound by its terms, contact SBI within 10 days for authorization to return the unused Product containing WPRE and to receive a full credit.

The WPRE technology is covered by patents issued to The Salk Institute for Biological Studies.

SBI grants you a non-exclusive license to use the enclosed Product containing WPRE in its entirety for its intended use. The Product containing WPRE is being transferred to you in furtherance of, and reliance on, such license. Any use of WPRE outside of SBI's Product or the Product's intended use, requires a license from the Salk Institute for Biological Studies.

This license agreement is effective until terminated. You may terminate it at any time by destroying all Products containing WPRE in your control. It will also terminate automatically if you fail to comply with the terms and conditions of the license agreement. You shall, upon termination of the license agreement, destroy all Products containing WPRE in you control, and so notify SBI in writing.

This License shall be governed in its interpretation and enforcement by the laws of California.

Contact for WPRE Licensing: The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037; Attn: Office for Technology Management; Phone: (858) 435-4100 extension 1275; Fax: (858) 450-0509.

#### CMV Promoter

The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.

#### CopGFP Reporter

This product contains a proprietary nucleic acid coding for a proprietary fluorescent protein(s) intended to be used for research purposes only. Any use of the proprietary nucleic acids other than for research use is strictly prohibited. USE IN ANY OTHER APPLICATION REQUIRES A LICENSE FROM EVROGEN. To obtain such a license, please contact Evrogen at license@evrogen.com.

SBI has pending patent applications on various features and components of the Product. For information concerning licenses for commercial use, contact SBI.

Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein.

#### Limited Warranty

SBI warrants that the Product meets the specifications described in the accompanying Product Analysis Certificate. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2006 System Biosciences (SBI).